AESCULAP®

D Service-Manual
ElanEC GA835, GA830

GB Service-Manual
ElanEC GA835, GA830



B|BRAUN


CE 0123

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Seite 1-17
Page 18-34

Legende

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Legend

- 1 Flexible shaft connection
 - 2 Flexible shaft sensor
 - 3 LED "CCW operation"
 - 4 LED "CW operation"
 - 5 "Mains Off" - Switch
 - 6 "Mains On" - Switch
 - 7 LED "Mains On"
 - 8 LED "Pump On" (GA835 only)
 - 9 Flow rate control (GA835 only)
 - 10 Lever (GA835 only)
 - 11 Coolant pump (GA835 only)
 - 12 Connection for sterile fluid bottle support bracket (GA835 only)
 - 13 Support bracket for sterile fluid bottle (GA835 only)
 - 14 Control cable socket
 - 15 Ventilation grille
 - 16 Equipotential socket
 - 17 Fuse holder including voltage selection insert
 - 18 Power plug socket
 - 19 Flexible shaft support connection
-
- A Motherboard GA835810
 - B Pump control card GA835801
 - C Display board GA835820
 - D Filter board GA835802
 - E Motor
 - F Fan
 - G End shield

Copyright © 2001

Aesculap AG & CO KG®

All rights reserved

June 2001

The information contained in this service manual is protected by copyright laws. All rights are protected.

The information contained in this service manual may not be copied or otherwise reproduced either in whole or in part without the express permission of Aesculap AG & CO. KG.

The technical information, illustrations and dimensions in this manual should not be construed as a commitment by Aesculap. No claims may be made on the basis of the information contained herein. We reserve the right to make improvements without altering this documentation.

We reserve the right to make technical changes without prior notice.!

- For service, maintenance and repairs, contact the AESCULAP/B. Braun representative in your country.
- All guarantees and warranty rights as well as licenses will be voided if modifications on the unit are performed.

Service address:

AESCULAP Technischer Service
 Am Aesculap-Platz
 D-78532 Tuttlingen
 Germany
 Tel.: +49 (7461) 95 27 00
 Fax: +49 (7461) 16 28 87
 E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Contents

1	Measuring and test equipment.....	20
2	Opening and closing the device.....	20
2.1	Opening the device.....	20
2.2	Closing the device.....	20
3	Visual inspection.....	21
3.1	Damage caused by disinfecting agents, UV radiation, etc.	21
3.1.1	Wiring	21
3.1.2	Electronic boards.....	21
3.1.3	Housing.....	21
3.2	Mechanical damage	21
4	Putting the device into operation	21
5	Selection of voltage.....	22
6	Basic function test.....	22
6.1	Preparation.....	22
6.2	Implementation of test.....	22
7	Description of technical function.....	23
7.1	Description of motor system.....	23
7.2	Mode of operation.....	23
7.3	Block diagram	24
8	Technical description of subassemblies.....	25
8.1	Motherboard GA835810.....	25
8.1.1	Component mounting diagram.....	25
8.1.2	Connector reference.....	26
8.1.3	Functions of motherboard.....	27
8.2	Pump controls GA835801.....	29
8.2.1	Component mounting diagram.....	29
8.2.2	Connector reference.....	29
8.2.3	Function.....	30
9	Changing fuses	30
10	Wiring diagram.....	31
11	Repairing the coolant pump	32
11.1	Replacing the pump springs.....	32
11.2	Replacing the pump motor.....	32
12	Final test.....	33
12.1	Measuring the protective conductor resistance....	33
12.2	Measuring the equivalent unit leakage current....	33
13	Spare parts list	34

Responsibility of manufacturer in case of repair

The manufacturer shall only be responsible for the effects on the safety, reliability and performance of the device, if

- the installation, extension, readjustment, alteration or repair is carried out by technically experienced and trained staff and
- the electrical installation of the respective room complies with the requirements of the IEC guidelines
- the device is used as specified in the Instructions for Use.
- Where any form of repair, i.e. change of printed circuit board etc. has been carried out on an opened device, the final test, described in section 12 must be carried out. Only original AESCULAP spare parts may be used for any repairs.

AESCULAP offers respective training courses, after which trained members of staff are authorized to repair the devices covered during the training course.

1 Measuring and test equipment


- | | |
|---------------------------------|--|
| 1. Oscilloscope: | 50MHz / two-channel, where possible, digital storage with respective key buttons |
| 2. Universal digital multimeter | For measuring voltage and resistance |
| 3. Safety tester | i.e. Rigel 233, GERB GM-200, Bender Safety Tester Mod. 601, etc. |

2 Opening and closing the device

Warning:

- Disconnect from mains prior to opening the device as appliance sockets, cables or other contact parts inside the device may still contain hazardous mains voltages.

2.1 Opening the device

- A) Remove the 4 screws on the base
- B) Fold back the housing cover
Caution: the earthing cable of the cover may not be excessively tensioned!
- C) Disconnect the protective conductor from the cover

2.2 Closing the device

Note:

Prior to replacing the cover:

- ensure the correct connection between the cover and the earthing cable
 - inspect the device as specified in section 3. (visual inspection)
- Pull bottom edges of device cover slightly apart and position over device.
 - Ensure that no leads are caught in between the cover and the device.
 - Replace the 4 screws on the base.

3 Visual inspection

Caution!

- **Disconnect from mains prior to opening the device as appliance sockets, cables or other contact parts inside the device may still contain hazardous mains voltages.**

Prior to carrying out the testing described below, the inside and outside of the device must be visually checked for any defects.

It must also be visually checked that all subassemblies removed during the repair are correctly secured and that the movement of moveable parts is not impeded.

When visually checking the safe condition of the device, it must be ensured that

- the fuses comply with the values specified by the manufacturer,
- any labelling and stickers on the device are legible,
- the mechanical condition allows the further secure operation,
- there is no contamination that could adversely affect the safety of the device
- no solder residue, cut-off wires, loose screws etc. are left inside the device,
- the protective earth conductors are not damaged,
- any screw connections providing the electrical connection between two housing parts, are securely fixed,
- the connecting wires are not adversely strained.

3.1 Damage caused by disinfecting agents, UV radiation, etc.

The disinfecting agents used in the operating theatre as well as the ozone released during the use of UV radiation can change the surface of the device subassemblies.

3.1.1 Wiring

Cables must be checked for their correct condition and seat. In case where inadequate insulation, breaks and inadequate contacts are suspected, the respective cables must be replaced.

3.1.2 Electronic boards

In case of any signs of corrosion or other damage to the contacts or solder points of the boards, the boards must

be replaced to ensure the safe operation of the device. After replacing the board, all plug-in connections must be carefully re-connected and the board catch must be closed again.

3.1.3 Housing

The condition of the front panel and housing must be checked. No damage that could adversely affect the function and safety of the device may be detected. The correct functioning of the control panels and LEDs must be checked.

3.2 Mechanical damage

Where serious mechanical damage has occurred during transportation or due to incorrect use, the device must be returned to the manufacturer for repair.

4 Putting the device into operation

 **Carefully read the Instructions for Use of the device prior to the first start-up. This will prevent any damage caused by incorrect installation or operation that is not covered by the guarantee.**

- The device may only be operated with an AC current of 50-60 Hz.
- The mains voltage must
 - correspond to the voltage specified on the type sign
 - correspond to the voltage specified on the appliance socket **(18)** (chpt. 5. Selection of voltage).
- The mains cable may only be plugged into a correctly installed socket.

Additional general information, such as technical data, is specified in the Instructions for Use.

5 Selection of voltage


Caution!

- The following work may only be carried out after the device has been disconnected from the mains, as contact parts may contain hazardous voltages.

The mains voltage must correspond to the mains voltage specified on the type sign.

A window (appliance socket) at the back of the control unit shows the voltage to which the device has been set. The standard pre-set voltage for European operation is: 230V.

If the specified value does not correspond to the actual mains voltage, the voltage selection must be changed.

- Release catch on fuse holder with suitable screwdriver.
- Pull out fuse holder.
- Pull out voltage selection insert.
- Replace voltage selection insert in such a way that the desired voltage range is readable.
- Insert fuse holder until audible engagement of catch.

6 Basic function test

6.1 Preparation

See also Instructions for Use

Ensure that the accessories used for the test function correctly.

- Plug cable of the foot controls into the control cable socket (14) at the back of the device.
- Place hose arrangement in coolant pump (11) (GA835 only).
- Connect device to mains.

6.2 Implementation of test

- Switch on the device with the mains ON (6) switch
 - The mains ON (7) display is illuminated.
- Push down foot controls up to stop.
 - The motor now runs at maximum speed.

Select flow rate of coolant pump (GA835 only)

Note

- The coolant pump (11) can only be switched on and off via the foot controls.
 - Upon switch on of the pump, also the pump display (8) is activated.
 - The coolant pump is only started by activating the motor.
 - After releasing the foot control pedal, the pump briefly reverses the direction of rotation (approx. ½ revolution), to prevent any further dripping of the coolant.
- To increase flow rate: turn button (9) clockwise.
 - To decrease flow rate: turn button (9) anticlockwise
 - The pump output is infinitely variable.

7 Description of technical function

General information

The functional description explains the functions of the individual circuit subassemblies and their cooperation.

As far as possible, complete functional units have always been positioned on the same printed board. The block diagram is mainly based on these functional units, so that there is always a direct connection between a block of the block diagram and the individual printed boards.

The aim of this technical description is to introduce service technicians authorized by the manufacturer to the technical aspects and functional principles of this device and to enable them to carry out certain repairs by themselves, where required. It is, however, not intended to provide such in-depth information to allow the repair of subassemblies. Instead, a detailed overall description and specific measuring instructions should allow technicians to recognize and replace defective subassemblies or reset them, where necessary.

7.1 Description of motor system

- | | |
|-------------------|--|
| a) Control unit | GA830 version w/o pump
GA835 version with pump |
| b) Foot pedal | GA136, GA139, GA188, GD411 |
| c) Flexible shaft | Macro: GA156, GA172, GB400, GB401
Micro: GA173, GA176 |

7.2 Mode of operation

The electric surgical motor GA835 (GA830), with clockwise and counterclockwise torque, transfers the motor output to the different AESCULAP handpieces via the flexible shaft. With the aid of the foot pedal (i.e. GA188) the control electronics allow an infinitely variable speed control of up to max. 20000 rpm. The speed can be precisely adjusted and maintained at a constant level even under extreme loading.

The electric surgical motor GA835 (GA830) contains an flexible shaft sensor.

Where a macro-flexible shaft is used, the motor automatically switches to clockwise torque and at the same time limits the motor speed to a maximum of 15000 rpm.

The motor cannot be run with counterclockwise torque.

Where an micro-flexible shaft is used, the motor can be operated with both clockwise and counterclockwise torque at a maximum speed of 20000 rpm. An acoustic signal sounds when the motor is run with counterclockwise torque.

Note:

The micro-flexible shaft must contain an flexible shaft sensor for the aforementioned work, i.e. flexible shafts (GA173 to SN: 06984; GA176 to SN: 13373) must be equipped with a new sleeve (for uncoupling the flexible shaft from the device). The new sleeve has a circumferential notch in the first third of the sleeve. The sleeve can be ordered as item number GA183R.

In case of micro-flexible shafts without a sensor, the motor automatically switches to clockwise torque and limits the maximum speed to 15000 rpm – counterclockwise operation is then blocked.


Flexible shafts GA173 are fitted as standard from SN: 06985 onwards and GA176 from SN: 13373 onwards.

Upon starting the unit, a cooling fan is activated.

The electric surgical motor GA835 contains an integral cooling/irrigation system that can be activated/deactivated via the foot controls. The flow rate of the irrigation fluid can be infinitely varied with rotary button (9). The pump does, however, only run when the motor has been activated.

The device can be operated at voltages of 100-120V~ and 220-240V~ (+/- 10%).


7.3 Block diagram

8 Technical description of subassemblies

8.1 Motherboard GA835810

8.1.1 Component mounting diagram

Mains potential assigned

Mains potential not assigned

8.1.2 Connector reference

Connector P1:

PIN	Signal name
1	MAINS 10
2	MAINS 6
3	MAINS 9
4	MAINS 7
5	MAINS 8
6	MAINS 5

Connector P2:

PIN	Signal name
1	MAINS 3
2	MAINS 4
3	NS14
4	NS44

Connector P3:

Isolating connector

Connector P4:

PIN	Signal name
1	+22V
2	PUMP1/0
3	PUMPSTART
4	+12V
5	0V

Connector P5:

PIN	Signal name
1	U-REF
2	MP22
3	F_OUT
4	MP24
5	MP25
6	MP26
7	MP27
8	MP28
9	MP29
10	MP30

Connector P6:

PIN	Signal name
1	MP31
2	MP32

3	MP33
4	MP34
5	POT. S
6	BREAK
7	MP37
8	MP38
9	MP39
10	START

Connector P7:

PIN	Signal name
1	0V
2	PUMP1/0
3	RE/LI
4	FOOT POT. S
5	FOOT POT. E

Connector P8:

PIN	Signal name
1	+12V
2	LEDP-K
3	LEDL-K
4	LEDR-K

Connector P9:

PIN	Signal name
1	VC
2	GND
3	C
4	B
5	A

Connector P10:

PIN	Signal name
1	+9V
2	HAL1
3	HAL2
4	0V

Connector P11:

PIN	Signal name
1	MAINS LED A
2	GND

Connector P12:

PIN	Signal name
1	0V
2	+22V
3	MP 3

4	+12V
5	+9V
6	RE/LI
7	FOOT POT. E
8	MP 8
9	MP 9
10	MP10

Connector P13:

PIN	Signal name
1	MP11
2	MP12
3	R/L-K
4	ENABLE A
5	MP15
6	MP16
7	MP17
8	MP18
9	MP19
10	NC

Connector P14:

PIN	Signal name
1	GND (TP21)
2	U _{in} (TP22)
3	U _{out} (TP23)
4	+12V-N (TP24)

Connector P15:

PIN	Signal name
1	+12V
2	0V

Connector P16:

PIN	Signal name
1	PHASE A
2	NC
3	PHASE C
4	NC
5	PHASE B
6	NC
7	TEMP
8	GND

8.1.3 Functions of motherboard

The motherboard is divided into two physical and potential areas.

1.) **Non-live part** incl. the following circuit components:

- a.) Power supply (Tr1, F1, GL2, U1)
- b.) Control of "Clockwise", "Counterclockwise" and "Pump on" displays (U5)
- c.) Circuit components for conditioning of foot controls signals (U24, U3, U4, U6)
- d.) Circuit for flexible shaft detection (U7, U8)
- e.) Optical transmission of control signals to the live side (OK1 - OK4)

re. a) **Power supply to non-live part of motherboard**

The rectified secondary voltage from Tr1 serves to supply the pump and is forwarded to the switching controller U1 for generating a stabilized direct voltage of +12V.

The secondary circuit from Tr1 is protected with fuse F1 (T0.63A).

Caution!

Electric shock hazard!

In case of measuring being carried out on an open device, the device must be operated with an isolating transformer.

Checking the mains part:

Measuring

points: P12-1: reference potential 0V
P12-2 : approx. 32V - 36V
P12-3: square wave voltage U1
P12-4: +12V

2.) **Live part** incl. the following circuit components:

- a) Mains input filter (C1 - C6, L1)
- b) Power supply of controls and regulating unit (TR2, F2, GL3, U2)
- c) Power supply of motor (RT1, GL1, C7, C8)
- d) Control part (U16, U17, U18, U19, U20, U22, U23)
- e) Regulating part (U10, U9)
- f) Driver stage (U11, U12, U13 - U15)
- g) Power stage (Q11A - Q11F, D29 - D40)

re. b) Power supply of controls and regulating unit

The rectified secondary voltage from Tr2 is forwarded to the switching controller U2 in order to generate a stabilized direct voltage of 12V.

The secondary circuit from Tr1 is protected with fuse F2 (T0.2A).

FETs, a value of approx. $2M\Omega$ or higher should be measured.

The resistance of failed Fets is considerably lower ($<100\Omega$).

 **Caution!****Electric shock hazard!**

In case of measuring being carried out on an open device, the device must be operated with an isolating transformer.

Checking the mains part

Measuring points: P14-1: reference potential GND
 P14-2 : approx. 30V - 35V
 P14-3: square wave voltage U2
 P14-4: +12V

re. c) Power supply of motor:

RT1 limits the starting current. Depending on the setting of the voltage setting device in the combined net input element, the mains voltage is rectified (set to 230V) or rectified and doubled (set to 110V). The motor will thus have a constant supply of approx. 200-250V (DC).

 **Caution!****Electric shock hazard!**

In case of measuring being carried out on an open device, the device must be operated with an isolating transformer.

Checking the power supply to the motor

Measuring points: P14-1: reference potential Gnd
 D - Q11A: 250 - 330V (DC)
 Measurement taken from stopped motor.
 (Q11A = left MOSFET on heat sink
 D = middle pin of Q11A)


re. g) Power stageChecking the power stages - FETs

Prior to measuring, the device must be disconnected from the mains!

For a quick check of the functioning of the FET, the resistance between the drain and source connections can be measured, using an ohmmeter. For correctly working

8.2 Pump controls GA835801

8.2.1 Component mounting diagram

8.2.2 Connector reference

Connector P1:

PIN	Signal name
1	+22V
2	PUMP 1/0
3	PUMPSTART
4	+12V
5	0V

Connector P2:

PIN	Signal name
1	TP1
2	+22V-S
3	TP3
4	TP4
5	TP5
6	TP6
7	TP7
8	PPOT. S
9	TP9

10	TP10
----	------

Connector P3:

PIN	Signal name
1	PPOT. E
2	PPOT. S
3	0V
4	MOTOR 1
5	MOTOR 2

8.2.3 Function

The coolant pump is controlled and regulated by card GA835801.

The entire board is activated by the "Pump on" switch on the foot controller **AND** the pedal activation for controlling the EC motor.

Caution!

Electric shock hazard!

In case of measuring being carried out on an open device, the device must be operated with an isolating transformer.

Functional test:

Measuring points: P12-1 (on motherboard): reference potential 0V
P2-1 (GA835801): approx. 12V when pump is activated

When the pump is activated, the pump output is infinitely variable with the potentiometer on the front of the device.

Functional test:

Measuring points: P12-1 (on motherboard): reference potential 0V
P2-9 (GA835801): 0V to 21V (depending on potentiometer setting)


Functional test of pump motor:

Remove connector P3
Connect motor to separate laboratory power supply.
Voltage: max. 24V

9 Changing fuses

Caution!

- The following work may only be carried out after the device has been disconnected from the mains, as contact parts may contain hazardous voltages.

Note

Always change both fuses at the same time.


Prescribed fuses:

2 off T4.0A breaking capacity L.

These fuses can be ordered from the manufacturer as item number TA020436.


- Release catch on fuse holder with suitable screwdriver.
- Pull out fuse holder.
- Replace both fuses
- Insert fuse holder until audible engagement of catch.

10 Wiring diagram


11 Repairing the coolant pump

11.1 Replacing the pump springs


- Remove caps.
- Remove the screws beneath the caps.
- Pull upper part of pump towards you and off.
- Pull pump rotor towards you and off as well.

- Replace springs
The springs can be ordered from the manufacturer as article no. GD412212.
- To reassemble the pump, follow the procedure above in the reverse order.

11.2 Replacing the pump motor

- Remove upper part of pump and rotor, see Replacing the pump springs.

- Remove caps.
- Position lever as shown in the illustration above.
- Remove screws
- Pull motor away from you.
- Unsolder pump cable (can be ordered from the manufacturer as article no. GA835825).
- Solder pump cable to replacement motor GD412221, see Wiring diagram.
- To mount replacement motor and pump, follow procedure above in reverse order. In doing this, make certain that the washer and seal between the pump housing and pump motor are seated properly.

12 Final test


Required test devices: safety tester
 i.e.. RIGEL Safety Tester 233
 GERB Safety Tester GM-200

Where any repair work such as replacement of printed boards, etc. was carried out on an opened device, the below final test must be carried out.

The device has been high-voltage tested by the manufacturer. This test may not be repeated by the service technician.

12.1 Measuring the protective conductor resistance


Measurement of the protective conductor resistance
 Maximum value: 0.2 Ohm

Circuit for testing the protective conductor for units of protection class I

12.2 Measuring the equivalent unit leakage current

Measurement of the equivalent unit leakage current
 Maximum value: 0.75 mA

Circuit for testing the equivalent unit leakage current for units of protection class I

When carrying out the electrical safety test according to VDE0751 part 1 10/90, the first measured value of the discharge current of the replacement device must be recorded. During each subsequent measurement, the first measured value may not be exceeded by 10µA or 50%. In case of higher values, a non-permissible change must be suspected. In this case, the device must be returned to the manufacturer for testing.

The safety evaluation must be carried out in accordance with IEC60601-1 requirements. The defined measuring equipment must be used for this purpose. Alternatively, other measuring arrangements, also specified by IEC60601-1 may be used.

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