EN  
Konelab™ / T Series  
**GLUCOSE (GOD-POD)**

REF	981379	10 x 60 ml
	981780	12 x 20 ml

**THIS PACKAGE INSERT IS APPLICABLE FOR USE OUTSIDE THE US. ANY REFERENCE TO THE KONELAB SYSTEMS ALSO REFERS TO THE T SERIES.**

**INTENDED USE**

For the *in vitro* quantitative determination of the glucose concentration in human serum or plasma on Konelab analyzers. All test results must be interpreted with regard to the clinical context.


**SUMMARY (1, 2)**

Glucose is the primary energy source for the human body. It is derived from the breakdown of carbohydrates in the diet and in body stores. When energy intake exceeds expenditure, the excess is converted to fat and glycogen for storage in adipose tissue and liver or muscle. When energy expenditure exceeds calorie intake, endogenous glucose formation occurs from the breakdown of carbohydrate stores and from noncarbohydrate sources.

The glucose level in the blood is maintained within a narrow range under diverse condition by regulatory hormones such as insulin, glucagon or epinephrine. The most frequently encountered disorder of carbohydrate metabolism is high blood glucose due to diabetes mellitus. Other applications are the detection of hypoglycemia in the neonate, in childhood, and due to extrapancreatic tumors.

**PRINCIPLE OF THE PROCEDURE**

This method employs glucose oxidase (GOD) and a modified Trinder colour reaction, catalysed by the enzyme peroxidase (POD). Glucose is oxidised to D-gluconate by glucose oxidase with the formation of an equimolar amount of hydrogen peroxide. In the presence of peroxidase, 4-aminoantipyrine and phenol are oxidatively coupled by hydrogen peroxide to form a quinoneimine dye, coloured in red. The intensity of colour in the reaction is measured at 510 nm and it is proportional to the glucose concentration in the sample.

**REAGENT INFORMATION**

10 x 60 ml  
12 x 20 ml

**Concentrations**

Phosphate buffer (pH 7.5)	250 mmol/l
Phenol	5 mmol/l
4-aminoantipyrine	0.5 mmol/l
Glucose oxidase	≥ 15 KU/l
Peroxidase	≥ 1 KU/l
NaN <sub>3</sub>	< 0.1 %

**Precautions**

For *in vitro* diagnostic use only. Exercise the normal precautions required for handling all laboratory reagents. The reagent contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.

**Preparation**

The reagent is ready for use.

**Note:** Check that there are no bubbles in the bottleneck or on the surface of the reagent when you insert the reagent vials or vessels in the Konelab analyzer.

**Storage and Stability**

Reagents in unopened vials are stable at 2...8 °C until the expiration date printed on the label. Keep away from sunlight. Refer to the Application Notes of your Konelab analyzer for the on board stability of reagents.

**SPECIMEN COLLECTION****Sample Type**

Unhemolyzed serum or heparin or EDTA plasma can be used. The sample should be separated from the cells as soon as possible after collection in order to avoid glycolysis. If the sample is not separated or analyzed without delay a glycolytic inhibitor should be used (3).

**Precautions**

Human samples should be handled and disposed of as if they were potentially infectious.

**Storage**

The separated, nonhemolyzed sterile sample can be stored for 8 hours at 25 °C or for 3 days at 4 °C (1). The stabilized sample can be stored for 1 day at 20...25 °C or for 7 days at 4...8 °C (4).

**TEST PROCEDURE**

Refer to the Reference Manual and Application Notes for an automated procedure on your Konelab analyzer. Any application which has not been validated by Thermo Fisher Scientific Oy cannot be performance guaranteed and therefore must be evaluated by the user.

**Materials provided**

Reagent as described above.

**Materials required but not provided**

Calibrator and controls as indicated below.

**Calibration**

Use sCal, code 981831, according to the instructions provided for your Konelab analyzer.

**Traceability:**

Refer to the package insert of sCal.

**Quality Control**

Use quality control samples at least once a day and after each calibration and every time a new bottle of reagent is used.

**Available controls:**


Nortrol, code 981043

Abtrol, code 981044

The Control intervals and limits must be adapted to the individual laboratory requirements. The results of the quality control samples should fall within the limits pre-set by the laboratory.

**CALCULATION OF RESULTS**

The results are calculated automatically by the Konelab analyzer using a calibration curve.

**Calibration Curve (example)**

Konelab 20/30/60. The calibration curve is lot dependent.

**LIMITATIONS OF THE PROCEDURE****Interference**

Criterion: Recovery within ± 10% of initial values.

Bilirubin: No interference found up to 1000 µmol/l (58 mg/dl).

Hemolysate: No interference found up to 12 g/l of hemoglobin.

Lipemia: No interference found up to 7 g/l of Intralipid® (trademark of Fresenius Kabi AB).

There is a poor correlation between turbidity and triglycerides concentration.

For other interfering substances, please refer to the reference 5.

**EXPECTED VALUES (1)**

Serum (fasting):

Adults: 4.1 – 5.9 mmol/l (74 – 106 mg/dl)

The quoted values should serve as a guide only. It is recommended that each laboratory verify this range or derives a reference interval for the population that it serves.

**MEASURING RANGE**

0.3 – 20.0 mmol/l (5 – 360 mg/dl)

Extended measuring range after secondary dilution:

0.3 – 60.0 mmol/l (5 – 1080 mg/dl).

**PERFORMANCE CHARACTERISTICS****Detection limit**

0.1 mmol/l (2 mg/dl)

The detection limit represents the lowest measurable concentration/activity that can be distinguished from zero. It is calculated as the concentration of zero sample + 3 SD (within run, n=24).

**Imprecision (result unit mmol/l)**

	Mean 5.5 mmol/l		Mean 14.9 mmol/l	
	SD	CV%	SD	CV%
Within run	0.07	1.3	0.14	0.9
Between day	0.10	1.8	0.28	1.9
Total	0.12	2.2	0.35	2.3

A precision study was performed using the NCCLS Document EP5-A as a guideline and Konelab 30 during 10 days, with the number of measurements being n = 40

**Method comparison**

A comparison study was performed using the NCCLS Document EP9-A as a guideline and the Konelab glucose hexokinase method (981304) as a reference.

Linear regression (result unit mmol/l):

$$y = 1.00x - 0.06$$

$$r = 0.999$$

$$n = 75$$

The sample concentrations were between 4.1 and 17.4 mmol/l.

The results obtained in individual laboratories may differ from the given performance data.

**BIBLIOGRAPHY**

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5<sup>th</sup> edition, W B Saunders Company, Philadelphia, 2001, p. 428, 444, 983.
- Thomas, L. (ed.), Clinical Laboratory Diagnostics; Use and Assessment of Clinical Laboratory Results, 1<sup>st</sup> edition, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main, Germany, 1998, pp. 131-139.
- Väisänen S, Eskelinen S and Halonen T, Glukoosin säilyvyys näytteenottoptukissa, Kliin Lab, 3-4, 2002, pp. 48 – 50.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 3-349 – 3-367, 2000.

**MANUFACTURED BY**

Thermo Fisher Scientific Oy  
Clinical Diagnostics Finland  
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland  
Tel. +358 9 329 100, Fax +358 9 3291 0300  
www.thermo.com/konelab

**Date of revision (yyyy-mm-dd)**  
2007-08-09

**Changes from previous version**  
Company name updated.

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