EpiAccess Needle

Instructions For Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

EPIEP, INC
195 Church Street
17th Floor
New Haven, CT 06510
USA
Telephone: (434) 951-0444
Fax: (434) 951-0478

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For any questions related to the EpiAccess Needle, please contact EpiEP Customer Service Department.

Address: 195 Church Street
17th Floor
New Haven, CT 06510
USA
Office: 434-951-0444
Fax: 434-951-0478

MASTER TABLE OF CONTENTS

1.	EPIACCESS NEEDLE INSTRUCTIONS FOR USE	1
2.	AGO EPIACCESS INSTUZIONI PER L'USO	18
3.	AIGUILLE EIPACCESS MODE D'EMPLOI	34
4.	JEHLA EPIACCESS POKYNY PROPOUŽITÍ	50
5.	EPIACCESS TÚ HASZNÁLATI UTASÍTÁS	66
6.	AGULHA EPIACCESS INSTRUÇÕES DE UTILIZAÇÃO	82
7.	GEBRAUCHSANLEITUNG FÜR DIE EPIACCESS NADEL	98

TABLE OF CONTENTS

1.	INTRODUCTION	5
2.	INDICATIONS FOR USE	5
3.	CONTRAINDICATIONS	5
4.	SYMBOL LEGEND	6
5.	WARNINGS AND PRECAUTIONS	7
6.	POTENTIAL COMPLICATIONS	8
7.	SYSTEM DESCRIPTION	9
8.	EPIACCESS NEEDLE CONTENTS	11
9.	PREPARING THE SYSTEM FOR USE	12
10.	INSTRUCTION FOR USE	13
11.	DISPOSAL	17
12.	LIMITED WARRANTY	17

INTRODUCTION

The EpiAccess System has two components: the EpiAccess Control Unit (L0006 – A0002 Instructions For Use with L0014 EpiAccess System Quick Reference Guide or L0038 – A0005 Instructions for Use with L0040 EpiAccess System Quick Reference Guide), which is used in conjunction with the EpiAccess Needle. The EpiAccess Needle has an integrated fiber optic sensor housed within a stainless steel tube in the inner lumen of the needle. The Control Unit contains instrumentation that captures the pressure signal from the sensor and delivers that signal to a computer. The components of the EpiAccess Control Unit are housed in a mobile cart that allows the user to position the system as desired.

INDICATIONS FOR USE

The EpiAccess System with introducer needle and integrated needle tip pressure transducer is intended to access the epicardial surface of the heart via a subxiphoid approach to facilitate guidewire placement in electrophysiology procedures in adult patients.

DESCRIPTION

The EpiAccess Needle is placed under fluoroscopic imaging. The system provides pressure measurements taken at the needle tip and via a commercially available (not provided) A-line catheter with pressure transducer that is also connected to the EpiAccess Control Unit. These pressure measurements are displayed on the EpiEP Control Unit on a digital touchscreen display monitor, with a custom, ergonomic graphical user interface (GUI).

Changes and differences in pressure measurements within the human body are known to be indicative of specific anatomies and therefore can be used by the physician to identify the location of the device tip within the body during the electrophysiology procedure. The pressure measurements are adjunctive information provided for user convenience.

CONTRAINDICATIONS

Not intended for use in patients with any of the following conditions:

- Congenital absence of the pericardium
- Absence of a free pericardial space

SYMBOL LEGEND

	Manufacturer
	Catalog number
	Batch code
	Date of manufacture
	Use by date
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Do not re-sterilize
	Do not use if package is damaged
	Keep away from rain
	Sterilized using ethylene oxide
	Consult Instructions For Use
	Temperature limitation
	Non-pyrogenic
	Do not re-use
	General Warning Sign
	Caution

Type CF Applied Part

WARNINGS AND PRECAUTIONS

Read and follow all instructions and labeling. Not doing so can cause damage to equipment and injury to user or patient. Also, not following instructions may prevent the device from being able to perform its intended actions.

Warning! Warning statements describe conditions or actions that can result in personal injury or loss of life.

Caution! Caution statements describe conditions or actions that can result in damage to the equipment or software.

- The EpiAccess Needle is intended for single patient use only. DO NOT REUSE or RE-STERILIZE. Reuse or re-sterilization may adversely affect the functional and structural performance of the device.
- Prior to use, carefully examine the EpiAccess Needle components and verify they, and their packaging, have not been damaged during shipment. If the components show any sign of damage DO NOT USE.
- Due to the risk of exposure to HIV or other blood borne pathogens, health care workers should always use standard blood and body fluid precautions. Sterile techniques should be strictly adhered to during any handling of the device.
- Fluoroscopy should be used to confirm positioning throughout the procedure.
- Subxiphoid procedures should be performed only in facilities appropriately equipped and staffed to perform such procedures. Lab capabilities should include, but are not limited to:
 - Intracardiac pressure monitoring capabilities
 - Systemic pressure monitoring
 - Contrast media injection, and management of untoward reactions to contrast
 - Pericardiocentesis
 - Surgical backup
- The EpiAccess Needle is not designed for use in the coronary, cerebral, and pulmonary vasculature, or diseased and atherosclerotic arteries.
- During insertion, use caution not to create excessive bends in this device.
- All fluid infusion should be through the proximal needle port.

-
- Do not attempt to use a guidewire of a size larger than 0.032 inches in diameter.
 - If resistance is met when advancing or withdrawing the guidewire, determine the cause and correct before continuing with this procedure.
 - Aspirate slowly, only from the proximal needle port.
 - Inject or saline flush only from the proximal needle port.
 - Individual patient anatomy and physician technique may require procedural variations.
 - Do not use after the expiration date.
 - For specific details in the use of electrophysiology catheters and the techniques employed in an electrophysiology study, the physician should follow the respective device Instructions for Use provided with the devices/equipment, be referred to the medical literature, and rely on training and practical experience.

POTENTIAL COMPLICATIONS

Potential complications that might be associated with the EpiAccess Needle are similar to those associated with any interventional procedure; particularly those accessing the pericardium from a subxiphoid approach, and include the following:

- Bleeding from the site of access
- Vascular rupture and perforation
- Hematoma
- Embolization
- Pain
- Epicardial irritation
- Local and systemic infection
- Neurological deficits including stroke and death
- Urticaria or ulceration may occur at the site of injection
- Thrombus formation
- Cardiac perforation
- Cardiac tamponade
- Pericarditis
- Hemopericardium
- Esophageal injury
- Coronary artery injury
- Abdominal bleeding
- Pneumopericardium
- Atrial fibrillation (AF)
- Ventricular tachycardia (VT) requiring cardioversion
- Ventricular fibrillation (VF)

SYSTEM DESCRIPTION

This manual and the equipment it describes are for use only by qualified medical professionals trained in the particular technique and surgical procedure to be performed.

Warning!

Please read all information carefully. Failure to properly follow the instructions may lead to patient injury.

Note:

This manual is designed to provide instructions for use of the EpiAccess Needle with the EpiAccess Control Unit.

HOW SUPPLIED

Unit Box Content – 1 EpiAccess Needle and 1 IFU

The EpiAccess needle is supplied sterile (Ethylene Oxide), non-pyrogenic, for single-patient use disposable (short-term use). It is supplied in appropriate protective packaging, consisting of a pouched tray in a SBS carton (box), to maintain sterility and product integrity. Do not use if the package is opened or damaged.

The EpiAccess Needle is used in conjunction with the EpiAccess Control Unit.

- **EPIACCESS NEEDLE**

This system is comprised of an EpiAccess Needle, which is a Tuohy Needle with a fiber optic pressure sensor housed within a stainless steel tube, integrated inside the lumen of the access needle. The instrumentation and usage are described in this manual, the EpiAccess Needle Instructions for Use (L0007 for A0002 or L0038 for A0005).

- **EpiAccess Control Unit**

Auxiliary instrumentation includes instrumentation, which captures the pressure signal from the sensor and displays the signals on a monitor. The data is displayed in a manner, which provides guidance to the physician on the position of the needle. The instrumentation and usage are described in the EpiAccess Control Unit Instructions for Use (L0006 for A0002 or L0038 for A0005) and the EpiAccess System Quick Reference Guide (L0014 for A0002 or L0040 for A0005).

TECHNOLOGICAL CHARACTERISTICS

The EpiAccess System is composed an EpiAccess Needle (0399-13950) and an EpiAccess Control Unit (A0002 or A0005).

EpiAccess Needle

The EpiAccess Needle consists of a luer hub and a 14 gauge, 12.6 cm stainless steel cannula with a slightly curved tip similar to a Tuohy Needle. The needle configuration has a slightly larger diameter (.078 inch diameter vs. .058 inch diameter, respectively) to accommodate the integral needle tip pressure sensor. The EpiAccess Needle has an 8°-tip angle (as with Tuohy Needles).

The needle has an integrated, biocompatible fiber optic sensor housed within a stainless steel tube, which is welded inside its lumen (Figure 1). The EpiAccess Needle allows a 0.032” diameter guide wire to pass through the entire length of the needle. The EpiAccess Needle will be provided as single-patient use, sterile (EO), disposable (short-term use) device. It is supplied in appropriate protective packaging, consisting of a pouched tray in a SBS carton (box), to maintain sterility and product integrity. The EpiAccess Needle is used in conjunction with the EpiAccess Control Unit.

Needle Sensor Characteristics

Characteristic	Value
Pressure Range	-300 mmHg to 300 mmHg
Accuracy	1.5% Full Scale
Operating Temperature	10°C to 50°C
Resolution	<0.3 mmHg

EpiAccess Needle Contents

The EpiAccess Needle is comprised of one (1) 2 mm needle assembly:


Figure - 1 EpiAccess Needle

Preparing the System for Use

1. The instrumentation and usage are described in the EpiAccess Control Unit Instructions for Use:
 - a. L0006 for the A0002 EpiAccess Control Unit
 - b. L0038 for the A0005 EpiAccess Control Unit
2. When using the A0002 EpiAccess Control Unit, follow the instructions listed below under Figure 2 (Steps 4 through 6).
3. When using the A0005 EpiAccess Control Unit, follow the instructions listed below under Figure 3.

Figure 2 –Signal Conditioner Box – Front View

4. Turn the Power switch on the Signal Conditioner Box to the “ON” position.

Warning! The indicator on the front of the ISO Single Channel High Resolution Signal Conditioner MUST be Green in color prior to proceeding. Failure to ensure that the Indicator light is green may place the patient at risk of injury. IF the Indicator light is not green, replace the EpiAccess Needle before proceeding.

5. Remove the protective cover from the EpiAccess Needle Port. Retain cover for replacement after use.
6. Carefully remove the protective end cap from the EpiAccess Needle.
7. Plug in EpiAccess Needle into EpiAccess Needle Port by gently inserting into the EpiAccess Needle Port.
8. Verify that the indicator light, located on the front of the FISO Single Channel High Resolution Signal Conditioner is green. Should the light not be green, remove the EpiAccess Needle and turn the power switch to the “OFF” position. Repeat process, starting at Step 1 above.


Figure 1 - A0005 EpiAccess Control Unit

1. Remove the protective cover from the EpiAccess Needle Port on the front of the A0005 Control Unit. Retain for replacement after use.
2. Carefully remove the protective end cap from the EpiAccess Needle.
3. Plug the EpiAccess Needle into the EpiAccess Needle Port by gently inserting into the EpiAccess Needle Port.

Warning! The Needle Indicator Light on the front of the EpiAccess Control Unit MUST be Green in color prior to proceeding. Failure to ensure that the Indicator light is green may place the patient at risk of injury. If the Indicator light is not green, replace the EpiAccess Needle before proceeding.

4. Verify that the Needle Indicator Light, located above the “Start/Stop” switch is green. Should the light not be green, remove the EpiAccess Needle and replace with another EpiAccess Needle, starting with Step 1 above.

Other Installation Aspects

The EpiAccess Control Unit has the ability to replicate the installed monitor to a secondary monitor should you desire to utilize a second monitor for easier viewing during the procedure. Take note of the installation steps regarding the secondary monitor in the EpiAccess Control Unit Instructions for Use (L0006 for A0002 or L0038 for A0005).

INSTRUCTIONS FOR USE

Warning! Should the EpiAccess System fail at any time during the following procedure, revert to fluoroscopy to complete the procedure.

Note: Variations may occur within these steps, depending on available capabilities, patient morphology/condition, and operator preference.

1. Flush the needle with sterile heparinized saline.
2. Prepare and drape the selected access site as required per hospital procedures.
3. Administer sedation and general anesthesia may be considered.

Caution! Failure to make the incision as described in Step 4 below could result in skin accumulation in the EpiAccess Needle as it is advanced, resulting in an occluded tip of the EpiAccess Needle.

Note: For ablation procedures, it is best to first position the intracardiac catheters as these may act as anatomical landmarks of the right ventricular (RV) apex anteriorly and coronary sinus posteriorly.

4. Make a small incision approximately 2 cm below the subxiphoid process.
5. Under fluoroscopic guidance, direct the EpiAccess Needle posteriorly toward the left shoulder (Figure 4).


Figure 4 - Subxiphoid Approach

6. After crossing the subcutaneous tissue, monitor the needle movement under fluoroscopy in the left anterior oblique view, 35-40°.
7. The needle should always point to the left shoulder, and it should be introduced more horizontally if the target is the anterior portion of the ventricles and more vertically if the diaphragmatic portion of the heart is the area of interest.
8. Carefully move the needle toward the heart silhouette until tactilely detecting the heart movement.

Caution! The EpiAccess Needle is sensitive to cardiac pulsations. Cardiac pulsations can propagate through tissue and bone and be present in the diaphragm muscle. The sensitivity to these propagated cardiac pulsations may result in a false positive reading (indication by the EpiAccess System that the needle has entered the pericardium when it has not). The clinician may observe elevated indication values while the needle passes through the diaphragm. These indications should be ignored if the clinician is sure of the current needle position. The needle position should always be verified under fluoroscopy prior to any needle advancement per standard of care. In a 25 patient study conducted at 5 EU centers, the rate of false positives was 2/25 or 8%.

9. Under fluoroscopic guidance slowly advance the needle until penetrating the pericardial sac.

Needle signal pressure is presented in the middle window. This signal can vary throughout the procedure. Figure 5 displays the graphical user interface of needle pressure; arterial pressure and heart beat segmentation. This information can be used to augment the fluoroscopic information to confirm needle tip location.

Warning! Heart beat segmentation may be impacted by irregular cardiac rhythms.


Figure 5 – User Interface

10. To confirm that the needle tip is within the pericardial space, inject a small volume of contrast media through the needle. A thin film will surround the cardiac silhouette if the needle is in the pericardial space or it will pool in the extra pericardial tissue until the needle enters the potential pericardial space.
11. Insert the “J” tip of a commercially available 0.032-inch guidewire through the needle into the pericardial space, under fluoroscopy, in the left anterior oblique view.
12. Ensure the guidewire travels freely within the pericardial space and does not exit the cardiac silhouette.
13. Maintain the left anterior oblique view and confirm the guidewire can be visualized along the left heart border.
14. Retract the needle and leave the guidewire to maintain access to the pericardial space.
15. Use only commercially available 0.032-inch diameter guidewires.
16. Dispose of the EpiAccess Needles per hospital protocol for biohazardous medical devices.

Storage

The EpiAccess Needle should be stored in a cool dry place in accordance with the parameters set forth on the Unit Box Label (L0011).

Disposal

EpiAccess Needle is a single use device and should be disposed of properly in the medical waste stream.

Limited Warranty

EpiEP, Inc. provides this limited warranty that the product has been manufactured and released for use in accordance with its specifications and tested using its established test methods. EpiEP makes no other express or implied warranties, including the implied warranties of merchantability and fitness for a particular purpose. EpiEP specifically does not make any representation or warranty about the use or performance of this device used in conjunction with any other intervention devices. Handling, storage, cleaning, and sterilization of this device, as well as other factors relating to the patient, diagnosis, treatment, surgical procedure, and other matters beyond EpiEP's control may directly affect the device and the results obtained from its use. EpiEP makes no warranty, expressed or implied, including but not limited to the warranties of merchantability or fitness for a particular use with respect to devices reused, reprocessed, or re-sterilized.

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EPIEP, INC
195 Church Street
17th Floor
New Haven, CT 06510
USA
Telefono: +1 (434) 951-0444
Fax: +1(434) 951-0478

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