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Originalversion

OPERATING MANUAL

MODE D'EMPLOI

MANUAL DE INSTRUCCIONES

CT 72/P

DURCHSICHT - THERMOSTAT
TRANSPARENT THERMOSTAT
THERMOSTAT TRANSPARENTS
TERMOSTATOS TRANSPARENTES

SI Analytics

a xylem brand

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Important notes: Before starting the Transparent Thermostat CT 72/P please read the operating instructions carefully. For safety reasons the Transparent Thermostat CT 72/P may only be used for the purposes described in this present operating instructions. Please also observe the operating instructions for the units to be connected. All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics GmbH may perform additions to the Transparent Thermostat CT 72/P without changing the described properties.

Mode d'emploi	99	144
---------------------	----------	-----

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Operating Instructions Transparent Thermostat CT 72/P

Table of contents

QUALITY ASSURANCE.....	53
Guarantee	53
YOUR CONTACTS AT SI ANALYTICS GMBH.....	53
UNPACKING AND CHECKING.....	53
INTENDED USE	54
Description	54
OPERATOR RESPONSIBILITY – SAFETY RECOMMENDATIONS	54
OPERATOR RESPONSIBILITY.....	55
SAFETY NOTES.....	56
Explanation of safety notes.....	56
Explanation of other notes	56
Safety recommendations.....	57
1. TECHNICAL DATA OF TRANSPARENT THERMOSTATS CT 72/P.....	60
1.1. Technical data of immersion thermostat CT 72	60
1.2. Technical data of transparent baths	61
2. OPERATING CONTROLS AND FUNCTIONAL ELEMENTS.....	62
2.1. Overall view of transparent thermostats CT 72/P	64
3. PREPARATIONS	65
3.1. Bath liquids.....	65
3.2. Filling / draining	65
3.3. Countercooling	65
4. OPERATING PROCEDURES	66
4.1. Power connection	66
4.2. Switching on / Start - Stop	66
5. T SETTING OF TEMPERATURES.....	67
5.1. Using the pre-settings in the T menu.....	67
5.2. Direct setting of temperatures	68
6.  SAFETY INSTALLATIONS, WARNING FUNCTIONS	68
6.1. Excess temperature protection	69
6.2. Low level protection.....	69
6.3. Switch-over from warning to shutdown function.....	70
6.4. Over and Sub temperature warning function	71

7. MENU MENU FUNCTIONS	72
7.1. MENU PROGRAM – START	73
7.2. MENU PROGRAM – creation, administration	76
7.3. MENU CONFIG - Configuration of unit	78
7.3.1. REMOTE	79
7.3.2. AUTOSTART	79
7.3.3. Setting of clock and date	80
7.3.4. RESET - Factory settings	80
7.4. MENU CONTROL – Control characteristics and parameters	81
7.4.1. Dynamik internal	81
7.4.2. Control parameters – XP-, TN-, TV- INTERNAL	82
7.5. MENU SERIAL - BAUDRATE, HANDSHAKE, PARITY	83
7.6. MENU ATC - Absolut Temperature Calibration	84
7.6.1. ATC STATUS - YES / NO	86
7.6.2. CALIBRATION TYPE: 1 -/ 2 -/ 3 POINT	86
7.6.3. Example: 3-point calibration for internal control	87
7.7. MENU LIMITS	89
7.7.1. Limits for internal control	89
8. TROUBLESHOOTING GUIDE / ERROR MESSAGES	90
9. ELECTRICAL CONNECTIONS	92
10. REMOTE CONTROL	93
10.1. Setup for remote control	93
10.2. Communication with a PC or a superordinated data system	93
10.3. List of commands	94
10.4. Status messages	96
10.5. Error messages	96
11. CLEANING / REPAIRING THE UNIT	97

Quality assurance

Dear customer,

SI Analytics GmbH implements a **Quality Management System** certified according to EN ISO 9001.

This guarantees the presence of organizational structures which are necessary to ensure that our products are developed, manufactured and managed according to our customers expectations. Internal and external audits are carried out on a regular basis to ensure that our QMS system is fully functional.

Please inform us if, despite our precautionary measures, you should find any product defects. You can thus help us to avoid such faults in future.

Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss.

Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.


Your contacts at SI Analytics GmbH

Please get in contact with us or the authorized agent who supplied you with the unit if you have any further questions.

SI Analytics GmbH

Hattenbergstrasse 10
D-55122 Mainz

Telefon +49 61 31 / 66 51 11
Telefax +49 61 31 / 66 50 01
E-Mail: si-analytics@xyleminc.com
www.si-analytics.com

The following specifications should be given when product enquiries are made:

- **Unit name** printed on the front of the unit and as specified on the name plate.

Unpacking and checking

Unpack the circulator and accessories and check for damages incurred during transit. These should be reported to the responsible carrier, railway, or postal authority, and a request for a damage report should be made. These instructions must be followed fully for us to guarantee our full support of your claim for protecting against loss from concealed damage. The form required for filing such a claim will be provided by the carrier.

Intended use

The transparent thermostats CT 72/P are intended to temper water in a bath tank. Precise information about the bath media and their operating temperature can be found in Section 3.1 ("bath liquids") in this operating manual.

The transparent thermostats are not suitable for direct temperature control of foods, semi-luxury foods and tobacco, or pharmaceutical and medical products. Direct temperature control means unprotected contact of the object with the bath medium (bath fluid).

Description

The main functional elements are the heater, circulation pump, and control electronics.

The transparent thermostats CT72 are operated via the splash-proof keypad. The microprocessor technology allows different values to be set, stored, and displayed on the VFD -DISPLAY.

Temperature- and time-dependent processes can be stored and executed using the integrated programmer.

The adjustable PID cascade temperature control automatically adjusts the heat supply to the thermal requirements of the bath.

Absolute Temperature Calibration (ATC3) provides high temperature stability at all points in the bath. With the 3-point calibration, an offset is adjusted at three temperatures to ensure an accurate temperature pattern at the selected spot in the bath over the entire temperature range.

Electrical connections:

RS232 interface for modern process technology without an additional interface.

Alarm output for external alarm message or control of refrigerating baths.

The excess temperature protection according to IEC 61010-2-010 is a safety device independent from the control circuit. Its safety value can be displayed and adjusted on the VFD -DISPLAY

Operator responsibility – Safety recommendations

The products of **SI Analytics GmbH** ensure safe operation when installed, operated, and maintained according to common safety regulations. This section explains the potential dangers that may arise when operating the circulator and also specifies the most important safety precautions to preclude these dangers as far as possible.

Operator responsibility

Make sure all persons tasked with operating, installing, and maintaining the unit have read and understand the safety information and operating instructions.

Set up, maintenance, repair

Only qualified personnel is authorized to perform configuration, installation, maintenance and repairs of the thermostat.

Appropriate operation:

Routine operation can also be carried out by untrained personnel who should however be instructed by trained personnel.

The personnel operating the units should be regularly instructed about the dangers involved with their job activities as well as measures to avert these dangers.

When using hazardous materials or materials that could become hazardous, the circulator may be operated only by persons who are absolutely familiar with these materials and the circulator. These persons must be fully aware of possible risks.

Handling

Make sure the product is checked for proper condition regularly (depending on the conditions of use).

Regularly check (at least every 2 years) the proper condition of the mandatory, warning, prohibition and safety labels.

Make sure that the mains power supply has low impedance to avoid any negative effects on instruments being operated on the same mains.

Magnetic radiation may affect other devices with components sensitive to magnetic fields (e.g., monitors). We recommend maintaining a minimum distance of 1 m.

Do not expose the unit to sunlight.

Safety notes

These notes are intended to enhance your own safety consciousness.

We have set the highest quality standards for ourselves and this unit during development and production. Every unit meets relevant safety regulations.

The correct unit usage and proper handling is however solely your responsibility.

Explanation of safety notes

In addition to the safety warnings listed, warnings are posted throughout the operating manual. These warnings are designated by an exclamation mark inside an equilateral triangle. "Warning of a dangerous situation (Attention! Please follow the documentation)." The danger is classified using a signal word. Read and follow these important instructions for averting dangers.

Warning

Describes a **possibly** highly dangerous situation. If these instructions are not followed, serious injury and danger to life could result.

Caution

Describes a **possibly** dangerous situation. If this is not avoided, slight or minor injuries could result. A warning of possible property damage may also be contained in the text.

Notice

Describes a **possibly** harmful situation. If this is not avoided, the product or anything in its surroundings can be damaged.

Explanation of other notes

Note!

Draws attention to something special.

Important!

Indicates usage tips and other useful information.

This icon is used in the operating instructions to indicate flashing values or parameters which have to be set or confirmed.

Safety recommendations

Caution:

Handling chemicals, compounds and other

During the use of thermostats, viscosimeters are tempered directly inside the thermostat, e.g. submerged. You can also temper test tubes, Erlenmeyer flasks and similar items.

We do not know which substances are contained within these vessels.
Many substances are:

- inflammable, easily ignited or explosive
- hazardous to health
- environmentally unsafe

The user alone is responsible for the handling of these substances!

Our advice:

- If in doubt, consult a safety specialist.
- Read the product manufacturer's or supplier's "EC – SAFETY DATA SHEET"
- Read relevant regulations concerning dangerous materials.
- Observe relevant guidelines for laboratories in your country.

The following measures were taken for the protection of the operator:

- Protection Class I according to VDE 0106 T1 (IEC536) i.e. protection against electric shocks by grounding all parts which carry the risk of electric contact.

Notice:

Check the safety installations at least twice a year!

- Excess temperature protection according to IEC 61010-2-010.
With a screwdriver turn back the adjustable excess temperature protection until the shut-down point (actual temperature).
- Low level protection according to IEC 61010-2-010.
To check the function of the float, it can be manually lowered with a screwdriver for example.

The following measures were taken for the protection of the operator:

- Protection IP 20 according to EN 60529, i.e. regarding the protection against accidentally touching live parts and damage by foreign matter, it has been ensured that foreign bodies with a thickness or diameter of more than 12 mm cannot penetrate.

Notice:

No special precautions were taken against the penetration of water and dust.

- The device should therefore not be used in a dusty atmosphere or in the neighborhood of spray water.
- Do not insert wires or tools in any of the openings.

Notice:

Burn hazard

- Some parts of the device may become extremely hot during continuous operation.
- Always protect yourself against accidental contact.

Before initial operation

These operating instructions must be carefully studied! It contains important information on the connection to the local mains supply, correct unit usage and safe handling.

Check for transportation damage during unpacking.

Do not try to start up a damaged unit before the damage has been repaired or you have ascertained the effect of the damage. Get in contact with supplier and/or carrier for settlement of damage claims

Repairs, alterations or modifications must only be carried out by specialist personnel. Considerable damage can be caused by improper repairs. The SI Analytics GmbH service department is at your disposal for repair work.

Never remove warning labels!

Initial operation

Only connect the unit to a power socket with an earthing contact (PE – protective earth). Electrical connections and any other work must be performed by qualified personnel only. The power supply plug serves as a safe disconnecting device from the line and must always be easily accessible. Never operate units with damaged mains power cables.

Place the unit on an even surface on a base made of nonflammable material.

Set the excess temperature safety installation at least 25 °C below the fire point of the bath fluid.

Operation

Only use this unit solely for the intended application.

Never operate the unit without bath fluid in the bath! Only use the heat transfer liquids recommended by SI Analytics GmbH

Do not operate the unit with wet or oily hands.

Do not expose the unit to spray water or immerse it in water.

Do not drain the bath fluid while it is hot.

Cleaning

Do not clean the unit with solvents (fire risk!) – a wet cloth soaked in household detergent is normally sufficient.

Transportation

Do not move the instrument while it is running or hot. Always empty the bath before moving the unit. Transport the unit with care. Sudden jolts or drops may cause damage in the interior of the unit.

Notice:

Shock hazard

Complete separation from mains is required when:

- all dangers caused by this device are to be avoided,
- cleaning is carried out;
- repairs or maintenance work is about to be carried out

Complete separation means:

Pull out the mains plug!

Disposal

Valid in EU countries

See the current official journal of the European Union – WEEE directive.

Directive of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE).

This directive requires electrical and electronic equipment marked with a crossed-out trash can to be disposed of separately in an environmentally friendly manner.

Contact an authorized waste management company in your country.

Disposal with household waste (unsorted waste) or similar collections of municipal waste is not permitted!

1. Technical data of Transparent Thermostats CT 72/P

Translation of the legally binding German version

(Version as of 10. May 2010)

Country of origin: Germany / Made in Germany

1.1. Technical data of immersion thermostat CT 72

Temperature ranges:

Operating temperature range			
with additional cooling	°C		+ 10 ... + 60
Working temperature range			
(20 °C ambient temperature)	°C		+ 30 ... + 60
with tap water cooling*	°C		+ 20 ... + 60
with flow-through cooler CK300/CK310 (option)	°C		+ 10 ... + 60
Temperature stability	K		± 0,02
Temperature drift at changing ambient temperature			
between 5 and 40 °C:			max. 0,002 K/1°C
Temperature control			PID3 cascade temperature control
<u>Heater capacity:</u>			
at 230 V	W		2000
at 115 V	W		1000

<u>Pump</u> :	Pressure max. at 0 liter	bar	0.35
	Flow rate at 0 bar	l/min	15

Connections:

electrical:			
connector: control cable to the flow-through cooler CK300/CK310	Vdc/mA		24-0 / max. 25
Computer interface			RS232
mechanical			
nozzles for flow-through cooler. CK300 / CK310			2 x Ø 12 mm
nozzles for tap water cooling			2 x Ø 8 mm
Overall dimensions	(WxDxH)	cm	13x16x33
Weight		kg	5.0
Mains power connection	230 V/50-60 Hz	V/ Hz	190 ... 253 / 50-60
Current draw	(at 230 V)	A	9
Mains power connection	115 V/60 Hz	V/ Hz	90 ... 127 / 60
Current draw	(at 115 V)	A	9

* If tap water is used, an operating temperature of approx. 3 °C above the given tap water temperature can be maintained.

Safety installations according to IEC 61010-2-010:

Excess temperature protection	adjustable from 0 °C ... 230 °C
Low liquid level protection	float switch
Classification according to DIN 12876-1	class III

Supplementary safety installations

Early warning system for low level	float switch
High temperature warning function	optical + audible (in intervals)
Low temperature warning function	optical + audible (in intervals)
Supervision of working sensor	plausibility control
Reciprocal sensor monitoring between working and safety sensors	difference >35 °C
Alarm message	optical + audible (permanent)
Warning message	optical + audible (in intervals)

Environmental conditions:

- Use only indoor.
- Altitude up to 2000 m - normal zero.
- Ambient temperature: +5 ... +40 °C
- Air humidity:
- Max. rel. humidity 80 % for temperatures up to +31 °C,
- linear decrease down to 50 % relative humidity at a temperature of +40 °C
- Max. mains fluctuations of ±10 % are permissible.

Protection class according to IEC 60 529	IP20
The unit corresponds to Class I	
Overtoltage category	II
Pollution degree	2

Caution:

The unit is not for use in explosive environment

Standards for interference resistance according to EN 61326-1

This unit is an ISM device classified in Group 1 (using high frequency for internal purposes) Class A (industrial and commercial range).



1.2. Technical data of transparent baths

<u>Filling volume:</u>	L	18
<u>Overall dimensions</u>	(WxDxH)	cm 38 x 24 x 32
<u>Weight:</u>	kg	Approx. 4.6 (without bath liquid)
<u>Bath liquid:</u>		Water
<u>Temperature range:</u>	°C	+10 ... +60
<u>Material:</u>		Acryl, high-grade steel
<u>Country of origin:</u>		Germany / Made in Germany

2. Operating controls and functional elements

Front view CT 72

Rear view CT 72


1

Mains power switch, illuminated

2

VFD-DISPLAY

Header: Control indicators see sections 2.1 and 2.2

- Line 1: Actual value internal
- Line 2: Working temp. setpoint, constantly S xxx.xx
- Line 3: Actual value (I = internal)
Alternating with the display in line 1

2.1

Control indicators in the header:

Heating / Cooling / Alarm /
Remote control

2.2

Control indicators in the header:

Temperature indication **Internal** actual value (**Ext** not possible on this unit)
Temperature indication in °C (°F not possible on this unit)

3

Adjustable excess temperature protection according to IEC 61010-2-010

4

Navigation keys

- 4.1  1. Key: > ENTER < Start / Stop (pump / heater)
 2. > ENTER < in the menu Menu item / select submenu for setting
 Save set value
 Save selected parameter

A beep signals the end of setting

After the actions Start, Stop and change from VFD Display to standard display the key  is locked for a short time. The above graph "front side" shows an example for standard display.

4.2

ESC

1. Key: >Return< Stop (pump / heater)
 2. >Return< in the menu one menu level down
 Correction function for parameters or values (prior to ENTER)

immediately back to standard display

 -  icon for „keep key pressed down“.

4.3

1. Key: >Up / Down <temperature – increase/decrease setpoint
 Push key quickly for single steps,
 Keep key pressed for fast change.
 2. >Up/Down< in the menu selection of menu items / parameters

Menu keys

4.4

Key: start the menu > warning and safety values<

4.5

T

Key: start the menu >temperature setpoints<

4.6

MENU

Key: display of MENU structure

5

RS232

Interface RS232: remote control via personal computer

6

5-pin connecting socket for flow-through cooler CK300/CK310 via adapter cable

7

Mains power cable with plug

8


Mains circuit breakers (resettable) 15 A


2.1. Overall view of transparent thermostats CT 72/P

Basic unit und transparent bath are supplied ready for operation.

Front view CT 72/P


Rear view CT 72/P

10a / 10b Cooling coil nozzles for tap water cooling

11 / 12 Nozzles for flow-through cooler CK300 / CK310

11 – IN - Return

12 – OUT – Feed

3. Preparations

3.1. Bath liquids

Caution

Water:

The transparent thermostat CT 52 is designed for operation with water **only!** SI Analytics GmbH recommends the mixing ratio of 70 % deionized water and 30 % tap water. The quality of water depends on local conditions.

- Due to the high concentration of lime, hard water is not suitable for temperature control because it leads to calcification in the bath.
- Ferrous water can cause corrosion - even on stainless steel.
- Chloric water can cause pitting corrosion.
- Distilled and deionized water is unsuitable. Their special properties cause corrosion in the bath, even in stainless steel.

3.2. Filling / draining

Notice

- Prevent humidity from entering into the immersion thermostat CT72.
- Do not drain the bath fluid while it is hot!
Check the temperature of the bath fluid prior to draining (by switching the unit on for a short moment, for example).
- Store and dispose the used bath fluid according to the laws for environmental protection.
- Acryl- baths are not resistant to solvents.
Avoid contamination of the bath fluid with solvents.

Filling

- Recommended maximum filling level with water as bath liquid is 33 mm below the cover plate.

Draining

- Press the mains switch to turn the transparent thermostat CT 52 off.
- Disconnect the power plug from the mains socket!
- Empty the bath tank (with water dispense pump, for example).

3.3. Countercooling

Notice

- Observe the laws and regulations of the water distribution company valid in the location where the unit is operated.
- Make sure that the tubes are securely attached.
- Employ suitable connecting tubing (temperature range).
- Avoid mechanical work (torsion, bending) to be applied to tubing.
- Regularly check the tubing for material defects (e.g. for cracks).

For applications near the ambient temperature, the cooling coil must be connected to the water mains.

Using tubing, connect the cooling coil (10a / 10b) to the tap water supply, and lead the tap water in a sink through the return connector.

4. Operating procedures

4.1. Power connection

Caution

- Only connect the unit to a power socket with earthing contact (PE – protective earth)!
- The power supply plug serves as safe disconnecting device from the line and must be always easily accessible.
- Never operate equipment with damaged mains power cables.
- Regularly check the mains power cables for material defects (e.g. for cracks).
- We disclaim all liability for damage caused by incorrect line voltages!

Check to make sure that the line voltage matches the supply voltage specified on the identification plate. Deviations of $\pm 10\%$ are permissible.

4.2. Switching on / Start - Stop

Switching on:

- Turn on the mains power switch (1).

① The unit performs a self-test. Then the software version appears. The display „OFF“ or „R OFF“ indicates the unit is ready to operate.

① The circulator enters the operating mode activated before switching the circulator off:

keypad control mode (manual operation)

or

remote control mode (operation via personal computer).

Start:

- Press  key.


The actual bath temperature is displayed on the VFD COMFORT-DISPLAY. The circulating pump starts with a slight delay.

Stop:

Press  key.

or

Keep  key pressed.

The unit also enters the safe operating state "OFF" or "r OFF" after a mains power interruption. The temperature values entered via the keypad remain in memory. With the circulator in keypad control mode, press the start/stop key to restart operation. With the circulator in remote control mode, the personal computer must first resend the parameters set via the interface before the circulator may be restarted.

5. T Setting of temperatures

Notice

Observe the limited working temperature range when using plastic bath tanks.
Working temperature range see technical data page 60.

5.1. Using the pre-settings in the T menu

Press the **T** key to call up the menu for temperature selection.

3 different working temperatures can be adjusted. Their values are freely selectable within the operating temperature range.

- ① The temperatures can be set in start or stop mode.
- ① Press **ESC** key if a value is to be retained

Setting of working temperature in the T menu



1. Press the key **T**. The value flashes .
2. Select SETPOINT 1 or 2 or 3 using the key  or .
3. Confirm by pressing the  key.

- ① The circulator uses the new working temperature value for temperature control.

Factory setting:

SETPNT 1	25 °C
SETPNT 2	37 °C
SETPNT 3	70 °C

Example: Adjustment/modification of the pre-setting of "SETPOINT 3"

1. Press the **T** key.
2. Select SETPOINT 3 by pressing the  key.
Example: SETPNT 3 / 70.00 °C
3. Keep the  key pressed until the integer digits flash .
- (example: <70>)
4. Adjust value by pressing the  key and the  key to 85.00 °C and confirm by pressing the  key.
The decimal digits flash  and can be adjusted if desired.
Confirm once more by pressing the  key.
Example on the left: SETPNT 3 / 85.00.

- ① If the active setpoint (SETPNT) is changed, the new value is immediately used for the control of the working temperature.
The heater control indicator flashes.
- ① If the other two setpoints (not activated for control) are changed the **T** MENU has to be left by pressing the **ESC** key after the decimal digits have been confirmed

Notice: Refer to SETPOINT MAX / MIN in chapter 7.7. MENU LIMITS

5.2. Direct setting of temperatures

The circulator uses the setpoint of SETPNT 1 or 2 or 3 for temperature control

The indicated setpoint temperature can be changed directly any time.

Example: change 25.00 °C to 45.00 °C

1. By pressing the key the circulator switches to the active SETPOINT< example on the left: >SETPNT / 1 25.00°C<.

The integer digits flash (example: <25>).

2. By pressing the keys and the value is changed to 45.00 °C and is confirmed by pressing the key.

The decimal digits flash and can be adjusted if desired.

Confirm once more by pressing the key.

① The circulator uses the new working temperature value for temperature control.

① The temperatures can be set in start or stop mode.

6. Safety installations, warning functions

Check the safety installations at least twice a year! Refer to (page **Fehler!** **Textmarke nicht definiert.**)

SECVAL
(Security Values)

- SAFETMP
- AL-TYPE
- OVERTMP
- SUBTEMP

Settings for the excess temperature protection > **SAFETMP**<

and for the warning functions for high > **OVERTMP**< and low > **SUBTEMP**< temperature are made in a menu which is called up by pressing the key .

Menu item > **AL-TYPE**< allows choosing between a warning and

an alarm cut-off for the menu items > **OVERTMP**< and > **SUBTEMP**<.

6.1. Excess temperature protection

Warning

The excess temperature protection **must** be set at least 25 °C below the fire point of the bath fluid used!

In case of wrong setting there is a fire hazard!

We disclaim all liability for damage caused by wrong settings!

-OFF-
ALARM
CODE 14

This excess temperature protection is independent of the control circuit. When activated heater and circulating pump are completely shut down.

The alarm is indicated by optical and audible signals (continuous tone) and the error message "ALARM-CODE 14" appears on the VFD COMFORT-DISPLAY together with the ticker:

> EXCESS TEMPERATURE PROTECTOR ALARM-CHECK ADJUSTMENT <

Setting range: 20 °C ... 230 °C

① Rough setting can be effected by using the temperature scale.

Exact setting:


1. Press the key  to display menu >SAFETMP<.

2. Press the  key and the set shutdown value is indicated.

Set the new shutdown value within 30 seconds using a screwdriver. The value is indicated on the VFD COMFORT-DISPLAY

Example: SAFETMP / 50 °C

SECVAL
SAFETMP

Recommendation:

Set the excess temperature protection at 5 °C to 10 °C above the working temperature setpoint.

6.2. Low level protection

This low level protection is independent of the control circuit.

If the low liquid level protection device is triggered, a complete shutdown of the heater and circulating pump is effected.

A continuous alarm sounds and together with the ticker: > LOW LEVEL ALARM-FILL MEDIUM < a message appears on the VFD -DISPLAY: "ALARM-CODE 1"

-OFF-
ALARM
CODE 1

Turn off the unit with the mains switch, refill bath fluid and turn the unit on again!

Notice

When adding bath fluid, always use the type of fluid which is identical with the fluid in the bath.

6.3. Switch-over from warning to shutdown function

```
SECVAL
AL-TYPE
```

If a shutdown of functional elements (e.g. heater, circulating pump) is required when the limit values are exceeded or undercut the circulator can be changed over from warning function >WARNING< to shutdown function >ALARM<.

Factory setting:
>WARNING<

```
AL-TYPE
WARNING ✱
```

```
AL-TYPE
ALARM
```

1. Press the key .
2. Select the menu >SECVAL -AL-TYPE< by pressing the  key.
3. Press the  key and the set parameter will flash ✱.
(Example: WARNING)
4. Change the parameter by pressing the  key and confirm by pressing the  key.
or
press the **ESC** key if the parameter is to be retained.

i Setting >WARNING<

A mere warning function with optical and audible warning signal (interval tone)
A message appears on the VFD COMFORT-DISPLAY:

```

  88888  88888
  WARNING WARNING
  CODE 03 CODE 04
  OVERTMP SUBTEMP
```

• Setting >ALARM<

Temperature limit with shutdown of heater and circulating pump.
An audible alarm sounds (continuous tone) and a message appears on the VFD COMFORT-DISPLAY:

```

  -OFF-  -OFF-
  ALARM  ALARM
  CODE 03 CODE 04
  OVERTMP SUBTEMP
```

6.4. Over and Sub temperature warning function


Over temperature

OVERTMP
200.00

Sub temperature

SUBTMP
-99.00

If the observance of a working temperature value >SETP< has to be supervised for a sensitive temperature application, then set over and sub temperature warning values. In the example below the SETPOINT 45 °C is surrounded by the values OVERTMP 487 °C and SUBTMP 43 °C. The electronics immediately register if the actual temperature breaches one of the set limit values. The resulting reaction is defined in a further menu item. (See chapter 6.3. Switch-over from warning to shutdown function)

1. Press the key .
2. By pressing the  or  key select the menu >OVERTMP< or >SUBTMP<.
3. Press the  key. The integer digits flash
4. Change the values to 87. °C and/or 83. °C by pressing the  and  key and confirm with the  key.
The decimal digits flash and can be adjusted if desired.
Confirm once more by pressing the  key.
See above examples.

 The warning functions are only activated if the actual bath temperature remains within the set limit values for 3 seconds after switch-on.

Recommendation:

Set the over temperature warning value >OVERTMP< 5 °C to 10 °C above the working temperature setpoint.

Set the sub temperature warning value >SUBTMP< 5 °C to 10 °C below the working temperature setpoint.

7. MENU Menu functions

The term „Menu functions“ refers to settings such as

MENU ↓

Menu level 1

1. Open the menu by pressing the **MENU** key.
2. Use the   keys to scroll in menu level 1.
3. Press the  key to change to menu level 2.
Press the **ESC** key if settings are to be retained.

MENU
P-START 
 

Start program

Page 73

P-START
STEP

MENU
PROGRAM

Administration and creation of programs

Page 76

MENU
CONFIG

Configuration of the unit
REMOTE – on / off (remote control via RS232)
AUTOST – AUTOSTART on / off
TIME / DATE – setting time and date
RESET – factory settings

Page 78

MENU
CONTROL

Control characteristics and parameters
C-TYPE – Internal or external control

DYNAMICS - internal
Control parameter - XP-, TN-, TV- INTERNAL

Page 81

MENU
SERIAL

Adjustable interface parameters
BAUDRATE, PARITY, HANDSHAKE

Page 83

MENU
ATC

ATC - Absolute Temperature Calibration,
Sensor calibration INTERNAL SENSOR,
3-point calibration

Page 84

MENU
LIMITS

Limitations of temperature and capacity
SETPOINT MAX / MIN - Maximum and minimum setpoint
HEAT MAX – Set maximum heating
COOLING MAX – Set maximum cooling
INTERN MAX / MIN – Limitation of the temperature range
BAND HIGH / LOW – Band limit

Page 89

7.1. MENU PROGRAM – START

Start-Menu

This menu will start a previously set program.

i Requirements:

1. Create a program. (refer to next chapter)
2. Return to the Start-MENU and confirm the desired setting of each MENU item with the key
3. Set a start time (>TIME< >DATE< >YEAR<) if the program is to be started by the internal timer.

Menu level 1

MENU
P-START

- > STEP< Program start at section 1 ... 10
- > RUNS < Number of repetitions 1 ... 99
- > END< Status at end of program (STDBY/SETPNT)
Standby or last setpoint
- > GO < Time of start (NOW/TIMER)

Level 2

Parameter level

i Press the **ESC** key if a parameter is to be retained.
Correction function for parameters or values (prior to ENTER)

P-START
STEP

STEP

- Set program step with and
- example: STEP 1

P-START
RUNS

RUNS

- Set number of runs with and
- example: 1 run

P-START
END

END
STDBY

or

END
SETPNT

- Set desired parameters with and .
- i** (STDBY / SETPNT)
Parameter **STanDBY**: the circulator switches to -OFF-.
Parameter **SETPoiNT**: the circulator constantly keeps the temperature at the value of the last step.

i

P-START
GO

GO
NOW



or

- Confirm >NOW< with the key and the program will start immediately

- i** or
start at the set time under parameter (TIMER).
Set time in the example below:
09. August 2009, 11:15 hrs

GO
TIMER

Menu functions

 <p>60 TIMER</p> <p>Submenu TIMER</p>	<p>Parameter level</p>	<p>① set the time for the start of the program in the submenu >TIMER<.</p>
 <p>TIMER TIME</p>	<p>TIME 11:15</p>	<p>>TIME< hours/minutes (hh:mm), set both values one after the other and confirm</p> <ul style="list-style-type: none"> hours flash, set by pressing   +  minutes flash, set by pressing   + 
 <p>TIMER DATE</p>	<p>DATE 09:08</p>	<p>>DATE< day/months (TT/MM), set both values one after the other and confirm.</p> <ul style="list-style-type: none"> day flashes, set by pressing   +  month flashes, set by pressing   + 
 <p>TIMER YEAR</p>	<p>YEAR 2009</p>	<p>>YEAR< year</p> <ul style="list-style-type: none"> Set the year with   and .
 <p>TIMER START</p>	<p>START OK</p>	<p>① The program starts at the set time.</p>
 <p>-OFF- 5 25:00 TIMER</p>		<p>① Display of time until start: In line 3 the notice >TIMER< and the set values for „TIME“ and „DATE/YEAR“ are alternately indicated</p> <p>① Check the correct setting of the internal real time clock if required (see MENU CONFIG)</p>

The started program

After the start the program will indicate the currently calculated setpoint in line 2 S XX.XX. The value increases within the time period >TSLICE< until the target temperature >SETPNT< of the section is reached.

If the time period in a section is set to „0“, the next section will not begin until the target temperature has been reached.

```

42.36
S  XX.XX
STEP XX
    
```

A

```
S  XX.XX
XXhXX.XX
```

B

```
S  XX.XX
XXhXX.XX
```

C1

```
S  XX.XX
I  XX.XX
```

C2

```
S  XX.XX
E  XX.XX
```

D1

```
S  XX.XX
RUN
```

D2

```
S  XX.XX
PAUSE
```

Use the edit keys   to scroll to line 3. The display changes approximately every 4 seconds between the current section (STEP XX) and the

- A remaining time of the section
- B remaining time of the program
- C current bath temperature
I xxx.xx - internal actual value or
E xxx.xx – external actual value
- D RUN – the program has started or
PAUSE – the progress of the program has been interrupted by pressing the **ESC** key. While the time is stopped the temperature will constantly remain at the last calculated setpoint
Continue with the **ESC** key.

Termination / Interruption of a program

```

-OFF-
S  XX.XX
I  XX.XX
    
```

- ① The program can be terminated any time by pressing the key **ESC** .
- ① In case of power failure the program is interrupted.
The circulator switches to –OFF- .
- ① If the AUTOSTART-function is activated the programmer starts again at the point in time approx. 5 minutes prior to the interruption. However, an uncontrolled change of the bath temperature has occurred.

7.2. MENU PROGRAM – creation, administration

Menu level 1


MENU
PROGRAM

1 program

10 sections

The integrated programmer permits fast and easy programming of setpoint temperature sequences. This temperature sequence is called program. A program is composed of individual sections (STEP). The sections are defined by duration (TSLICE) and target temperature. The target temperature is the setpoint (SETPNT), which is achieved at the end of a section. The programmer calculates the temperature ramp from the difference in time and temperature.

STEP	(Nr.)	1	2	4	6	8
SETPNT	(°C)	100	180	180	75	75
TIME	(hh.mm)	00:20	00.10	00:20	00:50	00:20

Graph1

① Sections without set value and time are skipped. They can be defined retroactively and the integrated into the program.

Menu level 1

▼ ▲
MENU
PROGRAM
← →

>EDIT< Create, administer program
>STEP< Program step (1 ... 10)
>SETPNT < Temperature setpoint of step ...
>TSLICE< Duration of step ...
>DELETE< delete program step (01 ... 10, ALL)

① Press **ESC** key, if a parameter is to be retained.
Correction function for parameters or values (prior to ENTER)


Level 2	Level 3	Parameter level
PROGRAM EDIT ▼	EDIT STEP ← → EDIT 01 STEP ← ▼ ▲	STEP ⏏ (STEP 1 ... 10) • Set program step with ▼ ▲ and ← → ← (Example: EDIT STEP 01) ① For STEP 01 the values for SETPOINT 01 and TSLICE 01 are set one after the other

Menu functions

Level 2	Level 3	Parameter level
	<pre>EDIT 01 SETPNT</pre>    	<pre>SETPNT 100.00</pre>  (values within working temp. range) <ul style="list-style-type: none"> Integer digits flash, set by pressing    Decimal digits flash, set by pressing   
	<pre>EDIT 01 TSLICE</pre>   (time slice)	<pre>TSLICE 00.10</pre>  <ul style="list-style-type: none"> Set duration by pressing   and 
 <pre>PROGRAM DELETE</pre>  	<pre>DELETE STEP --</pre> 	<ul style="list-style-type: none"> Delete program Program steps can be deleted individually or entirely. (STEP 01, 02,... 10, ALL). Set parameters by pressing   and 

7.3. MENU CONFIG - Configuration of unit

<p>Menu level 1</p> <p>▼ ▲</p> <p>MENU CONFIG</p> <p>← ↵ ↓</p>	<p></p> <p>ⓘ A RESET can be effected only in the >OFF< mode. Switch off circulator by pressing the ← ↵ key and call up the menu CONFIGURATION.</p>	
<p>Level 2</p>	<p>Parameter level</p>	<p>ⓘ Press the key ESC if a parameter is to be retained. Correction function for parameters and values (prior to ENTER).</p>
<p>CONFIG REMOTE</p> <p>▼ ▲</p> <p>← ↵ →</p>	<p>REMOTE OFF ✖</p> <p>or</p> <p>REMOTE ON</p>	<ul style="list-style-type: none"> • Switch on and off remote control by pressing ▼ ▲ and ← ↵ ⓘ Control display in the topline R for Remote ⓘ For remote control refer to page 93 ⓘ ○ ○ Connect RS232 with PC.
<p>CONFIG AUTOST</p> <p>▼ ▲</p> <p>← ↵ →</p>	<p>AUTOST OFF ✖</p> <p>or</p> <p>AUTOST ON </p>	<ul style="list-style-type: none"> • Switch on and off Autostart by pressing ▼ ▲ and ← ↵ AUTOSTART on = on AUTOSTART off = off See WARNING page 79
<p>Level 3</p>	<p>Parameter level</p>	<p>Parameter level</p>
<p>CONFIG TIME/DT</p> <p>▼ ▲</p> <p>← ↵ →</p>	<p>TIME/DT TIME</p> <p>← ↵ →</p>	<p>TIME 11.15 ✖</p> <ul style="list-style-type: none"> • Hours flash, set by pressing ▼ ▲ + ← ↵ • Minutes flash, set by pressing ▼ ▲ + ← ↵
<p>TIME/DT DATE</p> <p>← ↵ →</p>	<p>DATE</p> <p>← ↵ →</p>	<p>DATE 08.09 ✖</p> <ul style="list-style-type: none"> • Day flashes, set by pressing ▼ ▲ + ← ↵ • Month flashes, set by pressing ▼ ▲ + ← ↵
<p>TIME/DT YEAR</p> <p>← ↵ →</p>	<p>YEAR</p> <p>← ↵ →</p>	<p>YEAR 09 ✖</p> <ul style="list-style-type: none"> • Year flashes, set by pressing ▼ ▲ + ← ↵

	Level 3	Parameter level
	 	<ul style="list-style-type: none"> Return to factory settings by pressing  <p>RESET returns all set values to the factory setting except for date and time.</p> <ul style="list-style-type: none"> i A RESET can be effected only in the –OFF- mode. i During the message –RUN- all parameters are reset to factory settings.


7.3.1. REMOTE

Factory setting: OFF

The control electronics offer two ways of adjusting a setpoint.

- Adjustment of setpoint using the keypad or the integrated programmer.
- Adjustment of setpoint via the serial interface RS232 using a PC or a superordinated process control system.

i The topline of the VFD-DISPLAY shows a bright „R“ for remote control; - remote control discontinued.

 RS232	<p>IMPORTANT: additional measures for remote control</p> <ul style="list-style-type: none"> i Connect the circulator to the PC using an interface cable. i Check the interface parameters of both interfaces (Thermostat and PC) and make sure they match. (refer to Setup for remote control page 93)
--	---

7.3.2. AUTOSTART

Warning

For supervised or unsupervised operation with the “AUTOSTART“ function avoid any hazardous situation to persons or property

Take care to fully observe the safety and warning functions of the circulator.

Factory settings: OFF

Notice:

The circulator has been configured and delivered by JULABO in accordance with the NAMUR recommendations. This means for the start mode that the unit must enter a safe operating status after a power failure. This safe operating status is indicated by the message „OFF“ or „R OFF“ on the VFD COMFORT-DISPLAY.

A complete, all-pole shutdown of the main functional elements such as heater and pump motor is effected.

The values set on the circulator remain saved and the unit is restarted by pressing the start/stop key in manual control.

In remote control mode the values need to be reset by the PC via the interface.

If such a safety standard is not required, the NAMUR recommendations can be bypassed with the AUTOSTART function thus allowing a direct start of the circulator by pressing the mains switch or using a timer.

7.3.3. Setting of clock and date

The internal real time clock allows starting a program any time. The clock is set to the local mean time (CET) at the factory.

- ① If the unit is operated in a different time zone, the clock can be adjusted in this menu.
- ① Change summer/winter time in this menu

7.3.4. RESET - Factory settings

A Reset will return all values to factory setting except for date and time.



- ① A RESET can be effected in the >OFF< mode only.
Switch off the circulator by pressing the key  and call up the menu CONFIGURATION.

7.4. MENU CONTROL – Control characteristics and parameters

<p>Menu level 1</p> <p>▼ ▲</p> <p>MENU CONTROL</p> <p>← →</p>	<p>The Thermostat is qualified for internal temperature control</p> <p>ⓘ Press the ESC key if a parameter is to be retained. Correction function for parameters or values (prior to ENTER)</p>	
<p>Level 2</p> <p>CONTROL DYNAMIC</p> <p>▼ ▲</p>	<p>Parameter level</p> <p>DYNAMIC APER ✖</p> <p>or</p> <p>DYNAMIC NORM</p>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ← →
<p>CONTROL xP INT</p> <p>▼ ▲</p>	<p>xP INT 16 ✖</p> <p>0.1 ... 99.9</p>	<ul style="list-style-type: none"> The parameter flashes, set by pressing ▼ ▲ + ← →
<p>CONTROL TN INT</p> <p>▼ ▲</p>	<p>TN INT 100 ✖</p> <p>3 ... 9999</p>	<ul style="list-style-type: none"> The parameter flashes, set by pressing ▼ ▲ + ← →
<p>CONTROL TV INT</p> <p>▼ ▲</p>	<p>TV INT 5 ✖</p> <p>0 ... 999</p>	<ul style="list-style-type: none"> The parameter flashes, set by pressing ▼ ▲ + ← →

7.4.1. Dynamik internal

CONTROL DYNAMIC

This parameter affects the temperature sequence only in case of internal control.

Factory setting: APER (aperiodic)

Possible parameters:

NORM Allows for reaching the setpoint faster – with setpoint change or ramp function – but overshooting of up to 5 % is possible.

APER Ramp function: the increase of temperature occurs temporally offset and achieves the target temperature without overshooting.
Setpoint change: The temperature increases at the same rate, the target temperature is achieved without overshooting.

ⓘ With both settings constant temperature is achieved after approximately the same time.

7.4.2. Control parameters – XP-, TN-, TV- INTERNAL

In most cases the control parameters preset in the factory are adequate for achieving an optimum temperature sequence. The control parameters allow adjustment to special control processes..

```
xP INT
  16
```

Setting range: 0.1 ... 99.9

Proportional range >Xp<

The proportional range is the range below the setpoint in which the control circuit reduces the heating capacity from 100 % to 0 %

```
TN INT
  100
```

Setting range: 3 ...9999

Reset time >Tn< (Integral component)

Compensation of the remaining control deviation due to proportional regulation. An insufficient reset time may cause instabilities. Excessive reset times will result in unnecessary prolongation of compensation of the control difference.

```
TV INT
  5
```

Setting range: 0 ... 999


Lead time >Tv< (Differential component)

The differential component reduces the transient time. An insufficient lead time will prolong the time required for compensation of disturbance effects and cause high overshooting during run-up. An excessive lead time could cause instabilities (oscillations)

Optimization instructions for the PID control parameters

Optimum setting


Control parameters XP-, TN-, TV- INTERNAL


The heat-up curve reveals possible faulty settings of the control parameter.

Inappropriate settings may produce the following heat-up curves:


Xp too low


Tv/Tn too low


Xp too high or Tv too high

Tv/Tn too high or Xp too high

7.5. MENU SERIAL - BAUDRATE, HANDSHAKE, PARITY

<p>Menu level 1</p> 	<p>For communication between circulator and a PC or a superordinated process control system the interface parameters of both units must be identical.</p> <p>❗ For remote control refer to page 93</p> <p>Factory settings: 4800 Baud even hardware handshake</p>	
<p>Level 2</p>	<p>Parameter level</p>	<p>❗ Press the ESC key if a parameter is to be retained.</p>
		<ul style="list-style-type: none"> The parameter flashes, switch by pressing 
		<ul style="list-style-type: none"> The parameter flashes, switch by pressing  <p>even: Data bits = 7; Stop bits = 1 odd: Data bits = 7; Stop bits = 1 no: Data bits = 8; Stop bits = 1</p>
		<ul style="list-style-type: none"> The parameter flashes, switch by pressing  <p>Xon/Xoff-protocol (Software handshake) Protocol RTS/CTS (Hardware handshake)</p>


7.6. MENU ATC - Absolut Temperature Calibration

MENU
ATC


ATC serves to compensate a temperature difference that might occur between circulator and a defined measuring point in the bath tank because of physical properties.

Example:

1-point calibration


3-point calibration

T_T 1 = Original curve

Principle:

For ATC calibration, in steady state the bath temperature at the location of the temperature sensor (CT) is determined at the respective adjusted working temperature. This value is then set on the circulator in the menu >ATCalibration< under menu item >CTEMP X<. This can be a 1-point, 2-point or 3-point calibration.

- M = Temperature measuring instrument with temperature sensor
- B = Bath tank (INTernal)
- T = circulator
- CT = Temperature on measuring point
- TT = Temperature on Thermostat

Menu-Ebene 1

▼ ▲
MENU
ATC
◀ ▶

ⓘ Press the **ESC** key if parameter is to be retained. Correction function for parameters or values (prior to ENTER).

Level 2	Parameter level	
<p>ATC (I)</p> <p>STATUS</p> <p>▼ ▲</p>	<p>STATUS</p> <p>YES ✕</p> <p>or</p> <p>STATUS</p> <p>NO</p>	<ul style="list-style-type: none"> • The parameter flashes, switch by pressing ▼ ▲ and ▶ ⓘ >NO< Carry out an ATC calibration ⓘ >YES< return to standard operation after calibration.

Menu functions

Level 2	Parameter level	
<pre> ATE <I> TYPE ▼ ▲ </pre>	<pre> TYPE I-POINT ※ TYPE 2-POINT TYPE 3-POINT </pre>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ↵ ⓘ A >1-point<, >2-point< or >3-point< calibration can be carried out. <p>The selected calibration is indicated on level 2 by 1 or 2 or 3.</p>
<pre> ATE <I> TMPVAL1 ▼ ▲ ATE <I> CALVAL1 ▼ ▲ </pre>	<pre> TMPVAL1 8000 ESC CALVAL1 79.70 ※ </pre>	<p>The value >TMPVAL< is only indicated</p> <ul style="list-style-type: none"> ⓘ In addition the measured temperature value >CALVAL X< is saved during the next step. Integer digits flash, set by pressing ▼ ▲ + ↵ Decimal digits flash, set by pressing ▼ ▲ + ↵
<p>ⓘ If only a 1-point calibration is carried out, the following menu items are not indicated anymore</p>		
<pre> ATE <I> TMPVAL2 ▼ ▲ </pre>	<pre> TMPVAL2 12000 ESC </pre>	<p>The value is only indicated</p>
<pre> ATE <I> CALVAL2 ▼ ▲ </pre>	<pre> CALVAL2 119.50 ※ </pre>	<ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ + ↵ Decimal digits flash, set by pressing ▼ ▲ + ↵
<p>ⓘ If only a 2-point calibration is carried out, the following menu items are not indicated anymore</p>		
<pre> ATE <I> TMPVAL3 ▼ ▲ </pre>	<pre> TMPVAL3 16000 ESC </pre>	<p>The value is only indicated</p>
<pre> ATE <I> CALVAL3 ▼ ▲ </pre>	<pre> CALVAL3 159.30 ※ </pre>	<ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ + ↵ Decimal digits flash, set by pressing ▼ ▲ + ↵

7.6.1. ATC STATUS - YES / NO

```
ATC (I)
STATUS
```

```
STATUS
YES
```

```
STATUS
NO
```

In the second submenu the ATC function for the temperature sensor selected above is activated >YES< or deactivated >NO<.

>YES< (factory setting) The controller of the circulator uses the original curve of the temperature sensor or the new curve measured during the ATC calibration.

Important: Set to >NO< during the calibration process

>NO< An ATC calibration is to be carried out.

Important: Set to >YES< after calibration.

ⓘ In the > ATC STATUS < >YES< the ATC calibration always affects the current working temperature; also the one set via interface.

7.6.2. CALIBRATION TYPE: 1 -/2 -/3 POINT

```
ATC (I)
TYPE
```

```
TYPE
1-POINT
```

```
TYPE
2-POINT
```

```
TYPE
3-POINT
```

```
TMPVAL1 CALVAL1
8000 79.70
```

```
TMPVAL2 CALVAL2
12000 119.50
```

```
TMPVAL3 CALVAL3
16000 159.30
```

A >1-point<, >2-point< or >3-point< calibration can be carried out.

First geometrically define the location for calibration (measuring point CT), then determine the temperature values of the calibration points.

The type of calibrations also determines the number of the following pairs of values indicated on the VFD-DISPLAY.

Pairs of values:

TMPVAL X: Circulator temperature 1 or 2 or 3 (actual value TT)

The actual temperature of the bath is simultaneously saved with the "calibration value" >CALVAL< and can be indicated for control purposes (value does not flash).

CALVAL X: Calibration temperature 1 or 2 or 3 (actual value CT)

The „calibration value“ is determined with a temperature measuring device and saved under menu item >CALVAL<.

(value flashes \overline{X})

7.6.3. Example: 3-point calibration for internal control.

In the temperature range from 80 °C to 160 °C the calibration curve of the temperature sensor (TT) is to be adjusted to the actual temperatures at measuring point (CT).

This is only an example. The transparent thermostat CT 72 / P can only be set to maximum 60.00 °C.


1. Set working temperature setpoint – SETPNT:

Refer to „Direct temperature setting “ page 68

- By pressing the key the circulator switches to the active >SETPOINT< see example on the left: >SETPNT / 1 25.00°C<. The integer digits flash (Example: <25>).
- Change the value to 80.00 °C by pressing the keys and and confirm by pressing the key . The decimal digits flash.
- Confirm once more by pressing the key .
- The bath is heated up.
Wait for approx. 5 minutes until the temperature is constant.

2. Reading of temperature measuring device

Read the value of measuring point CT on the device and enter under menu item >CALVAL X< by using the keypad.

- >CALVAL 1< (79.70 °C)
- >CALVAL 2< (119.50 °C)
- >CALVAL 3< (159.30 °C)

3. Calibration

Menu level 1






- Press the **ESC** key if parameter is to be retained. Correction function for parameters or values (prior to ENTER).
- Setting is required only for the first calibration point.

Level 2	Parameter level	
		<p>To prepare a ATC calibration. adjust >NO<</p> <ul style="list-style-type: none"> The parameter flash, switch over with + .

Level 2	Parameter level	
<pre> ATE <I> TYPE ▼ ▲ </pre>	<pre> TYPE 3-POINT ✱ </pre> 	<ul style="list-style-type: none"> The parameter flash, switch over with ▼ ▲ + ←. <p>A >3. Point< calibration is performed.</p>
<pre> ATE <I> TMPVAL1 ▼ ▲ ATE <I> CALVAL1 ▼ ▲ </pre>	<pre> TMPVAL1 8000 ESC </pre> <pre> CALVAL1 79.70 ✱ </pre>	<p>The value >TMPVAL< is only indicated In addition the measured value >CALVAL X< is saved during the following step</p> <ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ (79) + ← Decimal digits flash, set by pressing ▼ ▲ (70) + ← <p>The first of 3 points is calibrated.</p>
Return to 2. Set working temperature value SETPNT: 120.00 °C		
<pre> ATE <I> TMPVAL2 ▼ ▲ </pre>	<pre> TMPVAL2 12000 ESC </pre>	<p>The value is only indicated</p>
<pre> ATE <I> CALVAL2 ▼ ▲ </pre>	<pre> CALVAL2 119.50 ✱ </pre>	<ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ (119) + ← Decimal digits flash, set by pressing ▼ ▲ (50) + ← <p>The second of 3 points is calibrated.</p>
Return to 2. set working temperature value SETPNT: 160.00 °C		
<pre> ATE <I> TMPVAL3 ▼ ▲ </pre>	<pre> TMPVAL3 16000 ESC </pre>	<p>The value is only indicated.</p>
<pre> ATE <I> CALVAL3 ▼ ▲ </pre>	<pre> CALVAL3 159.30 ✱ </pre>	<ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ (159) + ← Decimal digits flash, set by pressing ▼ ▲ (30) + ← <p>The 3-point calibration is completed</p>
<h4>4. Return to standard operation</h4>		
<pre> ATE <I> STATUS </pre>	<pre> STATUS YES ✱ OK </pre>	<ul style="list-style-type: none"> Set >YES< after calibration. (Standard operation)

7.7. MENU LIMITS

<p>Menu- Level 1</p> 	<p>The transparent thermostat CT 72 / P can only be set to maximum 60.00 °C. The configuration does not permit higher values.</p> <p>ⓘ Press the ESC key if parameter is to be retained. Correction function for parameters or values (prior to ENTER).</p>	
<p>Level 2</p>	<p>Parameter level</p>	
		<ul style="list-style-type: none"> Integer digits flash, set by pressing  Decimal digits flash, set by pressing 
		<ul style="list-style-type: none"> Integer digits flash, set by pressing  Decimal digits flash, set by pressing 
		<ul style="list-style-type: none"> The value flashes, set by pressing 
		<ul style="list-style-type: none"> The value flashes, set by pressing mit 

7.7.1. Limits for internal control

C-TYPE
INT

SETPOINT MAX / MIN – Maximum and minimum setpoint
Restriction of the adjustable temperature range

The limitation of the operating temperature range effects the temperature setting in the menu with the key **T**.

SETMAX
200.00

Only setting of working temperatures which lie within the determined limits is possible

SETMIN
-50.00

Existing settings for SETPNT 1, -2, -3, as well as those for >OVERTMP< and >SUBTMP < (refer to page 70), are automatically deferred into the limit range.

Setting range: -94,90 °C ... +200,0 °C

SET MAX > SET MIN

Interchange of values is not possible.

Set maximum heating / cooling

The heating and cooling capacity of the unit are adjustable. 100 % corresponds to the technical specification of the equipment.

Setting range:

HEAT MAX – 0 to 100 % in 1 % steps

COOLING MAX – 0 to 100 % in 1 % steps

HEATMAX
100

COOLMAX
0

8. Troubleshooting guide / error messages

-OFF-
ALARM
CODE 01

Alarm with complete shutdown:

If one of the following failures occur a complete, all-pole shutdown of the heater and circulating pump is effected.

„“ lights up and a continuous signal sounds.

The code for the cause of alarm is indicated on the VFD COMFORT-DISPLAY.

 + 
WARNING
CODE 03

Alarm without shutdown:

The code for the cause of alarm is indicated on the VFD COMFORT-DISPLAY. The warning signal sounds in regular intervals. The messages appear every 10 seconds.

Press the key  to stop the signal

ALARM
CODE 01

Error message with ticker: >LOW LEVEL ALARM-FILL MEDIUM <

Low level alarm

The circulator is operated without or insufficient bath fluid.

Switch the unit off with the mains switch, refill bath fluid and switch on!

Tube breakage has occurred (insufficient filling level of bath fluid caused by pumping-out)

Replace the tubing and refill bath liquid.

The float is defect (e.g. transport damage).

Repair by authorized JULABO service personnel.

ALARM
CODE 02

Error message with ticker:

> REFRIGERATOR ALARM-CHECK CONNECTION <

During the self-test after switch-on a short –circuit is detected between pin 2 and pin 4 of the control line or the control line was disconnected during operation.

Reconnect the control line or repair short-circuit.

WARNING
CODE 03
ALARM
CODE 03

Error message with ticker:

> EXCESS TEMPERATURE WARNING-CHECK LIMITS <

Excess temperature warning

or

Excess temperature alarm

Type of warning: set to >warning< or >alarm< (refer to page 70)

WARNING
CODE 04
ALARM
CODE 04

Error message with ticker:

> LOW TEMPERATURE WARNING-CHECK LIMITS <

Low temperature warning

or

Low temperature alarm.

Typ of warning: set to >warning< or >alarm< (refer to page 70)

ALARM
CODE 05

Error message with ticker:

> WORKING SENSOR ALARM-CALL SERVICE <

Cable of working temperature sensor is disconnected or short-circuited.

ALARM
CODE 06

Error message with ticker:

>SENSOR DIFFERENCE ALARM-CHECK VISCOSITY AND PUMP STAGE<

Defect of working or excess temperature protector.

Working temperature sensor and excess temperature protector report a temperature difference of more than 35°C.

ALARM
CODE 07

Error message with ticker:
> INTERNAL HARDWARE ERROR-CALL SERVICE <
Other errors

ALARM
CODE 12


Error in A/D converter

ALARM
CODE 14

Error message with ticker:
> EXCESS TEMPERATURE PROTECTOR ALARM-CHECK ADJUSTMENT <
Excess temperature protector defect.
The protection temperature is below the set working temperature setpoint.
Set the protection temperature to a higher value.

ALARM
CODE 33

Error message with ticker:
> SAFETY SENSOR ALARM-CALL SERVICE <
The cable of the excess temperature protector has been disconnected or short-circuited

By quickly switching off and restarting the unit the alarm is cancelled.
If the error occurs once more after the restart, a remote diagnosis is required.


[Err
PRESS
OK

Error message with ticker:
> CONFIGURATION ERROR-PRESS OK<
The configuration of the circulator does not correspond with its current application.
Press the  key for a non-recurring, automatic change of the configuration.

In this case please call the Technical Service or an authorized dealer.

Disturbances that are not indicated.

The electronic pump motor is overload-protected by an electronic current limiter. If viscosity of the bath fluid is or becomes too high, the motor stops running.

Mains circuit breakers (resettable) 15 A


9. Electrical connections

Notice

Use shielded cables only.

The shield of the connecting cable is electrically connected to the plug housing.

RS232 serial interface

This port can be used to connect a computer with an RS232 cable for remote control of the circulator.

Designation:

Pin 2	RxD	Receive Data
Pin 3	TxD	Transmit Data
Pin 5	0 V	Signal GND
Pin 7	RTS	Request to send
Pin 8	CTS	Clear to send


Pin 1; 4; 6, 9 reserved - do not use!

Interface connecting cable RS232:

CT72 (9pol)			PC (9pol)	
Pin 2	RxD	⇔	Pin 3	TxD
Pin 3	TxD	⇔	Pin 2	RxD
Pin 5	GND	⇔	Pin 5	GND
Pin 7	RTS	⇔	Pin 8	CTS
Pin 8	CTS	⇔	Pin 7	RTS

Connection of flow-through cooler CK300 / CK310

The flow-through cooler **CK300 / CK310** is connected to pins 1 and 3. For the connection use the shielded control cable supplied with the cooler.


Circuit:	Operation	= relay powered
	Alarm	= relay not powered

Pin assignment:

Pin	Signal
1	+24 V (I max. current 25 mA)
2	0 V
3	Alarm relay
4	Reserved - do not use!
5	Cooling pulse

10. Remote control

10.1. Setup for remote control

- Check the interface parameters for both interfaces (on circulator and PC) and make sure they match. (Serial interface refer to page 83)
- In the menu > MENU CONFIG < set the menu item > REMOTE < to > ON < .
- Connect both units with an interface cable..

Like all parameters which can be entered through the keypad, interface parameters are stored in memory even after the thermostat is turned off.

10.2. Communication with a PC or a superordinated data system

If the thermostat is put into remote control mode via the configuration level, the VFD COMFORT-DISPLAY will read "R -OFF-" = REMOTE STOP. The circulator is now operated via the computer.

In general, the computer (master) sends commands to the circulator (slave).

The thermostat sends data (including error messages) only when the computer sends a query.

In remote control mode, the start command and all values to be set must be resent by the PC via the interface in case of a power interruption.

AUTOSTART is not possible.

A transfer sequence consists of:

- command
- space (↔; Hex: 20)
- parameter (decimal separation with a period)
- end of file (↵; Hex: 0D)

The commands are divided into **in** and **out** commands.

in commands: retrieve parameters

out commands: set parameters

Important times for a command transmission:

To ensure a safe data transfer, the time gap between two commands should be at least 250 ms.

The circulator automatically responds to an **in** command with a data string followed by a LF (Line Feed). The next command should only be sent after 10 ms.

The out commands are valid only in remote control mode.

Command to set the working temperature > SETPNT 1< to 55.5 °C
 out_sp_00 ⇔ 55.5.↓

Command to retrieve the working temperature > SETPNT 1< in_sp_00.↓

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in commands: Asking for parameters or temperature values to be displayed.

Command	Parameter	Response of circulator
version	none	Number of software version (V X.xx)
status	none	Status message, error message (see page 96)
in_pv_00	none	Actual bath temperature.
in_pv_01	none	Heating power being used (%).
in_pv_03	none	Temperature value registered by the safety sensor.
in_pv_04	none	Setpoint temperature of the excess temperature protection
in_sp_01	none	Working temperature „SETPNT 1“
in_sp_02	none	High temperature warning limit „OverTemp“
in_sp_03	none	Low temperature warning limit „SubTemp“
in_par_01	none	Te - Time constant of the external bath.
in_par_02	none	Si - Internal slope
in_par_03	none	Ti - Time constant of the internal bath.
in_par_05	none	Factor pk/ph0: Ratio of max. cooling capacity versus max. heating capacity
in_par_06	none	Xp control parameter of the internal controller.
in_par_07	none	Tn control parameter of the internal controller.
in_par_08	none	Tv control parameter of the internal controller.
in_mode_01	none	Selected setpoint: 0 = SETPNT 1 1 = SETPNT 2 2 = SETPNT 3
in_mode_05	none	Circulator in Stop/Start condition: 0 = Stop 1 = Start
in_mode_08	none	Adjusted control dynamics 0 = aperiodic 1 = standard

10.4. Status messages

Status messages	Description
00 MANUAL STOP	Thermostat in „OFF“ state.
01 MANUAL START	Thermostat in keypad control mode.
02 REMOTE STOP	Thermostat in „r OFF“ state.
04 REMOTE START	Thermostat in remote control mode.

10.5. Error messages

Error messages	Description
-01 TEMP / LEVEL ALARM	Low liquid level alarm.
-02 REFRIGERATOR ALARM	Control cable of the flow-through cooler is short-circuited or interrupted.
-03 EXCESS TEMPERATURE WARNING	High temperature warning.
-04 LOW TEMPERATURE WARNING	Low temperature warning.
-05 TEMPERATURE MEASUREMENT ALARM	Working temperature sensor short-circuited or interrupted.
-06 SENSOR DIFFERENCE ALARM	Sensor difference alarm. Working temperature and safety sensors report a temperature difference of more than 35 °C.
-07 I²C-BUS ERROR	Internal error when reading or writing the I ² C bus.
-08 INVALID COMMAND	Invalid command.
-09 COMMAND NOT ALLOWED IN CURRENT OPERATING MODE	Invalid command in current operating mode.
-10 VALUE TOO SMALL	Entered value too small.
-11 VALUE TOO LARGE	Entered value too large.
-12 WARNING : VALUE EXCEEDS TEMPERATURE LIMITS	Error in A/D converter.
-14 EXCESS TEMPERATURE PROTECTOR ALARM	Excess temperature protector alarm
-30 CONFIGURATION ERROR: CONFIRM BY PRESSING <OK> ON CIRCULATOR	The configuration of the circulator does not conform to its present use. Press  to automatically perform a single modification of the configuration.
-33 SAFETY SENSOR ALARM	Excess temperature sensor short-circuited or interrupted.

11. Cleaning / repairing the unit

Caution

- Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit.
- Prevent humidity from entering into the circulator.
- Electrical connections and any other work must be performed by qualified personnel only.
- Acryl-baths are not resistant to solvents.
Do not use detergents containing alcohol or solvents.

Cleaning:

For cleaning the bath tank and the immersed parts of the circulator, use low surface tension water (e.g., soap suds).

Clean the outside of the unit using a wet cloth and low surface tension water.

The Thermostat is designed for continuous operation under normal conditions.

Periodic maintenance is not required.

The tank should be filled only with a bath fluid recommended by SI Analytics GmbH.

To avoid contamination, it is essential to change the bath fluid from time to time.

Repairs:

Before asking for a service technician or returning a circulator for repair, please contact our service station.

SI Analytics GmbH
Hattenbergstrasse 10
D-55122 Mainz
Telefon +49 61 31 / 66 51 11
Telefax +49 61 31 / 66 50 01
E-Mail: si-analytics@xyleminc.com
www.si-analytics.com

When returning a unit, take care of careful and adequate packing. SI Analytics GmbH is not responsible for damages that might occur from insufficient packing.

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