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An important property of liquids is their viscosity. A measure for the viscosity is the flow time, determined at a specified temperature, using a flow cup.

Application/application areas

- for paint and varnish manufacturers, its customers and test laboratories
- the test can be carried out in the lab and during production
- for paints, varnishes and similar liquids


Standard extent of delivery

- 1 Flow cup/Immersion flow cup with certificate of manufacturer

Options

- ACC 072 three-foot stand completely made of stainless steel, levelling, adjustable feet, with spirit level, height of 260 mm
- ZKZ 005 certificate of calibration

eingebaute Libelle /
spirit level

Handling Flow cups

- choose a flow cup that will give a flow time of between 30 s and 100 s (depending on the standard)
- put the flow cup into the levelled three-foot stand
- close the orifice of the cup by a finger and fill the cup with the liquid under test, avoiding the formation of air bubbles
- slide a glass plate over the rim of the cup so that no air bubbles will be formed between the glass plate and the surface of the sample
- draw the glass plate horizontally away, remove the finger from the orifice of the cup and simultaneously start the timing-device
- stop the timing-device as soon as the stream of the sample first breaks close to the orifice

Handling Immersion flow cups

- immerse the immersion flow cup entirely below the surface of the liquid under test
- withdraw the cup vertically and start the timing-device
- stop the timing-device as soon as the stream of the sample first breaks close to the orifice
- the average out of 3 measurements will give the flow time

Technical specification

ZFC 3010	Flow cup in accordance with ASTM D 1200 (Ford) with fixed jet* No. 1, 2, 3, 4 or 5.
ZFC 3011	Flow cup with fixed jet* No. 1, 2, 3, 4**, 5, 6, 7 or 8.
ZFC 3012	Immersion flow cup with fixed jet* No. 1, 2, 3, 4**, 5, 6, 7 or 8.
ZFC 3013	Flow cup in accordance with DIN EN ISO 2431 with fixed jet* No. 3, 4, 5 or 6.
ZFC 3014	Immersion flow cup in accordance with DIN EN ISO 2431 with fixed jet* No. 3, 4, 5 or 6.
ZFC 3015	Immersion flow cup in accordance with ASTM D 1200 (Ford) with fixed jet* No. 1, 2, 3, 4, 5 or 6.

* Please specify jet No. when ordering.

** Jet No. 4 in acc. with DIN 53211, withdrawn in October 1996

material: cup in anodized aluminium,
jet in stainless steel

warranty: 2 years

ZFC 3010 - 3015 Flow cups developed and manufactured in Switzerland by Zehntner GmbH Testing Instruments. Special customer versions are available on request. Subject to technical changes.