EN

# Konelab™/ T Series CO<sub>2</sub> (Bicarbonate)

REF 981371 2 x 100 ml

**THIS PACKAGE INSERT IS APPLICABLE FOR USE  
OUTSIDE THE US. ANY REFERENCE TO THE KONELAB  
SYSTEMS ALSO REFERS TO THE T SERIES.**

## INTENDED USE

For the *in vitro* quantitative determination of the CO<sub>2</sub> (bicarbonate) concentration in human serum or plasma on Konelab analyzers. All test results must be interpreted with regard to the clinical context.



## SUMMARY (1)

Total carbon dioxide in plasma consist of CO<sub>2</sub> dissolved in an aqueous solution, CO<sub>2</sub> bound loosely to amine groups in proteins (carbamino compounds), CO<sub>2</sub> that exists as HCO<sub>3</sub><sup>-</sup> or CO<sub>3</sub><sup>2-</sup> ions, and other undissociated bicarbonates as well as carbonic acid.

Approximately 90 % of carbon dioxide present in serum or plasma is in the form of bicarbonate. Alterations of HCO<sub>3</sub><sup>-</sup> and CO<sub>2</sub> dissolved in plasma are characteristic of acid-base imbalance. Bicarbonate's value has the most significance in the context of other electrolyte values and with blood gases and pH values.

## PRINCIPLE OF THE PROCEDURE

The CO<sub>2</sub> (in the form of bicarbonate ions) reacts with phosphoenolpyruvate to form oxaloacetate and phosphate. The reaction is catalyzed by phosphoenolpyruvate carboxylase (PEPC). Malate dehydrogenase (MDH) then catalyzes the reduction of oxaloacetate to malate and the oxidation of NADH to NAD<sup>+</sup>. The resulting decrease in absorbance at 380 nm is proportional to the amount of bicarbonate present in the sample. (2, 3)

## REAGENT INFORMATION

2 x 100 ml dry powder reagent

## Concentrations

### in reconstituted reagent:

Phosphoenolpyruvate	8.0 mmol/l
NADH	1.6 mmol/l
Phosphoenolpyruvate carboxylase (microbial)	> 1000 U/l
MDH (microbial)	> 200 U/l
Buffer	66 mmol/l
NaN <sub>3</sub>	< 0.1 %
Stabilizers	
pH 8.05 at 20 °C	

## Precautions

For *in vitro* diagnostic use only. Exercise the normal precautions required for handling all laboratory reagents.

The reagent contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes. Sodium azide may react with lead and copper in plumbing to form highly explosive metal azides. On disposal flush with large volumes of water.

See separate sheet for Risk and Safety Phrases (R21/22-52/53; S28-45).

## Preparation

Dissolve the dry powder reagent with 100 ml of freshly distilled or deionized water.

If necessary, dissolve the reagent by mixing gently – AVOID SHAKING.

Use the empty lot specific vial in the kit.

**Note:** Check that there are no bubbles in the bottleneck or on the surface of the reagent when you insert the reagent vials or vessels in the Konelab analyzer.

## Storage and Stability

Reagents in unopened vials are stable at 2...8 °C until the expiration date printed on the label.

The reconstituted reagent is stable for 6 months when stored capped at 2...8 °C.

Refer to the Application Notes of your Konelab analyzer for the on board stability of reagents.

Indications of reagent deterioration are turbidity, absorbance < 1.0 at 380 nm (1 cm), and/or failure to recover control values within the assigned range.

## SPECIMEN COLLECTION

### Sample Type

Serum or heparin plasma can be used. Heparin plasma gives slightly lower values than serum. Keep sample tubes closed.

### Precautions

Human samples should be handled and disposed of as if they were potentially infectious.

### Storage (4)

The sample can be stored for 1 day at 20...25 °C (1 hour after opening the tube), for 7 days at 4...8 °C or for 2 weeks at -20 °C.

## TEST PROCEDURE

Refer to the Reference Manual and Application Notes for an automated procedure on your Konelab analyzer. Any application which has not been validated by Thermo Fisher Scientific Oy cannot be performance guaranteed and therefore must be evaluated by the user.

## Materials provided

Reagent as described above. Two empty lot specific vials.

## Materials required but not provided

Calibrator as indicated below.

## Calibration

Use sCal, code 981831, according to the instructions provided for your Konelab analyzer.

## Traceability:

Refer to the package insert of sCal.

## Quality Control

Use quality control samples at least once a day and after each calibration and every time a new bottle of reagent is used.

Available controls:

Nortrol, code: 981043


Abtrol, code: 981044

The Control intervals and limits must be adapted to the individual laboratory requirements. The results of the quality control sample should fall within the limits pre-set by the laboratory.

## CALCULATION OF RESULTS

The results are calculated automatically by the Konelab analyzer using a calibration curve.

### Calibration Curve (example)

Konelab 20/30/60. The calibration curve is lot dependent.

## LIMITATIONS OF THE PROCEDURE

### Interference

Criterion: Recovery within ± 10% of initial values.

Bilirubin: No interference found up to 1000 µmol/l (58.5 mg/dl).

Hemolysis: No interference found up to 10 g/l of hemoglobin in hemolysate.

Lipemia: No interference found up to 10 g/l of Intralipid® (trademark of Fresenius Kabi AB).

There is a poor correlation between turbidity and triglycerides concentration.

For other interfering substances, please refer to the reference 5.

### EXPECTED VALUES (1)

Adult: 22 – 29 mmol/l (22 – 29 mEq/l)

The quoted values should serve as a guide only. It is recommended that each laboratory verify this range or derives a reference interval for the population that it serves.

### MEASURING RANGE

5 – 40 mmol/l (5 – 40 mEq/l)

Extended measuring range after secondary dilution:

5 – 120 mmol/l (5 – 120 mEq/l).

## PERFORMANCE CHARACTERISTICS

### Detection limit

1 mmol/l (1 mEq/l)

The detection limit represents the lowest measurable concentration/activity that can be distinguished from zero. It is calculated as the concentration of zero sample + 3 SD (within run, n=24).

### Imprecision

	Mean 20 mmol/l		Mean 25 mmol/l		Mean 33 mmol/l	
	SD	CV%	SD	CV%	SD	CV%
Within run	0.4	2.2	0.7	3.1	1.2	3.7
Between day	0.4	2.1	1.3	5.4	0.7	2.1
Total	0.7	3.7	1.6	6.7	1.5	4.6

A precision study was performed using the NCCLS Document EP5-A as a guideline and Konelab 30 during 10 days, with the number of measurements being n = 40.

### Method comparison

A comparison study was performed using the NCCLS Document EP9-A as a guideline and a commercially available enzymatic method as a reference.

Linear regression (result unit mmol/l):

$$y = 1.10x - 0.5$$

$$r = 0.946$$

$$n = 177$$

The sample concentrations were between 13 and 62 mmol/l.

The results obtained in individual laboratories may differ from the given performance data.

## BIBLIOGRAPHY

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCPress, Washington, D.C., 2000, 3-117 – 3-120.

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**LITERATURHINWEISE**

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, 5th Edition, AACC Press, Washington, D.C., 2000, 3-117 – 3-120.

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#### BIBLIOGRAPHIE

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.

3. Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
4. Guder WG, Narayanan S, Wisser H, Zawta B, List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
5. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACC Press, Washington, D.C., 2000, 3-117 – 3-120.

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#### SEZNAM LITERATURY

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, L.J, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCPress, Washington, D.C., 2000, 3-117 – 3-120.

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**BIBLIOGRAFIA**

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG. Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: >From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG. Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, 5th Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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**KIRJANDUS**

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analyses; Preanalytical Variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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**BIBLIOGRAFIA**

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, 5th Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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**BIBLIOGRAFIJA**

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, pp. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO2 in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. Analītiskās saraksts: Pirmsanalīzes vērtības. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG, Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCC Press, Washington, D.C., 2000, 3-117 – 3-120.

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## BIBLIOGRAFIA

- Burtis, CA and Ashwood, E R (ed.), Tietz Fundamentals of Clinical Chemistry, 5th edition, W B Saunders Company, Philadelphia, 2001, p. 732, 966.
- Norris, KA, Atkinson, AR and Smith, WG. Colorimetric Enzymatic Determination of Serum Total Carbon Dioxide, as Applied to the Vickers Multichannel 300 Discrete Analyzer, Clin. Chem. 21(8), 1975, pp. 1093-1101.
- Forrester, RL, Wataji, LJ, Silverman, DA and Pierre KJ, Enzymatic Method for Determination of CO<sub>2</sub> in Serum, Clin. Chem. 22(2), 1976, pp. 243-245.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Fifth Edition, AACCPress, Washington, D.C., 2000, 3-117 - 3-120.

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