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Instructions for Use EN

GOeasy Provi Temp - AUTOMIX 4:1

Self cure esthetic temporary crown and bridge material

GOeasy Provi Temp is an automatically mixed two component material in Automix-cartridges based on multifunctional methacrylic esters. The temporary material is a very easy and comfortable material for creating short- and long-term temporary crowns or bridges, onlays and veneers. Because of its flexibility, the material is especially suitable for longer bridge spans.

GOeasy Provi Temp is free of methylmethacrylate. Its temperature derived while curing is lower than 40 °C (104°F). As a temporary crown or bridge it protects the prepared teeth against external influences and preserves the occlusion. It shows increasing transverse strength, good abrasion resistance and low polymerization shrinkage. It fits perfect. Showing good polishability, good color stability and fluorescence it looks perfect.

Contraindications

If the patient is known to be allergic to acrylates.

Side effects

In singular cases, GOeasy Provi Temp may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Application

1. Impression taking

Before preparing of stumps or before a planned extraction of a tooth a situation impression is made by addition curing silicones (long storage stability!) or by alginates. For better stability of the temporary interdental areas are carved out. In molar areas with teeth absent it may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

Note:

In the silicone impression block out undercuts and if necessary cut grooves into the impression.

2. Preparing of the Automix-cartridge

Remove the cap of the Automix-cartridge and throw it away (**do not use it again!**). Attach one of the supplied 4:1 mixing cannulas. Make sure that the guidance of the Automix-cartridge is aligned with that of the mixing cannula and turn the cannula 90° clockwise until it locks in position. The application gun loaded with the prepared cartridge and is ready for application.

Note:

The initial extrusion from the mixing cannula (about

the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix.

Leave the used mixing cannula on the cartridge. It serves as a cap.

3. Forming of the temporary crowns / bridges

GOeasy Provi Temp is automatically mixed when dispensed with slight and even pressure directly into the situation impression made before. Filling should occur from bottom upward to prevent voids.

Working time of GOeasy Provi Temp is 50 seconds (at room temperature 23 °C/74 °F)

a. Application in the mouth

Load the situation impression with GOeasy Provi Temp. Seat the impression within the working time (50 seconds at 23 °C/74 °F) onto the prepared areas of the teeth.

After 1-2 minutes (setting time in mouth 37 °C/98 °F) the material shows hardened but still elastic condition and can be removed from the teeth together with the situation impression.

b. Application on the model

Load the situation impression with GOeasy Provi Temp. Seat the impression within the working time (50 seconds at 23 °C/74 °F) onto the prepared areas of the model.

After 3-4 minutes (setting time at 23 °C/74 °F) the material shows a hardened but still elastic condition and can be removed from the model together with the situation impression.

Note:

The setting reaction has to be checked with excess material intra-orally (e.g. with a scaler) resp. on the model. Temperature in mouth has a significant effect on setting reaction and the temporary can only be removed without destruction during the elastic state.

4. Post curing and finishing

If possible, leave the temporary in the situation impression during post curing.

Optimally the temporary is post cured in warm (45-55 °C /113-131 °F) water (e.g. in a hot cure polymerization device) until reaching its final hardness (approx. 4:30 min from start of mixing). Post curing at room temperature is completed after about 6 minutes from start of mixing.

After removal from the impression excess material and proximal undercuts are removed. Then the temporary acrylic can be worked out with rotary instruments and can be polished to high gloss.

Do not breathe polishing dust; use suitable mouth protective device, safety glasses and aspiration!

Note:

The oxygen inhibited smear layer on the surface caused by the acrylic system can not be fully avoided and should be removed before working out. It can easily be removed by alcohol or other suitable solvents.

5. Cementing of the temporary

GOeasy Provi Temp temporary crowns or bridges preferably should be cemented with an eugenol-free temporary cement (e.g. GOeasy Provi Temp Cem). In case of using eugenol-containing cements it should be considered that later used acrylate based cements could be hindered in hardening (inhibition of curing reaction by traces of remained eugenol).

6. Repairs

GOeasy Provi Temp temporaries show high mechanical strength. However, if a temporary breaks the following procedures are recommended:

a) Fracture shortly after production

Both ends of the fracture are refixed with freshly extruded GOeasy Provi Temp.

b) Fracture of longer existing temporary

The areas of fracture are cleaned and roughened and provided with some mechanical retentions. Thus prepared, the fractured areas can be joint with freshly extruded GOeasy Provi Temp. To enhance the curing place for some minutes in warm water.

7. Special hints

Non hardened GOeasy Provi Temp can easily be removed with alcohol or other suitable solvents.

8. Precautions

GOeasy Provi Temp is free of methylmethacrylate but contains other methacrylates.

With susceptible patients, sensitization to GOeasy Provi Temp cannot be excluded. GOeasy Provi Temp should not be used any more, if allergic reactions are observed. Do not use for patients with allergic reactions against acrylates.

Avoid contact with skin, mucous membrane and eyes. If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.

9. Storage

Do not store above 25 °C (77 °F). Do not use after expiry date.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no

other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law.

P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!


For dental professional use only!

The instructions for use has to be kept for the duration of the application.

Caution:

Federal law restricts this device to sale by or on the order of a dentist.

Time Table of application and curing of GOeasy Provi Temp

If the completion of setting takes place at room temperature, trimming / polishing can be done after 6 minutes (from start of mixing)