
Torsionsgerät 8405740

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07/06 JS

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Instruction sheet

07/06 JS

- 1 Torsion rods
- 2 Hemp string
- 3 Spring dynamometer
- 4 Weights
- 5 Counter-bearing with knurled screws
- 6 Pointer
- 7 Pulley
- 8 Circular scale

1. Description

The torsion apparatus allows for static or dynamic measurement of elastic properties of metal rods of different lengths and thickness and the corresponding torsional oscillations up to the limit of elasticity.

- 1 Pulley, 30 mm Ø, with 4 mm stem
- 1 Holder for pulley
- 2 Weights, 50 g, with 4 mm stem
- 1 Hemp string, 100 m
- 1 Spring dynamometer, 5 N

2. Scope of delivery

- 1 Basic torsion equipment, consisting of:
 - 1 Circular scale with bearing and chuck,
 - 1 pointer, 1 counter-bearing with chuck
- 1 Torsion rod, steel, 500 mm x 2 mm Ø
- 1 Torsion rod, brass, 500 mm x 2 mm Ø
- 1 Torsion rod, copper, 500 mm x 2 mm Ø
- 1 Torsion rod, Al, 300 mm x 2 mm Ø
- 1 Torsion rod, Al, 400 mm x 2 mm Ø
- 1 Torsion rod, Al, 500 mm x 2 mm Ø
- 1 Torsion rod, Al, 500 mm x 3 mm Ø
- 1 Torsion rod, Al, 500 mm x 4 mm Ø

3. Technical data

Circular scale

Dimensions: 5 mm x 230 mm Ø
Weight: 545 g

Counter-bearing with chuck

Length: 120 mm
Weight: 115 g

Pointer

Dimensions: 20 × 30 × 45 mm

4. Additionally required

For setting up:

1 H-shaped base	8611130
3 Stand rods, 750 mm, 12 mm Ø	8611340
1 Square clamp	8613180
1 Rod with hook	8615110

For dynamic measurements:

1 Stopwatch	8533135
-------------	---------

5. Sample experiments

5.1 Static measurement in the elastic range

- Set up experiment as in Fig. 1.
- Set up circular scale and chuck opposite its counter-bearing and at the same both bearings are aligned horizontally and that the orientation towards the opposing chuck is accurate.
- Secure torsion rod into both chucks without twisting it. Adjust the H-shaped base to the correct length for this.
- Calibrate the zero point for the angular scale by slightly loosening the knurled screw on the counter-bearing and rotating the rod. Tighten knurled screw again.
- Fasten pulley to stand rod. Make sure that the pulley is in the same vertical plane as the circular scale.
- Suspend the spring dynamometer over the pulley using the square clamp and the rod with a hook.
- Fasten hemp string to circular scale, thread through the groove and around the pulley, and fasten to spring scale.
- Increase the force acting on the circular scale (i.e., the torque) in 0.1 N steps by shifting the pulley along the stand rod and read the resulting angle of twist or torsional angle.
- Plot a graph of torsional angle and torque.

The experiment set-up guarantees that the force will always act vertically, since the string runs at a tangent to the circular scale. This means the simple equation

$$M = F \cdot r$$

can be used to calculate the effective torque.

It is imperative that the force is not increased too much in order not to exceed the elastic range of deformation (different from material to material). If the circular scale returns to zero after relieving the force, the elastic range has not been exceeded.

Within the elastic range, the effective torque M is proportional to the torsional angle φ . This relation can be described as follows:

$$M = D_w \cdot \varphi.$$

The proportionality factor D_w is called the torsional constant or torsion coefficient. This physical torsion coefficient depends on the material and the dimensions of the wire. It can be derived from the slope of the line in the diagram, or the quotient of torque M and the torsional angle φ , respectively:

$$D_w = M/\varphi$$

The torsion coefficient D_w of a material decreases with the length of the wire, as well as with its diameter.

5.1.1 Determining the limit of elasticity (yield point)

- Set up the torsion apparatus as in 5.1.
- By moving the pulley, increase the force in steps of 0.1 N and read off the resulting torsional angles.
- Each time the force is increased, release the force on the circular scale by taking away the dynamometer and check that the the scale returns to an angle of 0°.

If the circular scale does not return to an angle of 0° after relieving the force and some slight deformation remains, then the limit of elasticity, the yield point, has been passed.

The resulting torsional angles and torques depend on material and geometry (length and diameter of the wire investigated).

Once the limit of elasticity has been exceeded, the measured values do not follow a straight line anymore. Plastic (permanent) deformation of the wire results.

5.2 Static measurement in the plastic range – mechanical hysteresis

- Set up the torsion apparatus as in 5.1.
- By moving the pulley, increase the force in steps of 0.1 N and read off the resulting torsional angles.
- After obtaining a torsional angle of 180°, start to decrease the force, again in steps of 0.1 N.

Even when the force is removed entirely, an angle of twist is still observed. The wire has been physically deformed by a corresponding degree.

- Now reverse the direction of the string around the circular scale, so that the force is exerted by the dynamometer in the opposite direction.
- Once more, increase the force in steps of 0.1 N and read off the resulting torsional angles.

In order to get back to the 0° mark (original position) now, an extra counteracting force is required.

- Deform wire, step by step, up to a torsional angle of -180°.

- After reaching this angle, reverse the direction of the string and the direction of the force again.
- Increase force, step by step, up to a torsional angle of 180°.

Once this is done, a complete hysteresis curve should have been obtained. The characteristic values of this curve are the torsional angle in the relieved state (meaning with zero force acting, i.e. where the curve crosses the x-axis), and the extra force required to pull the wire back to the original value of 0° (where the curve crosses the y-axis).

5.3 Dynamic measurement (torsional pendulum)

- Set up torsion apparatus as in 5.1, but without the spring scale.
- Move the circular scale from the point of origin to an angle of approximately 25° and release it.
- Measure the time taken for 10 full, unimpeded torsional oscillations and take this as a base to calculate the oscillation period of the system.
- Conduct an initial measurement with the circular scale without adding any weight. In subsequent measurements, add the two weights with the 4 mm stems symmetrical either side the axis of the circular scale and repeat the measurement. Start at the innermost stem position, and move the extra weight one position outward for each measurement.
- Write down the resulting times, and the distance of the weights from the rotating axis.

The oscillation period T of the torsional pendulum can be calculated as:

$$T = 2\pi \cdot \sqrt{\frac{J}{D}} \quad (1)$$

where J is the moment of inertia of the oscillator, and D is the torsion coefficient. The moment of inertia J_{tot} is the sum of the moment of inertia J_K of the circular scale and the moment of inertia J_z of the added weights.

$$J_{\text{tot}} = J_K + J_z \quad (2)$$

If the additional weights are considered, their moment of inertia can be calculated as:

$$J_z = 2mr^2 \quad (3)$$

By squaring equations (1), (2) and (3), this results in:

$$T^2 = 4\pi^2 \frac{(J_K + 2mr^2)}{D} = \frac{4\pi^2}{D} J_K + \frac{8\pi^2 m}{D} r^2 \quad (4)$$

At first glance, the moment of inertia of the circular scale remains unknown. According the equation (4), the plot $T^2=f(r^2)$ must result in a straight line, which, due to the circular scale's moment of inertia, does not intersect with the origin.

Rearranging equation (1) for D and substituting (2) results in:

$$D = (J_K + J_z) \left(\frac{2\pi}{T} \right)^2 \quad (5)$$

J_K for the circular scale is always constant, J_z changes, in accordance with (3), corresponding to the distance of the added weights from the rotation axis. By taking two measurements for the oscillation period for different J_z , J_K can be eliminated from (5), and D can be calculated with the following equation:

$$D = \frac{(J_{z2} - J_{z1}) \cdot 4\pi^2}{T_2^2 - T_1^2} \quad (6)$$

Because of the subtraction, there needs to be a large difference between J_{z1} and J_{z2} , to keep the error small. An error analysis provides much more accurate values for dynamic measurements than for static measurements. In principle, dynamic measurement results in smaller values for D , since the friction in static measurements has the effect of making the torsion coefficient D appear bigger.


Fig. 1 Experiment set-up for static measurements



Fig. 2 Experiment set-up for dynamic measurements

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