
Olyps

TLIF
PEEK CAGE

SCIENT'X

YOUR EXPERTISE

OUR SUPPORT

The
th


Enhanced fusion

- Large graft windows to allow an optimum graft area and to improve graft vascularization
- Unique design to enhance anterior and posterior graft impaction

Pure PEEK Optima®

- Modulus of PEEK Optima® is close to the cortical bone's and provides micro-movements to enhance bone growth
- Radiolucent material to ensure accurate fusion follow-up

* From Invibio®

Anatomical shape

- Smooth design to preserve nervous structure during insertion
- Anatomical biconvex shape for a natural adaptation to the concavity of endplates

*easiest,
e safest*


Latest innovation in Scient'x range of PEEK cages, **Olys Transforaminal cage** is dedicated to the lumbar spine. Indications include :

- advanced discopathies
- extensive decompressions
- spondylolisthesis
- failed disc surgery
- recurrent disc herniation
- post-operative instability
- lumbar pseudarthrosis

Olys TLIF cage offers an optimized instrumentation and an innovative design restoring the natural lordosis and stabilizing the operated level.


Respect of the balance

- ➔ Complete range of cages to restore the intervertebral space and the lordosis

Safe approach

- ➔ Unilateral transforaminal approach to minimize dural exposure

Time saving


- ➔ Single cage insertion to reduce operative time

1 Approach and posterior fixation

The following procedure illustrates the insertion of a cage for a Transforaminal Lumbar Interbody Fusion (TLIF) in the L4-L5 intervertebral area.

The patient is in genu pectoral or in prone position. Standard posterior incision and approach are performed. A posterior fixation associated to a postero-lateral fusion is mandatory to ensure medium and long term stability to the segment. Therefore pedicular screws must be inserted before the Olys cage. Afterwards, distraction on these screws will promote exposure of the interbody space and will facilitate the insertion of instruments.


2 Arthrectomie

The lower facet of the cranial vertebra (L4) and the upper facet of the caudal vertebra (L5) are removed on one side with the **osteotome**. This maneuver can be completed by a nerve root decompression by removing the internal part of the upper facet on the caudal vertebra.


Osteotome
220ST02


3 Discectomy

The hemostasis of epidural veins is performed before any disc approach. The protruding nerve root is protected with the **nerve root retractor**. The spinal cord is slightly retracted with the second retractor. Incision of the L4-L5 disc is then performed.

Nerve root retractor
22ECT05-0x

Shavers are used to start the progressive restoration of disc height and facilitate disc removal. Available in 7 sizes (ranging from 6 to 12mm), they are chosen and inserted in accordance with intervertebral height and until resistance to instrument rotation occurs.

Shaver
22FRA01-xx

4 Completion of discectomy

Discectomy is completed using **straight, curved and angled disc rongeurs**. With those three types of disc rongeurs, optimal discectomy is obtained even in the contra-lateral disc space.


Straight disc rongeur
22PIN13

Angled disc rongeur
22PIN15

Curved disc rongeur
22PIN14

5 Distraction


Once discectomy has been achieved, the **lateral distractor** is inserted on the pedicle screws. Increased exposure of the intervertebral space is promoted by activating the toothed bar with the **wrench**.

Wrench
22CLE19


Lateral distractor
22DST17+22EMB02

6 Completion of the approach


The **osteotome tube** is impacted in the intervertebral space to calibrate the window. Four sizes of osteotome tube are available and chosen according to the intervertebral space and the cage size. The rounded extremities must be inserted before any impaction. The tube is gently impacted with a mallet to remove posterior osteophytes and adjust the edges of the window. The window thus obtained will ensure good visibility into the disc space and completion of the discectomy.

Osteotome tube
22OST01-xx

7 Endplates scraping

The remaining layers of annulus and cartilaginous endplates are removed with the **scrapers**. **Curved scrapers** allow to reach all parts of the vertebral endplates in the contra-lateral disc area.

Straight scraper
22CUR05


Right curved scraper
22CUR06D

Left curved scraper
22CUR06G

8 Trial instruments insertion

Inserting **trial cage instruments** allow to check that the discectomy is sufficient for cage insertion while confirming the cage size. To each cage size corresponds a trial cage instrument. At this step, graft can be inserted and impacted in the anterior part of the intervertebral disc space.


Trial cage instruments
22OLSxx-28F

Cage holder
22PRE14


Graft pusher
22COM09

Cage socket
22SOC05



9 Cage preparation

The cage is mounted and locked on the **cage holder**. The extremity of the cage holder is inserted in the cage groove. A 90° rotation of the cage holder with regards to the cage is performed and the ring on the instrument is locked. In order to respect the introduction direction of the cage, the horizontal marker located in the center of the cage must be in the superior part of the cage. Then the cage is inserted in the **cage socket** and filled in with graft (autologous bone, bone substitutes). The graft is impacted in the cage using the **graft pusher** and a mallet.

10 Cage insertion

The cage is first inserted straight into the intervertebral space (A).


Cage holder
22PRE14


The cage holder is then inclined laterally a few degrees (B) to start the cage inclination. The cage holder is unlocked by rotating the ring in order to recover its initial position (C). This maneuver straightens the cage holder without modifying the position of the cage in the intervertebral space. The cage holder is locked once again. The cage insertion in the intervertebral space goes on (D).

*** Warning :** When unlocking the cage holder, and until the cage holder is locked again, any rotation of the cage holder must be avoided. Otherwise, the cage will dissociate from the instrument. The flat parts on the instrument handle allow to control the relative position of the cage holder in relation to the cage.


B


C

D

Once the cage has reached an adequate position, the cage holder is unlocked by rotating the ring. After a 90° rotation of the instrument, the cage holder is released from the cage with a light traction and is removed from the disc space.

Final cage positioning

An X-ray must be performed to verify the proper position of the cage. The straight and curved pushers allow to adjust the position of the cage. The **straight pusher** is positioned at the cage extremity and the **curved pusher** at the center of the cage. Gentle impactions are performed to adjust the position of the cage. Graft can be inserted at the posterior part of the cage.

Straight pusher
22IMP13

Curved pusher
22IMP14

Completion of the posterior fixation

After final positioning of the Olys cage, the optimal stability of the segment will be ensured by adding a posterior fixation and a postero-lateral graft. The cage must be positioned in compression to maintain lordosis and to avoid any future movement of the cage in the intervertebral space.

Olyps

INSTRUCTIONS,

IMPLANTS AND INSTRUMENTS

Olys I N S T R U C T I O N S

Instructions for use

The Olys instrumentation is designed for the surgical treatment of spinal disc diseases. The treatment involves the fusion of two or more vertebrae. This fusion is realized between the vertebral bodies of two adjacent vertebrae. For best results, a detailed preoperative evaluation, a meticulous surgical technique and adequate post-operative care are mandatory. It is important that both the patient and surgeon are fully aware of the risks and possible complications associated with this type of surgery. Before attempting this technique, surgeon are advised to attend a training course with a surgeon already experienced with the use of the device.

DESCRIPTION :

The Olys system is composed of cages inserted into the discal space. The cages are introduced by posterior approach using special instruments. This instrumentation must be associated with a posterior fixation. Fusion is effected at the two vertebral end plates using bone grafts previously introduced into the cages. The implants used in the Olys system are in PEEK as indicated by the 🔗 symbol. It is essential to insert implants with instrumentation specifically designed for this purpose. For more description of the instrumentation it is necessary to read the technical documentation associated to the Olys product. The Olys implants must be assembled with new Olys components defined as being compatible with an other one.

INDICATIONS :

General criteria and principles related to instrumented spinal surgery are applied here :

- Degenerative disc pathologies
- Lumbar pseudarthrosis
- Grade 1 degenerative or isthmic spondylolisthesis (with or without reduction)

CONTRAINDICATIONS :

- Fracture, tumor
- Osteoporosis, calcium metabolism disorder
- Pregnancy
- Infection, recognized allergies to material used

SECONDARY AND POSSIBLE SIDE EFFECTS :

- Pseudarthrosis
- Implant penetration, migration or implant failure
- Infection
- Allergy to materials used

CAUTIONS OF USE :

Never reuse an implant, even in perfect state. Any implant which has been used, twisted, bent, implanted and then removed, even if appears intact, must be discarded. Use new implants routinely. This is indicate on the labelling of the implant package by the next symbol 🔗.

> Preoperatively :

The surgeon must be fully conversant with all aspects of the surgical technique and know the indications and contra-indications of this type of implant. The surgeon must have acquainted himself before the operation with the specific technique for insertion of the product which is available from the manufacturer. As part of the pre-operative examination, the surgeon must check that no biological, biomechanical or other factors will affect the correct conduct of the operation and the post-operative period. An appropriate range of sizes must be available at the time of the operation.

> Perioperatively :

The correct selection of the type of size of implant appropriate to the patient and the positioning of the implant are extremely important

> Postoperatively :

Patients must be informed of the precautions to be taken in their everyday life to guarantee a maximum implant service life. It is recommended that a regular postoperative follow-up is undertaken to detect early signs of failure of the implants and to consider the action to be taken. Deterioration of the device after bone consolidation cannot be considered to constitute a dysfunction or a deterioration in the characteristics of the implant. A suitable rehabilitation program must be designed and implemented.

HANDLING AND STORAGE :

The handling of the Olys material must be done as seldom as possible and always with the utmost care. The Olys implants (in their original packaging) must be stored with care in a clean and dry place. Do not expose the Olys implants to radiations or extreme temperatures. Should these requirements not be followed, reduced mechanical properties may occur which could lead to implant failure in some cases. Proper function of the surgical instruments specific to the Olys instrumentation should be verified prior to every surgical procedure.

CLEANING/DECONTAMINATION/STERILISATION :

The operations of cleaning, decontamination and sterilisation must be realized for implants and instruments before and after use. All packaging must be removed before the next steps :

> Decontamination :

Immerse the implants and the instruments in a bactericidal and fungicidal solution (didecyl ammonium chloride combine with proteolytic enzymes, diluted at 0.5% (5g for 1 liter of mild water). Time of dip : 30'. Rinse out

> Cleaning :

Clean the implants and the instruments in a machine (similar to MIELE machine) with products adapted, rinse out, and dry. All products which can alter the implants are forbidden. Lubricate all hinge and mechanism on instruments with paraffin oil (type Neodisher IP spray)

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We recommend the next sterilization procedure with an autoclave for implants and instruments :

- preheating 25' at 110°C (1 bar)
- vacuum 5' (0,8 bar under atmospheric pressure)
- heating 5' at 120°C (1 bar)
- vacuum 5' (0,8 bar)
- sterilization 18' at 134°C (2 bars)
- drying 20' return to surrounding temperature

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