Vercise™
Neural Navigator 4
Programming Manual

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sl

How to Use This Manual

This manual describes the usage of the Boston Scientific Vercise Neural Navigator Software. Read all instructions carefully before using the DBS Systems.

For indications for use, contraindications, warnings, precautions, adverse events, sterilization, component disposal, storage and handling, and warranty information, refer to the *Information for Prescribers DFU* for your Boston Scientific DBS System as listed in your *DBS Reference Guide*. For other device-specific information not included in this manual or labeling symbols, refer to the appropriate DFU for your Boston Scientific DBS System as listed in your *DBS Reference Guide*.

Guarantees

Boston Scientific Corporation reserves the right to modify, without prior notice, information relating to its products in order to improve their reliability or operating capacity.

Drawings are for illustration purposes only.

Trademarks

All trademarks are the property of their respective holders.

The **Bluetooth®** word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Boston Scientific Neuromodulation Corporation is under license. See also, Declaration ID D035363.

Contacting Boston Scientific

To contact Boston Scientific, see the “*Technical Support*” section of this manual.

Product Model Numbers

Model Number	Description
*DB-7164 and NM-7164	Vercise™ DBS Clinician Programmer
*DB-7164-R and NM-7164-R	Vercise™ DBS Clinician Programmer (Refurbished)
NM-7165	Keyboard
NM-7171	Power Supply for Tablet Computer
DB-7190 and NM-7190	Programming Wand
NM-4512	USB Splitter
NM-6316	International Power Converter
DB-5270	Vercise DBS Remote Control 4
DB-6386	DBS Pairing Magnet
**DB-7105-N4A	Vercise Neural Navigator Software 4.0 Installer

*Applicable after installation of Vercise Neural Navigator 4.0 (Software version 9028429-401).

** Both software part 9028429-401 and 9028602-100 must be installed prior to use.

Table of Contents

Introduction	1
Intended Use	1
Setting Up	2
Placing the Implanted Stimulator into CP Mode (Pairing Mode) (Vercise Genus Stimulators Only)	3
Placing ETS 3 into CP Mode (Pairing Mode)	4
Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)	4
Starting a Programming Session	5
Starting Vercise Neural Navigator	5
Programming the Stimulator	8
Configuring the Leads	8
Measuring Impedances	9
The Programming Screen	10
Creating or Modifying a Program	13
Selecting Stimulation Areas	13
Selecting Contacts	13
Steering Mode	14
Custom Mode	15
Turning Stimulation OFF for Individual Areas	15
Turning All Stimulation OFF	15
Increasing and Decreasing the Amplitude	15
Increasing and Decreasing the Pulse Width	16
Increasing and Decreasing the Rate	16
Programming Multiple Areas with Different Rates	16
Selecting Patient Amplitude Range	16
Viewing the Stimulation Field Model	17
Mapping the Patient's Clinical Effects of Stimulation	18
Ending a Programming Session	19
Magnet Pairing	19
Data	20
Export Database	21
Tools	22
Updates	22
Deleting Clinical Effects Data	23
Deleting Patient Data	23
Additional Information	24
Stimulator Programmable Characteristics	24
Charge Density	25
Energy Use Index	26
Vercise Genus P8, P16, and P32 Non-Rechargeable Stimulators	26
Vercise PC Non-Rechargeable Stimulator	26
Estimated Charge Time	27
Elective Replacement Indicator (ERI) Message	27
End of Service (EOS) Message	27
Clinician Programmer Management	28
Adjusting CP Date and Time	28
Resetting ClinicUser Password	28
Update Password	28
Lost/Forgotten Password	29
Software Installation, Uninstallation, and Removal	29

en

Software Licenses	30
Security.Cryptography.....	30
HIDSharp	30
Newtonsoft.Json.NET	32
SharpDX	32
vNext.....	32
Technical Support	33

Introduction

The Boston Scientific Vercise™ Neural Navigator software is used to program the Vercise™ PC, Vercise Gevia™, and Vercise Genus™ Deep Brain Stimulation (DBS) Systems.

A programming session may include the following activities:

1. Setting Up
2. Starting Vercise Neural Navigator
3. Connecting to the Stimulator
4. Configuring the Stimulator and Leads
5. Testing different Stimulation Settings

This manual will provide instructions on how to accomplish these steps and perform additional functions, such as exporting reports and backing up data.

Boston Scientific DBS Systems utilize MICC¹ technology which is designed to adapt for impedance changes and maintain consistent therapy over time. MICC technology enables current steering across Lead Contacts which is intended to provide precise positioning of stimulation.

If you have any issues, contact Boston Scientific Technical Support.

Note: Screens depicted in this manual may differ slightly from the screens on your Vercise Neural Navigator.

Intended Use

Vercise Neural Navigator is a software program that is used to set and adjust stimulation parameters for the Vercise PC, Vercise Gevia, and Vercise Genus DBS Systems.

¹ Multiple independent current control

Setting Up

The Clinician Programmer (CP) communicates with the Stimulator via wireless telemetry. Wireless Programming allows patients to move about, if instructed by the clinician, while the clinician adjusts parameters. Vercise Genus devices utilize Bluetooth technology for direct wireless communication between the CP and the Stimulator. Vercise PC and Vercise Gevia devices communicate via a Programming Wand. The Programming Wand uses an inductive telemetry radiofrequency (RF) link to communicate with the Stimulator.

Caution: Use only Vercise PC, Vercise Gevia, or Vercise Genus DBS System components with Vercise Neural Navigator. Failure to do so may result in the inability to program the Stimulator.

Caution: The CP is not considered as equipment for the patient environment as defined by IEC 60601-1. The CP and the person using the CP should not be in physical contact with the patient while programming.

1. Plug the CP into a power source.
2. Power ON the CP.
3. Login as ClinicUser. If this is the first login, enter password “bsn”. You will be prompted to setup a new ClinicUser password when logging into the CP for the first time. The new ClinicUser password must be 10 or more characters in length.

Note: In regions where software is installed by a third party, follow region-specific guidance for the ClinicUser password.

4. Determine the type of Stimulator to be programmed and follow the appropriate sections as listed in Table 1, then continue on to the “Starting a Programming Session” section of this manual.

Stimulator Type	Model Numbers	Reference Section
Vercise Genus IPG(s)	DB-1408, DB-1416, DB-1432, DB-1216, DB-1232	See “Placing the Implanted Stimulator into CP Mode (Pairing Mode) (Vercise Genus Stimulators Only)” on page 3 of this manual.
ETS 3	DB-5170	See “Placing ETS 3 into CP Mode (Pairing Mode)” on page 4 of this manual.
Vercise Gevia IPG	DB-1200	See “Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)” on page 4 of this manual.
Vercise PC IPG	DB-1140	See “Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)” on page 4 of this manual.
ETS 2	DB-5132	See “Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)” on page 4 of this manual.

Placing the Implanted Stimulator into CP Mode (Pairing Mode) (Vercise Genus Stimulators Only)

- Use the patient's Remote Control to enter the Stimulator into CP Mode (Pairing Mode). The Remote Control and Stimulator must be linked in order to place the Stimulator into CP Mode (Pairing Mode).
 - From the **Lock** screen, press and hold the **P** button until the screen displays "Stimulator Ready for CP Connection" (Figure 1 and Figure 2).


Figure 1. Programs Button on Remote Control


Figure 2. Stimulator Ready for CP Connection

- Alternatively, navigate to the **+** Clinician Menu on the patient Remote Control and select **+** CP Mode: Stimulator (Figure 3).

Figure 3. CP Mode: Stimulator

Note: If the Stimulator does not have a linked Remote Control, refer to the appropriate Remote Control manual as listed in your DBS Reference Guide for instructions on Stimulator Linking.

Note: If the linked Remote Control is not available, see the "Magnet Pairing" section of this manual.

Note: The Stimulator will automatically exit CP Mode (Pairing Mode) after two minutes if a CP connection is not established.

Placing ETS 3 into CP Mode (Pairing Mode)

1. Press and hold the  **Stimulation Button** on the side of the ETS until the lights flash (Figure 4). The ETS battery indicator light will alternate blinking green and yellow, indicating that the ETS is in CP Mode (Pairing Mode) and is available for connection.

Note: *The Stimulator will automatically exit CP Mode (Pairing Mode) after two minutes if a CP connection is not established.*

en

Figure 4. ETS Battery Indicator Light

Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)

1. Connect the Programming Wand to the CP using the USB cable provided with the Programming Wand (Figure 5).
 - (a). Plug the Mini USB end of the USB Cable into the USB Port on the side of the Programming Wand.
 - (b). Plug the Standard USB end of the USB Cable into the USB Port on the CP.



Figure 5. Vercise PC and Vercise Gevia: Clinician Programmer (CP) and Programming Wand

2. Wait for the Wand to perform a self-test. At the end of the self-test, the Wand will beep.
3. If the Power Light on the Wand is green, place the Wand over the Stimulator.
 - (a). If the Power Light on the Wand remains red, contact Technical Support.

Starting a Programming Session

Starting Vercise Neural Navigator

1. Select the Vercise Launcher icon  on the desktop.
2. Select  to start Vercise Neural Navigator.

Note: If Brainlab Elements is present on the CP, Vercise Neural Navigator may be launched from within Elements .

Note: Multiple software should not be run simultaneously on the same CP (except when launching Vercise Neural Navigator from within Elements).

Note: Vercise Neural Navigator can also be launched in Demo Mode using the Vercise Launcher. Demo Mode is used for demonstration purposes only (Figure 6).

Figure 6. Launch Screen with DEMO Option Mode

en

3. Upon starting Vercise Neural Navigator, the screen will show the **Connect Tab** and the software will automatically search for a Stimulator to connect to (Figure 7 and Figure 8).

(a). **Vercise Genus:** A Stimulator that is in CP Mode (Pairing Mode) and is within range will automatically appear on the **Connect Tab** (Figure 7). If no Stimulator is found and displayed on the **Connect Tab**, confirm that the Stimulator is in range and is in CP Mode (Pairing Mode) by viewing the Patient Remote Control Screen. Select the **Refresh** button on the CP (Figure 7).

Note: A Vercise Genus Stimulator must be in CP Mode (Pairing Mode) in order to establish a connection with the CP.


Figure 7. Vercise Genus Connect Tab

(b). **Vercise PC and Vercise Gevia:** The **Connect Tab** will show the Programming Wand (Figure 8). Press the **Scan** button. The CP will then scan for devices using the Programming Wand. If only one Stimulator is found in range, the CP will automatically connect to that Stimulator. If no Stimulator is found, move the Wand closer to the Stimulator that you are trying to connect to and select the **Scan** button.

Note: Press the Refresh button if the Programming Wand does not appear on the screen.


Figure 8. Vercise PC and Vercise Gevia Connect Tab

Note: The CP cannot connect to a Stimulator when the Stimulator is in MRI Mode. Exit MRI Mode using the Remote Control and Refresh to connect. For instructions on exiting MRI Mode, refer to the Remote Control Manual as listed in your DBS Reference Guide.

Table 2: Vercise Neural Navigator Tab Descriptions		
#	Feature	Description
1	Connect Tab	Displays the connection status between the CP and Stimulator.
2	Configure Tab	Configure Leads and edit patient profile. Note: If Brainlab Elements is present on the CP, some patient data may be imported from Elements.
3	Program Tab	Adjust the Stimulator program settings.
4	Data Tab	Generate, print and export reports, and export or delete selected patient(s) data.
5	Tools Tab	Manage data and updates.
6	End Session Tab	Disconnect from the Stimulator or exit the application.
7	Battery Indicator	Displays the battery status of the CP.

4. Press the **Connect** button next to the Stimulator.
5. Once a connection has been established between the CP and the Stimulator, the following screen will appear (Figure 9 and Figure 10).


Figure 9. Connection Established Between CP and Stimulator (Vercise Genus or ETS 3)





Figure 10. Connection Established Between CP and Stimulator (Vercise Gevia, Vercise PC, or ETS 2)

Table 3: Description of the Connect Tab														
#	Feature	Description												
1	Patient ID	Displays the Patient ID number.												
2	Connection Status	Displays the connection status between the CP and Stimulator along with the Model and Serial Number of each of the devices.												
		<table border="1"> <thead> <tr> <th>#</th> <th>Connection Status</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td>1</td> <td></td> <td>Connection between CP and Vercise Genus or ETS 3 Stimulator was established.</td> </tr> <tr> <td>2</td> <td></td> <td>Connection between CP, Programming Wand, and Vercise Gevia, Vercise PC, or ETS 2 Stimulator was established.</td> </tr> <tr> <td>3</td> <td></td> <td>Connection not established. Ensure the Vercise Genus or ETS 3 Stimulator is in CP Mode (Pairing Mode). For Vercise Gevia or Vercise PC, move the Programming Wand closer to the Stimulator and Scan again, and/or see the “Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)” section of this manual or refer to the Vercise DBS Programming Wand Manual for more information.</td> </tr> </tbody> </table>	#	Connection Status	Description	1		Connection between CP and Vercise Genus or ETS 3 Stimulator was established.	2		Connection between CP, Programming Wand, and Vercise Gevia, Vercise PC, or ETS 2 Stimulator was established.	3		Connection not established. Ensure the Vercise Genus or ETS 3 Stimulator is in CP Mode (Pairing Mode). For Vercise Gevia or Vercise PC, move the Programming Wand closer to the Stimulator and Scan again, and/or see the “Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)” section of this manual or refer to the Vercise DBS Programming Wand Manual for more information.
		#	Connection Status	Description										
1		Connection between CP and Vercise Genus or ETS 3 Stimulator was established.												
2		Connection between CP, Programming Wand, and Vercise Gevia, Vercise PC, or ETS 2 Stimulator was established.												
3		Connection not established. Ensure the Vercise Genus or ETS 3 Stimulator is in CP Mode (Pairing Mode). For Vercise Gevia or Vercise PC, move the Programming Wand closer to the Stimulator and Scan again, and/or see the “Connecting the Programming Wand to the CP (Vercise Gevia, Vercise PC, and ETS 2 Stimulators Only)” section of this manual or refer to the Vercise DBS Programming Wand Manual for more information.												
3	Refresh Button	Scan for available Stimulators or Programming Wand. Disabled if the CP is already connected to the Stimulator.												
4	Connection Status	Displays a green or blue filled circle if the Stimulator is connected to the CP.												
5	Connect or Disconnect Button	Connect or disconnect from a Stimulator. When a Stimulator is not connected, this button reads “Connect.” When a Stimulator is connected, this button reads “Disconnect.”												
6	Signal Strength	Displays the signal strength for communication between the CP and Stimulator (Vercise Genus or ETS 3 only).												

Programming the Stimulator

Configuring the Leads

en

Once a connection has been established between the CP and Stimulator, switch to the **Configure Tab** to configure the Leads that are connected to the Stimulator (Figure 11 and Figure 12). During an initial programming session, Lead Configuration must be completed prior to navigating to the **Program Tab**. Before proceeding to the **Program Tab**, verify the port assignment within the **Configure Tab** accurately reflects the physical connection of the Lead to Stimulator port(s).

Note: Once a Stimulator has been initially configured, you can directly switch to the Program Tab after connection has been established from the Connect Tab.

Note: If Brainlab Elements is available on the CP, you may import Lead information and objects into the Neural Navigator software from Elements. To import or remove Elements data, Link or Unlink from the Configure screen. This option will only be available if Neural Navigator is launched from within Elements.


Figure 11. Configuring a 16 Contact Stimulator

During the implantation procedure of a 32 Contact (4 Port) Stimulator, no more than two Leads should be connected to the Stimulator. Therefore, only two of the four ports are to be configured for programming. Only Ports L2 and R2 of the 32 Contact IPG are available for programming.


Figure 12. Configuring a 32 Contact Stimulator

Table 4: Description of the Configure Tab		
#	Feature	Description
1	Leads Configuration	For each Lead, select the Lead type, Stimulator Port to which the Lead is connected, and brain hemisphere. Enter the target Area. For Directional Leads, select the Directional Marker orientation. Note: When Lead information has been imported from Elements, changes made in the Configure Tab rather than in Elements may cause the Stimulator to be unlinked from Elements.
2	Measure Button	Measure Impedances. See the “Measuring Impedances” section of this manual for more information.
3	Link/ Relink/ Unlink Patient	Import or remove data from Elements. Note: This is only available if Brainlab Elements is installed on the CP and Neural Navigator is launched from within Elements.
4	Patient Amplitude Control	Turn ON/OFF the patient’s ability to change stimulation Amplitude. The range of Patient Amplitude Control is adjusted in the Program Tab .
5	Stimulator Information	Displays the Stimulator Information including Serial number, Model number, Firmware Version and type of Stimulator.
6	Implant Date	Displays the date on which a CP first connects to a new Stimulator. The Implant Date can be adjusted by selecting the Implant Date button.
7	Patient ID	The Patient ID automatically defaults to the Stimulator serial number. The Patient ID can be edited by typing into the Patient ID field.

Measuring Impedances

Impedances can be measured using the **Measure** button on the **Configure** or **Program Tabs**. The impedances of each Contact may be used to verify electrical integrity. When an impedance measurement is taken, impedances are assessed between a Contact and the Stimulator Case (monopolar), and between pairs of Contacts (bipolar). Impedances over 8000 Ω may be the result of open or unconnected wires and are displayed in yellow on the **Impedance Measurement** window. Impedances less than 200 Ω may be the result of short circuits and are displayed in orange. Contacts that have impedances outside the acceptable range are marked with a Ω symbol on the programming screen. The most recent set of impedance measurements are included in a report that can be printed or exported from the **Data Tab**.

The Programming Screen

Once the Leads have been configured, select the **Program** Tab to begin programming. The Programming screen is divided into the following sections and features as shown in Figure 13. Programming features specific to the Directional Lead and Programming a Directional System are shown in Figure 14. The Clinical Effects Map is shown in Figure 13 and can be viewed for both a Standard Lead and Directional Lead. STIMVIEW™, or the Stimulation Field Model (SFM), is shown in Figure 14 and can be viewed for both a Standard Lead and a Directional Lead.

en

Figure 13. Programming Screen


Figure 14. Directional Lead Programming Screen

#	Feature	Description				
1	Program Button	Select the Program that you would like to set up or adjust.				
2	Program Options Button	<ul style="list-style-type: none"> View Energy Use Index for non-rechargeable Stimulators. View battery recharge estimate for rechargeable Stimulators. Delete and copy Programs. Change ramp and cycle times for Programs. 				
3	Area Panel	Select the Area within a Program that you would like to set up or adjust.				
4	Area Options Button	Delete an Area within a Program or import simulations from Brainlab Elements (if Elements is present on the CP).				
5	+	Add an Area. Select from one of the Stimulator Ports defined in the Configure Tab .				
6	Stimulation ON/OFF Button	Turn stimulation OFF for the Area selected. Note: When the Amplitude is at 0 mA, increase the Amplitude to turn ON stimulation.				
7	Pulse Width Button	Adjust the Pulse Width. <table border="1"> <thead> <tr> <th>Default</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>60 μs</td> <td>20 μs to 450 μs</td> </tr> </tbody> </table>	Default	Range	60 μ s	20 μ s to 450 μ s
Default	Range					
60 μ s	20 μ s to 450 μ s					
8	Rate Button	Adjust the Rate. <table border="1"> <thead> <tr> <th>Default</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>130 Hz</td> <td>2 Hz to 255 Hz</td> </tr> </tbody> </table>	Default	Range	130 Hz	2 Hz to 255 Hz
Default	Range					
130 Hz	2 Hz to 255 Hz					
9	Units Button	Select the Units in which Amplitude is displayed on Contacts and Stimulator Case. <table border="1"> <thead> <tr> <th>Default</th> <th>Alternative</th> </tr> </thead> <tbody> <tr> <td>%</td> <td>mA</td> </tr> </tbody> </table>	Default	Alternative	%	mA
Default	Alternative					
%	mA					
10	Patient Amplitude Buttons	Adjust Maximum and Minimum Patient Amplitude. Note: The Patient Amplitude buttons are only displayed if the Patient Amplitude Control has been set to ON from the Configure Tab.				
11	Stop All Button	Turns off all stimulation.				
12	Contact and Stimulator Case Configuration	Displays percentage of anodic (+) or cathodic (-) energy assigned to the Lead Contacts and Stimulator Case for a given Area. See the "Selecting Contacts" section of this manual for more information.				
13	Measure Button	Measures impedances of the Contacts.				
14	Programming Modes	Select Steering or Custom Programming modes.				
15	Step Size	Select the step size for Amplitude adjustments: 0.1 mA or 0.5 mA.				
16	Level Up and Down Buttons	Steer stimulation focus along the Lead. Press the Level number to select a new Level.				
17	Total Amplitude	Increase or decrease the total Amplitude delivered for a given Area.				
18	Clinical Effects Panel	Make note of the Therapeutic Benefits and/or Side Effects for the current Stimulation settings.				
19	Text Notes	Capture text notes for a given Lead (up to 250 characters per Area).				
20	Clinical Effects Map	Graphical summary of assigned Therapeutic Benefits and/or Side Effects at a given position along the DBS Lead array and a stimulation Amplitude. Note: Clinical effects data is captured and listed in reports but not plotted on the Clinical Effects Map for configurations not possible in Steering Mode and for Directional Lead settings that are not 100% focused or spread.				
21	Display Drop-down	View control to switch between the Clinical Effects Map, 3D Overview, or 3D Split View of the stimulation field.				
22	Reference Head	The reference head demonstrates the relationship of the Lead currently being programmed to the position of the patient head.				
23	Clinical Effects Legend	The level of therapeutic benefit is indicated by the saturation of the dot.				

Table 5: Description of the Program Tab		
#	Feature	Description
24	STIMVIEW™ or Stimulation Field Model (SFM)	Visual representation of the estimated stimulation field for the currently programmed stimulation parameters. Note: <i>If Brainlab Elements is available on the CP, imported objects may be viewed in STIMVIEW, referred to as STIMVIEW XT when imported objects are displayed.</i>
25	Directional Presets	Select the one touch buttons to adjust the stimulation field. The directional presets will steer the fully focused stimulation field in one of four orthogonal directions or put the stimulation field into “ring mode.” Ring Mode generates, from a segmented Contact Level, stimulation fields equivalent to those generated by a standard “ring” or cylindrical Contact. Note: <i>Only applies to Directional Leads.</i>
26	Rotate Buttons	Steer the stimulation focus circumferentially around the Lead. Note: <i>Only applies to Directional Leads.</i>
27	Spread/Focus Buttons	Radially spread or focus the stimulation field. Note: <i>Only applies to Directional Leads.</i>
28	Directional Indicator	Visual indicator of the orientation of the radiopaque Directional Marker band on the Directional Lead. The orange line and dot correlate to the center of the radiopaque Directional Marker. Note: <i>Only applies to Directional Leads.</i>
29	STIMVIEW™ View Controls	Adjust the view of the SFM using Zoom, Rotate, Panning Control, or Reset to the original view. Both the Lateral and Axial views of the SFM will adjust in unison using these controls when in 3D Split View, but must be adjusted in the Lateral view.
30	Virtual Contact	Dotted ring illustrating the axial location of stimulation along the Lead. The arrow indicator illustrates the rotational orientation of stimulation around the Lead. The dotted ring and arrow indicator together form the Virtual Contact.

Creating or Modifying a Program

To create a new Program or modify an existing Program, select the **Program** button and choose one of the four Programs from the drop-down arrow. The system allows you to configure up to four Programs on a Stimulator.

For a given Program, you can view and/or adjust several options via the **Program Options** button. The Program Options include the following:

#	Feature	Description						
1	Battery	For a non-rechargeable Stimulator, the Energy Use Index for the current Program is displayed. This value is used to provide an estimate of battery longevity for the non-rechargeable Stimulator on the current Program. See the “ <i>Energy Use Index</i> ” section of this manual for more information. For a rechargeable Stimulator, an Estimated Charge Time for the current Program is displayed. This value provides an estimate of the duration and frequency of charging necessary to maintain stimulation.						
2	Ramp	Time to gradually increase the stimulation from zero to the programmed Amplitude when stimulation is turned ON. <table border="1"> <thead> <tr> <th>Default</th> <th>Options</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>ON</td> <td>ON/OFF</td> <td>1 sec to 10 sec</td> </tr> </tbody> </table>	Default	Options	Range	ON	ON/OFF	1 sec to 10 sec
Default	Options	Range						
ON	ON/OFF	1 sec to 10 sec						
3	Cycle	The cycled on and off duration of stimulation delivery. <table border="1"> <thead> <tr> <th>Default</th> <th>Options</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>OFF</td> <td>ON/OFF</td> <td>1 sec to 90 min</td> </tr> </tbody> </table>	Default	Options	Range	OFF	ON/OFF	1 sec to 90 min
Default	Options	Range						
OFF	ON/OFF	1 sec to 90 min						
4	Copy to	Copy the current Program settings to another Program.						
5	Delete	Delete the settings for the current Program.						

Selecting Stimulation Areas

For a given Program, you may configure up to four Areas. With a new Program, an Area will be automatically assigned to each brain hemisphere and named based on the defined Target and side of the brain selected in the **Configure Tab**. You may add an additional Area by selecting an empty Area (+) and choosing a configuration (e.g. Left STN). You may reassign an Area by first selecting **Options** for that Area and then choosing **Delete Area**. You may then select a different hemisphere configuration.

Note: *Predefined simulations may be imported from Elements by selecting Load Simulation from the Area Options button.*

Selecting Contacts

You can manually assign anodes and cathodes in **Custom Mode** or incrementally steer a stimulation field along the Lead in the **Steering Mode**. Steering Mode is limited to a monopolar configuration of either a single cathode or adjacent cathodes.

To provide anodic stimulation, you must use **Custom Mode**. You may assign the Stimulator Case and all the Contacts as anode or cathode individually in Custom Mode. The External Trial Stimulator (ETS) is limited to **Custom Mode** since the Stimulator Case cannot be assigned as a cathode or anode.

Note: *Switching from Custom Mode to Steering Mode will clear the Contact and Stimulator Case assignments.*

Note: *Stimulation using multiple independent current control (MICC) with a Directional Lead is referred to as Cartesia 3D.*

Note: *Stimulation with only Lead Contacts assigned as anode(s) and the IPG case assigned as cathode is referred to as Monopolar Anodic Stimulation (MAS).*

en

Steering Mode

Steering Mode is a simplified programming mode where the Contact(s) act as the cathode(s) and the Stimulator Case acts as the anode. This mode allows you to steer a monopolar cathode along the Lead, eliminating the need to turn ON and OFF individual Contacts. Steering Mode incrementally shifts a percentage of the cathodic current to the adjacent Contact(s) using current steering technology to create smooth transitions between Contacts.

The DB-2201 Standard Lead has eight ring Contacts per Lead, labeled 1 through 8 on each Lead.

To Steer along the DB-2201 Standard Lead:

1. Select **Steering Mode**.
2. Select a Contact to assign it as a 100% cathode.
3. Use the ▲ and ▼ buttons to incrementally steer the stimulation focus along the length of the Lead. The amount of cathodic current will shift in 10% increments.

Note: You may also adjust directly from one Contact (Level) or half level to another. The Amplitude for the selected Area will drop to 0 mA when another Contact is selected, but not when steering in 10% increments.

The DB-2202 Directional Lead has a total of eight Contacts per Lead, labeled 1 through 8 on each Lead. Contacts 1 and 8 are the distal tip and proximal ring Contacts, respectively, while Contacts 2 through 7 are the small directional Contacts (segments) grouped into two rows of three segments each.

To Steer along the DB-2202 Directional Lead:

1. Select **Steering Mode**.
2. Select a Contact to assign it as a cathode. You may create an equal spread of current across a Level of Contacts (“ring mode”) by selecting anywhere within that Level, then selecting the center button. To assign a single directional segment as the cathode, select anywhere within that Level, then select the corresponding button (Figure 15).


Figure 15. Directional Contact Selector

3. Use the ▲ and ▼ buttons to incrementally steer the stimulation focus along the length of the Lead.

Note: You may also adjust directly from one Contact (Level) or half level to another. The Amplitude for the selected Area will drop to 0 mA when another Contact is selected, but not when steering in 10% increments.

4. Select one of five preset directions (○, ↓, ←, ↑, →) for the stimulation field. The directional presets will steer the fully focused stimulation field in one of four orthogonal directions or put the stimulation field into “ring mode.”

The following steps can be used to refine the applied directional preset or selected directional segment.

5. Use the  and  buttons to rotate and steer the stimulation focus circumferentially around the Lead. Each rotation is 30 degrees.
6. Use the  and  buttons to radially spread or shrink the focus of the stimulation field.
7. To choose another starting point or to steer on another Contact, select another Contact. To select a segmented Directional Contact, select the Level, then select one of the three labeled segmented Contacts around the circumference of the center button on the Directional Contact Selector.

Note: The total Amplitude for the selected Area will drop to 0 mA when another Contact is selected.

Note: Stimulation using multiple independent current control and the Directional Lead is referred to as Cartesia 3D.

Custom Mode

Custom Mode allows you to assign a percentage of anodic or cathodic current to individual Contacts and the Stimulator Case.

To program in Custom Mode:

1. Select **Custom Mode**.
2. Select the Stimulator Case or Contact that you want to adjust. If it was selected, one tap will assign it as an anode (+). Another tap will reassign it as a cathode (-). Another tap will reassign it as OFF (blank). Tapping on a Contact will first select it without changing the polarity.

Note: Changing the Contact polarities will reset the Amplitude to zero.

3. Select the + and – buttons for the Contact to adjust the percentage of anodic or cathodic current assigned to the selected Contact.

Note: When using the External Trial Stimulator (ETS), monopolar configurations are not possible since the ETS “case” cannot be assigned as a cathode or anode.

Note: When using the ETS, Clinical Effects data is recorded but not plotted on the CEM.

Turning Stimulation OFF for Individual Areas

To turn stimulation OFF for Individual Areas:

1. Make sure the Area you wish to turn OFF is selected by clicking on the appropriate Area on the Area Panel.
2. Press the Stimulation OFF button to turn Stimulation OFF.

Note: When the Amplitude is at 0 mA, increase the Amplitude to turn ON stimulation.

Turning All Stimulation OFF

Selecting the  button will stop stimulation for all active Areas. This function is only meant for turning all Stimulation OFF. To turn Stimulation ON, select each Area that you want to turn ON and select the Stimulation ON/OFF switch.

Increasing and Decreasing the Amplitude

Amplitude is measured in milliamperes (mA). The default setting for Amplitude is 0 mA and the range is 0 to 20 mA. The maximum Amplitude for a single Contact is 12.7 mA.

To increase or decrease the Amplitude:

1. Use the + and – buttons labeled Total Amplitude to increase or decrease Amplitude.
2. The default step size for Amplitude changes is 0.1 mA. You can change the step size to 0.5 mA using the **Step Size** buttons.

Note: High stimulation levels can cause permanent tissue damage. A message will pop-up notifying you if you attempt to exceed a stimulation limit and settings that exceed this limit are not allowed.

Increasing and Decreasing the Pulse Width

The Pulse Width of the stimulation is the length of time a burst of energy is applied per pulse. The Pulse Width is measured in microseconds (μs). The default setting for Pulse Width is 60 μs and the range is 20 μs to 450 μs .

To increase or decrease the Pulse Width:

1. Select the **Pulse Width** button.
2. Select the desired Pulse Width from the options provided. Pulse Widths that exceed stimulation limits are greyed out.

Note: High stimulation levels can cause permanent tissue damage.

Note: Increasing the Pulse Width by more than 10 microseconds (μs) at a time will cause the total Amplitude to be reset to 0 mA.

Increasing and Decreasing the Rate

The pulse rate of the stimulation, often called the Rate or the Frequency, dictates how many stimulation pulses are delivered in a second, measured in Hertz (Hz) or pulses per second (pps). The default setting for Rate is 130 Hz and the range is 2 to 255 Hz.

To increase or decrease the Rate:

1. Select the **Rate** button.
2. Select the desired Rate from the table of available Rates. Incompatible Rates are greyed out.

Note: Areas that are assigned to the same Lead Port cannot have Rates that sum up to greater than 255 Hz.

Programming Multiple Areas with Different Rates

Areas are allowed to be programmed with different Rates. By default, the Multiple Rate option is disabled. When you enable Multiple Rates, only the Rates that are compatible with the Rates and Pulse Widths from the other active Areas are available.

Note: Modifying the Rate of an Area will alter the available Rates for the other Areas.

Note: If you disable Multiple Rates, the Rate for all Areas will be reset to the Rate selected for the current Area.

Selecting Patient Amplitude Range

By default, patients do not have the ability to adjust the Amplitude of their stimulation.

However, in some cases, you may want to give a patient the ability to adjust the Amplitude of their stimulation using the Remote Control. To give patients Amplitude Control, first turn ON the Patient Amplitude Control in the **Configure Tab**. Once the Patient Amplitude Control is turned ON, you can set the allowable Amplitude range in the **Program Tab** for each Area by setting a Minimum and Maximum.

Note: High stimulation levels can cause permanent tissue damage. A message will pop-up notifying you if you attempt to exceed a stimulation limit and settings that exceed this limit are not allowed.

Viewing the Stimulation Field Model

The Stimulation Field Model (SFM), called STIMVIEW, is a visual representation of the estimated stimulation field for the currently programmed stimulation parameters. The SFM includes both a visual representation of the DBS Lead as well as the approximated stimulation field shown in the color red (Figure 14 and Figure 16). As programming parameters are adjusted and the stimulation is steered along the Lead, the SFM will adjust accordingly.

You may switch between two different views by selecting either 3D Overview or 3D Split View from the Display drop-down (Figure 16). 3D Overview presents a three-dimensional view within which you can zoom, rotate, and pan. The 3D Split View provides a dual pane view centered on the Lead. The top pane is in-line with the Lead and the bottom pane is on an axis perpendicular to the Lead. If objects have been imported from Elements, you may show or hide individual objects using the Object Visibility drop-down.

Adjust the view of the SFM using  to Zoom,  to Rotate,  to Pan, or  to Reset the original view. When in 3D Split View, both the Lateral and Axial views of the SFM will adjust in unison using these controls. These controls will not affect or adjust any programming parameters.


Figure 16. Display Drop-down

Mapping the Patient's Clinical Effects of Stimulation


Figure 17. Clinical Effects Map

For a given stimulation setting, you may make note of a 0 to 4 rating for each therapeutic benefit and a 0 to 4 rating for each side effect by selecting the button labeled with the symptom or side effect, then selecting the appropriate numerical rating. If selection of a numerical rating is not desired, select anywhere outside of the Therapeutic Benefit and/or Side Effect rating box to close. To remove your Therapeutic Benefit and/or Side Effect selection, select the Therapeutic Benefit or Side Effect that you want to remove, and then select the Therapeutic Benefit or Side Effect from within the pop-up to remove the highlighted section. Each button that is selected is captured as data associated with that stimulation setting for that patient.

You may also select the **Text Notes** button to enter and save up to 250 characters of text associated with each Lead.

When capturing Clinical Effects in Steering Mode, a dot is plotted on the CEM at the Axial Lead position and Amplitude. When programming directionally, the CEM switches to a polar grid. A new CEM will be created for programming settings at different Levels (axial positions along the Lead). The rating scale of the therapeutic benefit determines the color saturation of the center of dot. A visual key indicating color saturation for a score appears at the bottom of the CEM when programming in ring mode (100% spread). If a side effect is selected, an orange ring is displayed. Selecting a dot will display a pop-up window containing the date and time at which the dot was captured along with the stimulation setting and effects details (Figure 17).

All of this data is saved on the Stimulator and available for export in the **Reports Tab**.

Note: Clinical effects data is captured and listed in reports but not plotted on the CEM for configurations that are not possible in Steering Mode.

Note: The reference head in the CEM view highlights the hemisphere in which you are currently programming.

Note: A Clinical Effects Map is only displayed at 100% focus or 100% spread (ring mode).

Ending a Programming Session

To end a Programming Session on the CP:

1. Select the **End Session Tab**.
2. Select **Exit Application** to end the programming session and close the application.
3. Alternatively, select **Disconnect from Stimulator** to end the programming session and disconnect from the patient's Stimulator. This will take you back to the **Connect Tab**.

All Programs and programming data are automatically saved in real-time during the programming session. No step to actively "save" is required. The Patient's Remote Control automatically syncs with the Stimulator to which it has been linked.

Caution: *Do not place a magnet over a Stimulator within 60 seconds of ending a programming session. If you need to reconnect to a Stimulator, either utilize the Remote Control to initiate CP Mode (Pairing Mode), or wait 60 seconds before utilizing a magnet to place the IPG into CP Mode (Pairing Mode).*

Magnet Pairing

If a Remote Control is not available to place the implanted Stimulator into CP Mode (Pairing Mode), you may utilize a magnet. Place the DBS Pairing Magnet directly over the Stimulator for at least 2 seconds and then remove it.

Note: *Once the DBS Pairing Magnet is removed, the Stimulator will automatically complete a device reset cycle. This device reset cycle will take several seconds to complete. During this reset cycle, the Stimulator will not provide Stimulation and cannot communicate. Once the reset is complete, the Stimulator will return to normal operation.*

For instructions on using the Remote Control to place the implanted Stimulator into CP Mode (Pairing Mode), see the "Placing the Implanted Stimulator into CP Mode (Pairing Mode) (Vercise Genus Stimulators Only)" section of this manual.

Data

Within the **Data Tab**, you may generate reports for the current programming session, or for patients that have been previously programmed using the same Clinician Programmer.

To generate a report for the current programming session, select the **Data Tab** (Figure 18). A Report may be printed and exported as a PDF or Excel file.


Figure 18. Data Tab

Select  on the **Data Tab** and select the desired information you want to include in the report by checking on any of the following check boxes:

- Programs
- Configuration
- Clinical Effects Maps
- Clinical Effects Details
- Anonymize Patient Data

You can also view reports for all Stimulators that were connected to the CP. Reports can be viewed when the CP is not connected to a Stimulator.

To view reports when the CP is not connected to a Stimulator (Figure 19):

1. Select the **Data Tab**.
2. Select the patient whose report you would like to view and select **View**.


Figure 19. Viewing Reports When CP Is Not Connected to Stimulator

Export Database

The Export feature allows you to backup a single patient's data or the entire patient database on the CP to a specified location. The backup location can be a folder on the CP or an external storage drive (for example, USB flash drive). This feature can be accessed from the **Data Tab**.


Figure 20. Backup Database

To create a backup of a single patient's data or the entire patient database (Figure 20):

Note: In order to export the data of multiple patients or another patient, you must disconnect from all Stimulators.

1. Select the **Data Tab**.
2. Select the patient record that you would like to export. To select multiple records, you may press **CTRL** (control) on the keyboard while selecting desired records or select **Select All**.
3. Select **Export**.
4. If desired, select **Anonymize Patient Data**.
5. Select **Browse** to choose a backup location.
6. Select **Backup** to perform the backup function.

Note: After the backup is completed, a pop-up window will confirm the location of the file and indicate if the backup was successful.

Tools

The **Tools Tab** allows you to activate Leads, delete clinical effects data, and delete patient data.

en

Updates

The **Updates Tab** allows you to enable supported Leads and features.

Note: You can only perform updates when disconnected from the Stimulator.

Note: Only features available in your region will be provided.


Figure 21. Enable a New Lead or Feature

To enable a new Lead or feature (Figure 21):

1. Disconnect from the patient's Stimulator by clicking on the **Disconnect** button on the **Connect Tab**.
2. Go to the **Tools Tab**.
3. Select the **Updates Tab**.
4. Select the row for the desired feature.
5. Select the column labeled **Key** and enter the key provided by your local representative.
6. Select **Verify Key(s)**.

Deleting Clinical Effects Data

All the Clinical Effects Data for a patient can be deleted from the **Tools Tab** under the **Clinical Effects Data Tab**.

Note: This feature is available only when the CP is linked to a Patient's Stimulator.


Figure 22. Delete All Clinical Effects Data

To delete the Clinical Effects Data (Figure 22):

1. Go to the **Tools Tab**.
2. Select the **Clinical Effects Data Tab** and select **Delete**.
3. Select **Continue**.

Deleting Patient Data

Patient data can be deleted after data deletion is enabled on the **Data Management Tab**.

To delete Patient Data for one or all patients:

1. Go to the **Tools Tab**.
2. Select the **Data Management Tab**.
3. Select the **Allow Data Deletion** checkbox (Figure 23) and select **Close**.
4. Navigate to the **Data Tab**. The **Delete** button is now displayed on the right-side panel.
5. Select the patient(s) data you would like to delete.
6. Select **Delete**.


Figure 23. Allow Data Deletion Checkbox

Additional Information

Stimulator Programmable Characteristics

The stimulation parameters are independent for each DBS Lead such that stimulation of different brain targets can have different Amplitudes, Pulse Widths, Stimulation Rates, and Contact configurations. It is possible to configure one Lead as monopolar, and one as multipolar. It is also possible to configure a single Lead with both monopolar and multipolar Areas.

The programmable parameter ranges for the Stimulator are shown below.

#	Parameter	Parameter Range
1	Waveform	Charge balanced, asymmetric biphasic
2	Pulse Shape	Rectangular
3	Current or Voltage Regulated	Current
4	Amplitude ²	0.1 mA to 20 mA
5	Rate ³	2 Hz to 255 Hz
6	Pulse Width ⁴	20 μ s to 450 μ s
7	Cycle On/Off	1 sec to 90 minutes
8	Ramp On	1 sec to 10 seconds
9	Contact Connections	16
10	Independent Areas of Stim (4 Programs with 4 Areas per Program)	16
11	Current Path Options	Unipolar, Bipolar, Multipolar

2 The programmable coverage for each individual Contact is limited to 12.7 mA. A programming interlock is enforced to limit the total output current to 20 mA or less per coverage Area. For example, a maximum current output of 12.7 mA on one Contact would limit the total summed current output on the remaining Contacts to 7.3 mA within one coverage Area.

3 The Rate is limited to 255 Hz for a given Area. The global Rate limit for each Lead is also 255 Hz.

4 The use of lower Pulse Widths than established (60 μ s to 450 μ s) is the sole responsibility of the user.

Charge Density


Figure 24. Charge Density Limits for Boston Scientific DBS Leads

Figure 24 displays the recommended maximum charge density for different combinations of Amplitude (mA) and Pulse Width (μ s). The solid black line (Limit: 6 mm²) refers to all Contacts on DB-2201 Standard Lead and the cylindrical (ring) Contacts of DB-2202 Directional Leads. The dashed black line (Limit: 1.5 mm²) refers to the small directional Contacts of DB-2202 Directional Leads. These estimates of charge density are only for the Boston Scientific DBS Leads.

WARNING: Patients may have the ability to change the Amplitude with the Remote Control. The physician should set and verify the maximum and minimum Amplitude levels allowed by the Remote Control to ensure that current levels remain safe.

Energy Use Index

The Energy Use Index applies to non-rechargeable Stimulators only. The Energy Use Index gives you an estimate of longevity of the battery life for the selected Program. After the optimal settings have been identified for a Program, from the **Program Tab**, select **Program Options** and then select **Battery** to obtain the Energy Use Index.

en

Use Figure 25 and Figure 26 to identify the longevity that corresponds to this Energy Use Index. The figures take into account nominal non-therapy power consumption, including shelf-life and patient Remote Control use. If the estimate for longevity obtained is below 12 months, consider evaluating a Boston Scientific rechargeable system.

Vercise Genus P8, P16, and P32 Non-Rechargeable Stimulators


Figure 25. Longevity Estimates Based on 24 Hour Per Day Usage

Vercise PC Non-Rechargeable Stimulator


Figure 26. Longevity Estimates Based on 24 Hour Per Day Usage

Estimated Charge Time

The Estimated Charge Time applies to rechargeable Stimulators only. The Estimated Charge Times provides an estimate of the duration and frequency of charging necessary to maintain stimulation for the selected Program. After the optimal settings have been identified for a Program, from the **Program Tab**, select **Program Options** and then select **Battery** to obtain the Estimated Charge Time.

Elective Replacement Indicator (ERI) Message

You will not be able to connect to a non-rechargeable Stimulator that is nearing the end of its battery life. The CP will display the Stimulator with an ERI message and the Stimulator battery voltage as shown in Figure 27 on the **Connect Tab**. During the ERI period, the Stimulator will continue to provide stimulation; however, no changes can be made to the Stimulator settings.

Note: The ERI Message applies to non-rechargeable Stimulators only.

Note: The battery voltage shown in Figure 27 is for demonstration purposes only. The ERI indicator varies per Stimulator type.


Figure 27. ERI Message Displayed on Connect Tab

End of Service (EOS) Message

When the Stimulator has reached its end of service, stimulation will no longer be provided. The CP will display the message shown in Figure 28 on the **Connect Tab**.

Note: The EOS Message applies to non-rechargeable Stimulators only.


Figure 28. EOS Message Displayed on Connect Tab

Clinician Programmer Management

Adjusting CP Date and Time

If the date and/or time are incorrect, select the **Adjust Date and Time** icon  on the desktop to launch the Date and Time adjustment window. Select the **Change date and time** button to modify the date and/or time as desired and select **OK** to confirm changes, as shown in Figure 29.


Figure 29. Date and Time Screens

Resetting ClinicUser Password

Update Password

If you know your ClinicUser password and would like to change it to a new password, follow the steps below:

1. Press and hold **CTRL + ALT + DEL** (control + alt + delete) on the keyboard and select **Change a password**.
2. Enter your old password in the “Old password” field, and enter your new password in the “New password” and “Confirm password” fields. The ClinicUser password must be 10 or more characters in length.
3. Ensure your new password is the same in each field. Then press **Enter**.

Lost/Forgotten Password

If you are unable to log in to the ClinicUser profile on the CP due to a lost or forgotten password, follow the steps below to reset your ClinicUser password:

1. From the CP login screen, select the PasswordReset user profile (Figure 30).


Figure 30. Login Screen with PasswordReset Profile

2. Refer to the “*Technical Support*” section of this manual and call your local Technical Support team. Technical Support will provide a password that will allow you to log in to the PasswordReset user profile.
3. Enter the password provided by Technical Support to log into the PasswordReset user profile.
4. When prompted, enter a new ClinicUser password in both the “New Password” and “Confirm Password” fields. The new ClinicUser password must be 10 or more characters in length. Ensure the ClinicUser password is the same in each field. Then select **Change Password** (Figure 31).

Note: If desired, you may exit the PasswordReset user profile without changing your ClinicUser password by selecting the  icon at the top right of the screen.


Figure 31. Update Password

5. Once the ClinicUser password has been changed successfully, you will see a message indicating the change was successful. Select **OK** to return to the login screen.
6. Using your new password, log in to the ClinicUser profile.

Software Installation, Uninstallation, and Removal

For instructions on how to install, uninstall, or remove programming software, refer to the *Software Installation Guide* for your Boston Scientific DBS System as listed in your *DBS Reference Guide*.

Software Licenses

Security.Cryptography

en

Copyright© 2017 Microsoft

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

HIDSharp

Copyright 2010-2019 James F. Bellinger

Licensed under the Apache License, Version 2.0 (the "License"); you may not use this file except in compliance with the License.

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

Apache License
Version 2.0, January 2004
<http://www.apache.org/licenses/>

TERMS AND CONDITIONS FOR USE, REPRODUCTION, AND DISTRIBUTION

1. Definitions.

"License" shall mean the terms and conditions for use, reproduction, and distribution as defined by Sections 1 through 9 of this document.

"Licensor" shall mean the copyright owner or entity authorized by the copyright owner that is granting the License.

"Legal Entity" shall mean the union of the acting entity and all other entities that control, are controlled by, or are under common control with that entity. For the purposes of this definition, "control" means (i) the power, direct or indirect, to cause the direction or management of such entity, whether by contract or otherwise, or (ii) ownership of fifty percent (50%) or more of the outstanding shares, or (iii) beneficial ownership of such entity.

"You" (or "Your") shall mean an individual or Legal Entity exercising permissions granted by this License.

"Source" form shall mean the preferred form for making modifications, including but not limited to software source code, documentation source, and configuration files.

"Object" form shall mean any form resulting from mechanical transformation or translation of a Source form, including but not limited to compiled object code, generated documentation, and conversions to other media types.

"Work" shall mean the work of authorship, whether in Source or Object form, made available under the License, as indicated by a copyright notice that is included in or attached to the work (an example is provided in the Appendix below).

"Derivative Works" shall mean any work, whether in Source or Object form, that is based on (or derived from) the Work and for which the editorial revisions, annotations, elaborations, or other modifications represent, as a whole, an original work of authorship. For the purposes of this License, Derivative Works shall not include works that remain separable from, or merely link (or bind by name) to the interfaces of, the Work and Derivative Works thereof.

"Contribution" shall mean any work of authorship, including the original version of the Work and any modifications or additions to that Work or Derivative Works thereof, that is intentionally submitted to Licensor for inclusion in the Work by the copyright owner or by an individual or Legal Entity authorized to submit on behalf of the copyright owner. For the purposes of this definition, "submitted" means any form of electronic, verbal, or written communication sent to the Licensor or its representatives, including but not limited to communication on electronic mailing lists, source code control systems, and issue tracking systems that are managed by, or on behalf of, the Licensor for the purpose of discussing and improving the Work, but excluding communication that is conspicuously marked or otherwise designated in writing by the copyright owner as "Not a Contribution."

"Contributor" shall mean Licensor and any individual or Legal Entity on behalf of whom a Contribution has been received by Licensor and subsequently incorporated within the Work.

2. **Grant of Copyright License.** Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable copyright license to reproduce, prepare Derivative Works of, publicly display, publicly perform, sublicense, and distribute the Work and such Derivative Works in Source or Object form.
3. **Grant of Patent License.** Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in this section) patent license to make, have made, use, offer to sell, sell, import, and otherwise transfer the Work, where such license applies only to those patent claims licensable by such Contributor that are necessarily infringed by their Contribution(s) alone or by combination of their Contribution(s) with the Work to which such Contribution(s) was submitted. If You institute patent litigation against any entity (including a cross-claim or counterclaim in a lawsuit) alleging that the Work or a Contribution incorporated within the Work constitutes direct or contributory patent infringement, then any patent licenses granted to You under this License for that Work shall terminate as of the date such litigation is filed.
4. **Redistribution.** You may reproduce and distribute copies of the Work or Derivative Works thereof in any medium, with or without modifications, and in Source or Object form, provided that You meet the following conditions:
 - (a). You must give any other recipients of the Work or Derivative Works a copy of this License; and
 - (b). You must cause any modified files to carry prominent notices stating that You changed the files; and
 - (c). You must retain, in the Source form of any Derivative Works that You distribute, all copyright, patent, trademark, and attribution notices from the Source form of the Work, excluding those notices that do not pertain to any part of the Derivative Works; and
 - (d). If the Work includes a "NOTICE" text file as part of its distribution, then any Derivative Works that You distribute must include a readable copy of the attribution notices contained within such NOTICE file, excluding those notices that do not pertain to any part of the Derivative Works, in at least one of the following places: within a NOTICE text file distributed as part of the Derivative Works; within the Source form or documentation, if provided along with the Derivative Works; or, within a display generated by the Derivative Works, if and wherever such third-party notices normally appear. The contents of the NOTICE file are for informational purposes only and do not modify the License. You may add Your own attribution notices within Derivative Works that You distribute, alongside or as an addendum to the NOTICE text from the Work, provided that such additional attribution notices cannot be construed as modifying the License.

You may add Your own copyright statement to Your modifications and may provide additional or different license terms and conditions for use, reproduction, or distribution of Your modifications, or for any such Derivative Works as a whole, provided Your use, reproduction, and distribution of the Work otherwise complies with the conditions stated in this License.

5. **Submission of Contributions.** Unless You explicitly state otherwise, any Contribution intentionally submitted for inclusion in the Work by You to the Licensor shall be under the terms and conditions of this License, without any additional terms or conditions. Notwithstanding the above, nothing herein shall supersede or modify the terms of any separate license agreement you may have executed with Licensor regarding such Contributions.
6. **Trademarks.** This License does not grant permission to use the trade names, trademarks, service marks, or product names of the Licensor, except as required for reasonable and customary use in describing the origin of the Work and reproducing the content of the NOTICE file.
7. **Disclaimer of Warranty.** Unless required by applicable law or agreed to in writing, Licensor provides the Work (and each Contributor provides its Contributions) on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied, including, without limitation, any warranties or conditions of TITLE, NON-INFRINGEMENT, MERCHANTABILITY, or FITNESS FOR A PARTICULAR PURPOSE. You are solely responsible for determining the appropriateness of using or redistributing the Work and assume any risks associated with Your exercise of permissions under this License.
8. **Limitation of Liability.** In no event and under no legal theory, whether in tort (including negligence), contract, or otherwise, unless required by applicable law (such as deliberate and grossly negligent acts) or agreed to in writing, shall any Contributor be liable to You for damages, including any direct, indirect, special, incidental, or consequential damages of any character arising as a result of this License or out of the use or inability to use the Work (including but not limited to damages for loss of goodwill, work stoppage, computer failure or malfunction, or any and all other commercial damages or losses), even if such Contributor has been advised of the possibility of such damages.
9. **Accepting Warranty or Additional Liability.** While redistributing the Work or Derivative Works thereof, You may choose to offer, and charge a fee for, acceptance of support, warranty, indemnity, or other liability obligations and/or rights consistent with this License. However, in accepting such obligations, You may act only on Your own behalf and on Your sole responsibility, not on behalf of any other Contributor, and only if You agree to indemnify, defend, and hold each Contributor harmless for any liability incurred by, or claims asserted against, such Contributor by reason of your accepting any such warranty or additional liability.

END OF TERMS AND CONDITIONS

Newtonsoft Json.NET

Copyright© 2007 James Newton-King

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

SharpDX

Copyright© 2010-2014 SharpDX - Alexandre Mutel

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

vNext

Copyright (c) 2015-2019, App vNext

All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of App vNext nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL <COPYRIGHT HOLDER> BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

Technical Support

Boston Scientific Corporation has highly trained service professionals to assist you. The Technical Support Department is available to provide technical consultation 24 hours a day. To speak to a representative, please choose your locality from the following list:

Argentina

T: +5411 4896 8556 F: +5411 4896 8550

Australia / New Zealand

T: 1800 676 133 F: 1800 836 666

Austria

T: +43 1 60 810 F: +43 1 60 810 60

Balkans

T: 0030 210 95 37 890 F: 0030 210 95 79 836

Belgium

T: 080094 494 F: 080093 343

Brazil

T: +55 11 5853 2244 F: +55 11 5853 2663

Bulgaria

T: +359 2 986 50 48 F: +359 2 986 57 09

Canada

T: +1 888 359 9691 F: +1 888 575 7396

Chile

T: +562 445 4904 F: +562 445 4915

China – Beijing

T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou

T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai

T: +86 21 6391 5600 F: +86 21 6391 5100

Colombia

T: +57 1 629 5045 F: +57 1 629 5082

Czech Republic

T: +420 2 3536 2911 F: +420 2 3536 4334

Denmark

T: 80 30 80 02 F: 80 30 80 05

Finland

T: 020 762 88 82 F: 020 762 88 83

France

T: +33 (0) 1 39 30 97 00 F: +33 (0) 1 39 30 97 99

Germany

T: 0800 072 3301 F: 0800 072 3319

Greece

T: +30 210 95 42401 F: +30 210 95 42420

Hong Kong

T: +852 2960 7100 F: +852 2563 5276

Hungary

T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore

T: +91 80 5112 1104/5 F: +91 80 5112 1106

India – Chennai

T: +91 44 2648 0318 F: +91 44 2641 4695

India – Delhi

T: +91 11 2618 0445/6 F: +91 11 2618 1024

India – Mumbai

T: +91 22 5677 8844 F: +91 22 2617 2783

Italy

T: +39 010 60 60 1 F: +39 010 60 60 200

Korea

T: +82 2 3476 2121 F: +82 2 3476 1776

Malaysia

T: +60 3 7957 4266 F: +60 3 7957 4866

Mexico

T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East / Gulf / North Africa

T: +961 1 805 282 F: +961 1 805 445

The Netherlands

T: +31 30 602 5555 F: +31 30 602 5560

Norway

T: 800 104 04 F: 800 101 90

Philippines

T: +63 2 687 3239 F: +63 2 687 3047

Poland

T: +48 22 435 1414 F: +48 22 435 1410

Portugal

T: +351 21 3801243 F: +351 21 3801240

Singapore

T: +65 6418 8888 F: +65 6418 8899

South Africa

T: +27 11 840 8600 F: +27 11 463 6077

Spain

T: +34 901 11 12 15 F: +34 902 26 78 66

Sweden

T: 020 65 25 30 F: 020 55 25 35

Switzerland

T: 0800 826 786 F: 0800 826 787

Taiwan

T: +886 2 2747 7278 F: +886 2 2747 7270

Thailand

T: +66 2 2654 3810 F: +66 2 2654 3818

Turkey – Istanbul

T: +90 216 464 3666 F: +90 216 464 3677

Uruguay

T: +59 82 900 6212 F: +59 82 900 6212

UK & Eire

T: +44 844 800 4512 F: +44 844 800 4513

Venezuela

T: +58 212 959 8106 F: +58 212 959 5328

Note: Phone numbers and fax numbers may change. For the most current contact information, please refer to our website at <http://www.bostonscientific-international.com/> or write to the following address:

Boston Scientific Neuromodulation Corporation

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Copyright© 2010-2014 SharpDX – Alexandre Mutel
Copyright© 2010–2014 SharpDX – Alexandre Mutel
Copyright© 2010–2014 SharpDX - Alexandre Mutel
Boston Scientific (Australia) Pty Ltd
Boston Scientific Limited
Corporation or its affiliates.