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## Mischbettwasservollentsalzer

### Operating instructions

Mixed-bed ion exchanger

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**Read these operating instructions carefully before installing and starting up the system!**

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# Preface

Dear Sir or Madam,

The stakpure team thanks you for the confidence you have placed in us.

With your selection of this mixed bed ion exchanger system you have chosen an innovative and high-quality product that will provide a long service life.

Please study the information that is given in these operating instructions on installation and operation before you start to install and start this water purification system.

Please also only appoint informed and trained personnel to operate the system and only use stakpure original accessories, replacement parts and consumables.

The above requests are particularly important as we, the manufacturer, cannot accept liability for any damages to the instrument, connected peripherals, building or persons that occur because of improper installation, operation or use of the system, or of the use of parts of any type of other manufacturers.

We wish you every success with and because of your new water purification system!

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## 12 General notes and safety precautions

- Heed these notes and precautions for your own well-being and proper functioning of the cartridge!
- Read these operating instructions carefully before installing and starting up the system and pay particular attention to the safety precautions.
- Improper use of this system or use of it in any other way than that described in this user manual invalidates the guarantee.
- When a constructional change is made to the system, or a product of a different manufacture is installed in it, the CE-mark is invalidated.
- The intended purpose of mixed-bed ion exchangers DS 1500 – DS 15000 is to deionize potable water.
- 32° C is the maximum permissible temperature for the potable water used as feedwater for mixed-bed ion exchangers DS 1500 – DS 15000. The guarantee does not cover damage to joints, gaskets and/or ion-exchanger material caused by water of a higher temperature (leakages).
- The room in which the water deionizer is installed must be frost-free and the temperature in it is never to go below +2° C.
- The installation area must have a floor drain with at least DN 50 pipe. When no floor drain, is available, a water watcher (article no.:13012900) is to be installed. The manufacturer is not liable for damages caused by overflow of water.
- For safety reasons, shut the water supply off when the system is not required for a longer period (weekends or holidays).
- Do not open the central cover lid (D) in pressure conditions on the cartridge. The opening should only be carried out by an authorized regeneration dealer.
- We recommend the additional installation of a back flow preventer and pipe cut-off combination to prevent return flow to the public potable water supply, as is required by the DIN EN 1717 part 4 regulations for potable water installations.
- A pressure vessel certificate is not required because of an appropriate change in the pressure vessels regulations of 21.04.89.
- Follow all currently appropriate laws and regulations, including those on accident prevention, which are in force at the location where the system is to be installed.
- The installation of a particle filter in the ion exchanger outlet is recommended to protect downstream users from possible harm.
- The validity term of the guarantee is 12 month.
- Under some circumstances, a defective cartridge could allow ion exchange resin to leak out.  
**CAUTION: Danger of slipping on leaked-out resin!**  
 As the resin material can cause slight irritation on contact with skin and/or eyes, immediately turn off the raw water tap to stop further leakage and proceed further as follows: On skin contact, immediately wash contacted skin areas with plenty of water. On contact with eyes, remove contact lenses if this is possible, then immediately and carefully wash eyes with plenty of water for about 10 minutes. When irritation still occurs, seek medical advice immediately. Wear protective gloves to collect up the leaked out resin, fill it into a clearly labelled plastic container and return it to the authorized specialist supplier for waste disposal.

## 13 Fields of application and process description

### Application

All natural waters contain salts that cause deposits which negatively affect numerous analyses or other applications. Mixed-bed ion exchangers remove such salts from water to a very great extent.

Fields of application include:

- Dispensing chemist's
- Printeries
- Laboratories
- Laboratory washing machines
- Humidifiers
- Medical technology
- Hospitals
- Battery water producers
- Reverse osmosis product water
- Pharmaceutical industry
- Industrial cleaning plants
- Steam generators
- etc.

### Process description

Stackpure mixed bed ion exchangers work according to the ion exchange principle. The dissolved salts that are present in the water are in the form of ions. The ion exchange resins exchange these for  $H^+$  and  $OH^-$  ions. The result is pure, demineralized water. The resins used have a certain ion exchange capacity, whereby the volume of raw water that can be processed depends on the total salt content of the raw water. A conductivity meter that is fitted to the outlet of the cartridge containing the resins shows the degree of resin exhaustion. When the reading reaches the limiting value, the ion exchange capacity has been reached. The cartridge containing the ion exchange resins is then to be sent to a central regeneration station for appropriate regeneration treatment and return.

## 14 Determination of the exchange capacity

### Feed water quality

With a feed-water quality of  $1^\circ \text{dH} = 1.25 \text{ eH}$  total dissolved solids (approx.  $30 \mu\text{S}/\text{cm}$ ), the exchange capacity is ten times the type number of the particular ion exchanger type. Using DS 4000 as example:

DS 4000 at  $30 \mu\text{S}/\text{cm}$  (approx.  $1.25^\circ \text{eH}$  total dissolved solids) = 40 000 liters  
 DS 4000 at  $300 \mu\text{S}/\text{cm}$  (approx.  $12.5^\circ \text{eH}$  total dissolved solids) = 4000 liters  
 DS 4000 at  $600 \mu\text{S}/\text{cm}$  (approx.  $25^\circ \text{eH}$  total dissolved solids) = 2000 liters

### Please note:

The exchange capacity is lowered by high amounts of carbonic acid and silicic acid in the feedwater.

### Pure water quality

The pure water taken from the ion exchanger is free of salts and carbon dioxide. It therefore has a pH of approximately pH 7. When a lower pH is measured in the pure water, this is usually the result of measurement error.

The reason for this is that pure water has a craving for carbon dioxide and can take up large amounts of it. Even when it only takes up extremely low amounts, the pH is lowered. Only flow-through pH measurement procedures that are carried out under the exclusion of air bring correct results. It is clear from this that pH measurements with indicator paper do not provide reliable results.

## 15 Note on the disposal of waste equipment

According to your state government requirements and the 2002/96 EC and 2006/66/EC directives, equipment that is to be scrapped can be brought to authorized collection points for recycling. Alternatively, it can be returned to us for proper recycling/waste disposal. In the case of return for regeneration, repair, incorrect delivery or double delivery, please use the original cardboard box/packaging whenever possible. Send fall and knock protected.

## 16 Installation and operation

### Site prerequisites:

Install the mixed bed ion exchange cartridge so that it is rigidly held. Ensure that it is freely accessible.

1. The installation area must be frost-free, minimum temperature +2° C
2. Potable water stopcock with R ¾" male thread
3. DN 50 waste pipe or water watcher with sensor
4. No source of heat in the vicinity of the ion exchange system or peripheral equipment
5. Use only the hoses supplied. The use of materials that are not suitable for pure water, (e.g. non-ferrous metal) could possibly lead to deterioration over time.
6. The cartridges are solely for connection to a potable water supply.

### Connection:

7. Remove the cartridge from the cardboard box and save the box for return of the exhausted cartridge for regeneration.
8. Remove the closing caps from the cartridge connectors and save them for fitting them to the exhausted cartridge.
9. Fit the dosing orifice (only supplied in the additional operating instructions consignment for cartridges DS 1500 to DS 4000) in the hose connecting nipple of the raw water hose.

**CAUTION: Proper system functioning and the high quality of the pure water can only be attained with inserted dosing orifice position.**

10. Screw the connecting nipple of the raw-water hose with inserted dosing orifice to the potable water shut off valve.
11. Connect the free end of the raw-water hose to the connector labelled **INLET (A)** of the mixed-bed ion exchange cartridge. Only hand tighten the hose nipple.
- 11a. When the conductivity meter has a built-in measuring cell: Screw the conductivity meter with built-in measuring cell to the cartridge connector

labelled **OUTLET (B)**. Connect the free connector of the conductivity meter with the pure water hose that is supplied.

- 11b. When the conductivity meter has a separate measuring cell: Screw the measuring cell to the connector of the mixed-bed ion exchanger labelled **OUTLET (B)**. Connect the free connector of the measuring cell to the pure water hose that is supplied. Only hand-tighten the hose nipple.
12. Connect the free end of the pure water hose to your user.

**CAUTION: Do not mix up the raw water hose INLET (A) with the pure water hose OUTLET (B) as the high quality of the product water would then not be attainable.**

### 16.1 Putting into operation

Air-vent the cartridge:

13. Open the vent-valve of the central cover lid of the cartridge
14. Carefully open the raw water tap and wait until water emerges from the vent-valve.

**CAUTION: Do not screw the vent valve out of the thread but just turn it anti-clockwise until air emerges.**

15. When water emergencies, close the vent-valve.
16. The mixed-bed ion exchanger is now ready for use.
17. The quality of the water that is produced can be read off from the conductivity meter.

**CAUTION: The values of the conductivity are irrelevant when no product water is drawn out. True values are only shown when water flows through the measuring cell.**

18. The mixed-bed ion exchanger cartridge must be regenerated when the maximum conductivity value for the particular application is reached.

## 17 Regeneration

We recommend to keep an additional cartridge in reserve to continually have pure water at your availability, even when sending a cartridge for regeneration. Please send the exhausted cartridge (in

its original packaging whenever possible) to the following address:

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\*On request, the cartridges can also be equipped with quick-connect couplings.

## 19 Specifications DS 7000 – DS 15000

<b>stakpure mixed bed ion exchanger</b>	<b>DS 7000</b>	<b>DS 11000</b>	<b>DS 15000</b>
Water demineralizer			
Stainless-steel cartridge			
Flow rate	2000 l/h	2500 l/h	3000 l/h
Capacity at 10°dH total dissolved solids	7000 l	11000 l	15000 l
Pure water quality			
Max. water temperature			
Material			
Dimensions mm (∅ x height)	363 x 600	363 x 850	363 x 1100
Connection*			
Weight, approx.	55 kg	70 kg	90 kg
<b>Order-No.</b>	<b>12700000</b>	<b>12110000</b>	<b>12150000</b>

\*On request, the cartridges can also be equipped with quick-connect couplings.

## 20 Accessories for ion exchangers

<b>Designation</b>	<b>Order no.</b>
Conversion set to quick connect nipple, nickel-plated	13150000
Conversion set to quick connect coupling hose, nickel-plated	13150600
Conversion set to R3/4" quick connect nipple with check valve	13150700
Conversion set to quick connect coupling – hose with check valve	13150800
2 x 1.5m hose, DN 10 R 3/4"	13169000
2 x 1.5m reinforced hose, DN 20 R 3/4"	13169100
Wall bracket for cartridge types DS 1500 – DS 4000	13140500
Water watcher with solenoid valve and sensor R 3/4"	13012900
Pressure reducer, 0 – 4 bar R 3/4"	13141400
Solenoid valve, nickel plated R 3/4"	13150300
Pure water tap, straight nozzle R 3/4"	13140000
Pure water tap, angled nozzle R 3/4"	13140100
Pure water shut off valve R 3/4"	13140300
Pure water two-way distributor R 3/4"	13142000
Pure water three-way distributor R 3/4"	13140200

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
## 22 EC-Declaration of Conformity

in the sense of **directive 97/23/EC** for pressure vessels, acc. to **article 3, paragraph 3 in 97/23/EC**.

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We hereby certify that the conception and form of the machine named below that we have brought to market are in accordance the determinations of the named article of the EC directive. This declaration is no longer valid when changes are made to the product.

**Product designation:** Mixed bed ion exchanger

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**Art. No.:** 12115000, 12200000, 12200102, 12280002, 12280102, 12400000, 12600000, 12700000, 12110000, 12150000

**The following harmonised standards were applied:**


DIN EN 287-1, DIN EN ISO 15614-1

**The following national standards, directives and specifications were applied:**

AD 2000 set of regulations

The technical documentation is completely available. The operating instructions for the product are available in the national language of the user.

Niederahr, den 24. 03. 2014

---

Manufacturer/CE authorised person: Leo Trumm

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**Manufacturer:** stakpure GmbH, Auf dem Kesseling 11, D-56414 Niederahr
**stakpure GmbH**