Akusto-optischer Modulator

Fortgeschrittenenpraktikum A 4.10

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Name:

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Datum:

Criterion for Bragg and Raman-Nath diffraction regimes

M. G. Moharam and L. Young

The idea is well entrenched in the literature that thin phase gratings (whether holographic or acoustically induced) should exhibit Raman-Nath behavior (and thus give several diffracted waves), and that thick phase gratings should show Bragg behavior (one diffracted beam and that only for Bragg angle incidence). The parameter Q of Klein and Cook, which is a normalized measure of grating thickness, has been extensively used as a criterion for deciding which regime will apply. It is perhaps not generally realized that Q is not a reliable parameter for this purpose but requires, as indeed Klein and Cook noted, a limitation on grating strength. This limitation is a matter of practical concern. For example, we have observed Raman-Nath behavior with Fe-doped LiNbO₃ even for very large values of Q . The purpose of the present paper is to note that a parameter ρ (first defined by Nath) is an effective replacement for Q , since ρ is reliable and Q is not. ρ is defined as $\lambda_0^2/\Lambda^2 n_0 n_1$, where λ_0 is the vacuum wavelength of the light, Λ is the grating spacing, n_0 is the mean refractive index, and n_1 is the amplitude of the sinusoidal modulation of the refractive index. The grating thickness does not enter ρ , so the terms thin and thick are, strictly speaking, irrelevant to the question of which regime is operative. However, thin enough gratings will tend to operate in the Raman-Nath regime because the index modulation must be large for a thin grating to produce appreciable diffraction.

Introduction

The phenomenon of light diffraction by periodic phase gratings, whether holographically or acoustically produced, has been treated extensively by many authors.¹⁻¹⁵ There is general agreement that it is convenient to define two regimes in which phase gratings operate. In the Raman-Nath regime, several diffracted waves are produced. In the Bragg regime, essentially only one diffracted wave is produced, and that only for near Bragg incidence. It has been customary to refer to gratings which operate in the Raman-Nath or the Bragg regimes as thin and thick gratings, respectively. Klein and Cook⁴ introduced a parameter Q , defined as $2\pi\lambda_0 L/\Lambda^2 n_0$ (where L is the grating thickness, λ_0 is the vacuum wavelength of light, Λ is the grating spacing, and n_0 is the mean refractive index) to distinguish between the two diffraction regimes: Values of $Q < 1$, i.e., thin gratings, were believed to give Raman-Nath operation. Large values of Q ($Q > 10$), i.e., thick gratings, were believed to give Bragg regime operation. The use of Q is well entrenched in the literature. However, as

Klein and Cook noted, their analysis was based on the assumption that $\nu = 2\pi n_1 L/\lambda_0$ remains less than six; and if ν exceeds this value, it is necessary to place more restrictive conditions on Q (which they proceeded to discuss). Actually, it can be shown that Q may fail with ν as low as three. Recently, Magnusson and Gaylord⁸ have shown theoretically that for a large modulation, higher order diffracted waves become important (Raman-Nath regime) even for large Q . Conversely, for small modulation, only a single wave is diffracted in spite of small values of Q . Kaspar⁷ came to a similar conclusion when he compared his theory for a grating with a complex dielectric constant with the coupled wave theory of Kogelnik.³ Alferness⁹ has discussed qualitatively the effect of the magnitude of the modulation of the refractive index on the validity of Q as a tool to predict the regime of diffraction. Various authors^{8,9,16} have observed several diffracted waves during hologram recording in lithium niobate and dichromated gelatine when Q was large. We have observed up to eight diffracted waves during hologram storage in 1-cm thick Fe-doped lithium niobate crystal. For this experiment, $Q = 55$.

These experimental observations, in addition to the theoretical predictions which show that Q does not work, led us to reconsider the basis of the use of Q and to seek to find a parameter which could be used to predict the regime of diffraction. We will show that a parameter ρ defined as $\lambda_0^2/\Lambda^2 n_0 n_1$ (where n_1 is the mod-

The authors are with University of British Columbia, Electrical Engineering Department, Vancouver, B.C., V6T 1W5.

Received 30 November 1977.

0003-6935/78/0601-1757\$0.50/0.

© 1978 Optical Society of America.

ulation of the refractive index and λ_0 , Λ , and n_0 are as defined before) can be reliably used to predict whether one is in the Raman-Nath regime or in the Bragg regime. The parameter ρ was defined first by Nath,¹² who was considering the case of normal incidence. Nath pointed out that, if ρ is very large, the diffraction effect will not be prominent, as is otherwise the case where ρ is nearly zero. Phariseau¹³ and later Bergstein and Kermisch⁵ and Chu and Tamir⁶ have obtained expressions for the intensities of the first few diffracted orders near the first Bragg incidence^{5,13} and higher order Bragg angle⁶ assuming that ρ is large (Bergstein and Kermisch⁵ and Chu and Tamir⁶ defined a parameter equal to $1/\rho$ called s and q in successive papers.) Those authors^{5,6,13} have used $\rho \gg 1$ (or $1/\rho \ll 1$) in mathematical approximations to obtain analytical expressions of the intensities. They did not use it to predict whether one is in the Raman-Nath or the Bragg regime. In fact, they did not address either the problem of determining the diffraction regime or the problem of the effectiveness of Q in distinguishing between the two regimes.

Theory

The problem of light diffraction by periodic phase gratings has, of course, been treated in detail in many papers.¹⁻¹⁵ We will only give a brief outline here.

The scalar wave equation is

$$\nabla^2 E + k^2 E = 0, \quad (1)$$

where $k^2 = \beta^2 - j\alpha\beta$ and $E(x,z,t)$ is the complex amplitude of the y component of the field. $\beta = 2\pi n/\lambda_0$ is the propagation constant, and n is the refractive index. A sinusoidal modulation of the refractive index is assumed.

$$n = n_0 + n_1 \cos(2\pi x/\Lambda). \quad (2)$$

Here, n_0 , n_1 , λ_0 , and Λ are as defined before. An unslanted grating is assumed for simplicity. The electric field E may be represented by its Fourier expansion as

$$E = \sum_{l=-\infty}^{\infty} \phi_l(z) \exp(-j\bar{\sigma}_l \cdot \bar{r}), \quad (3)$$

where $\bar{\sigma}_l = \bar{\sigma}_0 - l\bar{K}$, where \bar{K} is the grating vector and $K = 2\pi/\Lambda$, $\bar{\sigma}_0$ is the wave vector of the incident wave ($l = 0$), and $\bar{r} = (x,y,z)$.

Substituting Eqs. (2) and (3) in Eq. (1) and neglecting the second derivative with respect to z , we arrive at the well-known set of coupled wave equations:

$$\frac{\partial \phi_l}{\partial \xi} = j\rho l^2(1 - B_l)\phi_l + j(\phi_{l-1} + \phi_{l+1}), \quad (4)$$

where $B_l = 2\Lambda \sin\theta_0/\lambda_0 l$, θ_0 is the external angle of incidence, $\xi = \nu(z/L)$, $\nu = \pi n_1 L/\lambda_0 \cos\theta_1$, and θ_1 is the angle of diffraction of the zero-order mode. Thus, B_l is equal to 1 for Bragg diffraction forming the l th mode, and ν is a measure of the grating strength (and appears, e.g., in the coupled mode theory of Kogelnik for Bragg diffraction³). The absorption constant α may be neglected without loss of generality since it can be allowed for by the substitution $\phi_l' = \phi_l \exp(-\alpha z/\cos\theta_1)$.

Equation (4) shows that the l th mode is coupled to itself and to two adjacent modes ($l - 1$ th and $l + 1$ th modes). Effective energy transfer between modes requires essentially that the factor $\rho l^2(1 - B_l)$ be relatively small since, if this factor is much larger than 1, almost all the energy will be coupled back to the l th mode. [This can be seen from Eq. (1) by neglecting the second term of the right-hand side.] Therefore, if $\rho < 1$, appreciable energy may be transferred successfully to higher order modes up to some value of l provided the magnitude of B_l is appropriately limited. The number of higher order modes depends on how small ρ is. (The smaller ρ is, the larger the number of higher order modes.) If $\rho = 0$, Eq. (1) gives the well-known solution in terms of Bessel functions $\phi_l = j^l J_l(2\nu)$. This solution is often obtained by Fourier expansion of the transmitted wave with spatially sinusoidal phase modulations. A similar solution was shown by Klein and Cook for $Q = 0$. However, as the thickness L of the grating goes to zero ($Q \rightarrow 0$) the modulation of the refractive index must go to infinity to retain the finite phase shift ($n_1 L$); and as $n_1 \rightarrow \infty$, $\rho \rightarrow 0$. Clearly, the case where Q or $\rho = 0$ is a nonphysical situation. For $\rho < 1$, the diffraction process is in the Raman-Nath regime. If $\rho \gg 1$, appreciable energy may be transferred only into the mode for which the Bragg conditions holds or nearly holds, i.e., the mode with l such that $B_l \simeq 1$ for the given θ_0 . All the above applies to an infinite plane wave incident on an infinite grating, as in previous work.⁴

The above analysis was illustrated by solving the set


Fig. 1. The intensity of four modes vs the grating strength ν for $\rho = 1$ and $B_l = 1/l$ (i.e., Bragg incidence for the first-order mode). The sum of the intensities of the four modes shown was close to unity for this case, all the other modes being negligibly small.


Fig. 2. The intensity of the zero- and first-order modes vs the grating strength ν for $\rho = 50$ and $B_l = 1/l$. All the other modes are negligible. For this value of ρ , $Q = 100 \nu$.

of Eq. (1), using the UBC computer center, fourth-order Runge-Kutta routine with error control, for the boundary conditions $\phi_0(0) = 1$, $\phi_l(0) = 0$ for $l \neq 0$. The intensity of the l th mode is measured by $\phi_l^*(\nu)\phi_l(\nu)$. The number of modes included (i.e., the number of equations) was increased for each case until no further significant effect occurred for a given value of ρ , e.g., for $\rho = 50$, seven modes were included, and for $\rho = 0.1$, twenty-three modes were needed.

Case I: $\rho \leq 1$ (Raman-Nath regime)

Figure 1 shows the intensities of several modes for $\rho = 1$ as a function of the grating strength ν . Bragg incidence for the first mode was assumed (i.e., $B_l = 1$ for $l = 1$). Clearly, even for values of ν corresponding to large values of $Q = 2\nu\rho \cos\theta_1$, the Raman-Nath regime applies; and energy is transferred successively into a large number of modes.

Case II: $\rho \gg 1$ (Bragg regime)

Figure 2 shows a plot of the intensity of the zero- and first-order modes against ν (and by implication Q) for $\rho = 50$ and Bragg incidence for the first-order mode. The energy is transferred to and fro between these two modes as $\cos^2\nu$ and $\sin^2\nu$ in agreement with Kogelnik's coupled wave theory.³ Phariseau¹³ indeed showed that intensities of the higher order modes are of the order of $1/\rho^2$ or less of that of the incident mode (provided that $\rho \gg 1$). Figure 2 shows that the Bragg regime operates for all Q given that $\rho \gg 1$. As ρ becomes larger, the grating becomes more selective, i.e., B_l must be closer to 1 so that the product $\rho l^2 (1 - B_l)$ is small enough that the energy may be transferred to the l th mode.

Conclusions

The factor ρ has the advantage over Q as a criterion for whether the Raman-Nath or the Bragg regime will be observed with a given phase grating, that ρ is always reliable and Q is not. A large Q is not. A large ρ favors the Bragg regime. The relative light intensity going into higher order modes is of the order $1/\rho^2$, so a value of $\rho > 10$ indicates more or less ideal Bragg behavior. A further point is that the thickness L of the grating is irrelevant in itself, since it does not enter ρ . This may be explained as follows: to obtain a significant amount of diffraction, the grating strength ν has to be large enough in some sense. For ν to be large enough in this sense, then, if L is large, n_1 may be large or small; so ρ would accordingly be large or small, and we could be in either regime. On the other hand, if L is small, n_1 would have to be large to obtain a sufficiently large ν ; in which case ρ would be small, and the Raman-Nath regime would hold. Therefore, the distinction between thick and thin gratings or holograms as determined by the value of L (for given λ_0 , n_0 , and Λ) in $Q = 2\pi\lambda_0 L / \Lambda^2 n_0$ is invalid as a description of whether a single diffracted beam will be produced, and this only for light incident at or near the Bragg angle, or whether many diffracted beams will be produced regardless of the satisfaction of a Bragg angle, e.g., at normal incidence.

This work was supported by an operating grant from the National Research Council of Canada to L. Y.

References

1. C. V. Raman and N. S. N. Nath, Proc. Indian Acad. Sci. 4, 222 (1936).
2. C. B. Bruckhardt, J. Opt. Soc. Am. 56, 1502 (1966).
3. H. Kogelnik, Bell Syst. Tech. J. 48, 2909 (1969).
4. W. R. Klein and B. D. Cook, IEEE Trans. Sonics Ultrason. SU-14, 123 (1967).
5. L. Bergstein and D. Kermisch, Proc. Symp. Mod. Opt. 17, 655 (1967).
6. R. S. Chu and T. Tamir, IEEE Trans. Microwave Theory Tech. MTT-18, 486 (1970).
7. F. G. Kaspar, J. Opt. Soc. Am. 63, 37 (1973).
8. R. Magnusson and T. K. Gaylord, J. Opt. Soc. Am. 69, 1165 (1977).
9. R. Alferness, J. Opt. Soc. Am. 66, 353 (1976).
10. R. Alferness, Appl. Phys. 7, 29 (1975).
11. S. F. Su and T. K. Gaylord, J. Opt. Soc. Am. 65, 59 (1975).
12. N. S. N. Nath, Proc. Indian Acad. Sci. 8, 499 (1938).
13. P. Phariseau, Proc. Indian Acad. Sci. Sect. A 44A, 165 (1956).
14. D. Kermisch, J. Opt. Soc. Am. 59, 1409 (1969).
15. N. Uchida, J. Opt. Soc. Am. 63, 280 (1973).
16. V. E. Wood, N. F. Hartman, C. M. Verber, and R. P. Kenan, J. Appl. Phys. 46, 1214 (1975).

Double-pass acousto-optic modulator system

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Time and Frequency Division, National Institute of Standards and Technology, 325 Broadway, Boulder, Colorado 80305

(Received 2 March 2005; accepted 21 April 2005; published online 1 June 2005)

A practical problem that arises when using acousto-optic modulators (AOMs) to scan the laser frequency is the dependence of the beam diffraction angle on the modulation frequency. Alignment problems with AOM-modulated laser beams can be effectively eliminated by using the AOM in the double-pass configuration, which compensates for beam deflections. On a second pass through the AOM, the beam with its polarization rotated by 90° is deflected back such that it counterpropagates the incident laser beam and it can be separated from the input beam with a polarizing beam splitter. Here we present our design for a compact, stable, double-pass AOM with 75% double-pass diffraction efficiency and a tuning bandwidth of 68 MHz full width at half maximum for light transmitted through a single-mode fiber. The overall efficiency of the system (defined as the optical power out of the single-mode fiber divided by the optical power into the apparatus) is 60%.

[DOI: 10.1063/1.1930095]

I. INTRODUCTION

For many laser cooling and trapping experiments, it is necessary to shift and/or sweep the laser frequency at the end of a cooling cycle. This is done to either adiabatically cool the atoms to a lower temperature, or to manipulate them in some other way, such as launching them in an atomic fountain. Most often, the laser frequency must be switched on time scales shorter than 1 ms, and it may be that laser beams with several different frequencies must be produced from a single laser source.

Acousto-optic modulators (AOMs) are widely used to accomplish the frequency control in laser cooling experiments. When the laser frequency is scanned with an AOM, the angle of the first-order diffracted beam shifts as well, since the beam diffraction angle is a function of modulation frequency. Changes in beam diffraction angle may be desirable for some applications where spatially resolved diffracted beams are needed, but for many applications any change in the laser propagation direction is an unwanted side effect. Using an AOM in the double-pass configuration is a way to practically eliminate changes in beam steering during frequency sweeps and jumps within the frequency tuning bandwidth of the AOM. Here we present a detailed description of our modular double-pass system design and a study of the system's performance.

The article is organized as follows. In Sec. II, the laser frequency jump and ramp sequence that is used to cool and launch the atoms in the primary frequency standard at the National Institute of Standards and Technology (NIST) is presented in detail. This provides a real-world example of the level of laser frequency control that can be desirable for

atomic physics experiments. In Sec. III, background into the properties of acousto-optic modulators is presented. In Sec. IV, the double-pass configuration is introduced, our specific design is discussed in detail, and the performance of our design is presented.

II. ATOMIC FOUNTAIN LASER FREQUENCY CONTROL

For laser cooling and trapping experiments, the optimum laser intensity and detuning values that maximize the number of laser-cooled atoms are not the same as the optimum values for the lowest atom temperature. The laser intensity that maximizes the number of atoms can depend on the details of the atom source, but in general, a higher laser intensity is better for a thermal source. The optimum laser detuning is generally a few linewidths, but it can also depend on the details of the atom source and the laser intensity. In contrast, at low laser intensities the equilibrium atom temperature in an optical molasses is equal to a small constant term ($\sim 1 \mu\text{K}$ for ^{133}Cs) plus a term that is proportional to the laser intensity and inversely proportional to the laser detuning.^{1,2} To simultaneously maximize the atom number and minimize the atom temperature, often a brief "postcool" stage is applied as the molasses is turned off. During the postcool stage, the laser frequency is ramped further to the red of the atomic resonance while the intensity is ramped to zero.

The atom launch and postcool sequence for the vertical laser beams used in the NIST-F1 atomic fountain clock at NIST is shown in Fig. 1. NIST-F1 serves as the primary frequency standard for the United States, and a detailed description of the apparatus has been published previously.³ We accomplish the frequency sweeps and jumps shown in Fig. 1 with AOMs. The intensity ramp is produced by controlling the rf power into the AOM with a voltage-controlled attenuator. With this sequence, we are able to collect about 10^7

^{a)} Author to whom correspondence should be addressed; electronic mail: edonley@boulder.nist.gov

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
FIG. 1. Laser frequency and intensity ramps typically used for NIST-F1. The intensity is shown on the left vertical scale in units of the saturation intensity, I_{sat} . The laser detuning from resonance is indicated on the right vertical scale. The atoms are initially collected in an optical molasses for ~ 250 ms. To launch the atoms in a moving molasses, the frequencies of the up (down) beams are shifted from their initial frequencies by typically $+2.5$ MHz (-2.5 MHz) for 1 ms. Then the atoms are postcooled in the moving frame for 2 ms by ramping down the laser intensity while linearly shifting the frequencies of both beams further to the red by 44 MHz. Note that the frequency and intensity curves for the horizontal molasses beams are not shown. Those beams also have a postcool sequence, but there is no frequency step when the atoms are launched.

atoms per cycle, cool them to $\sim 0.5 \mu\text{K}$, and launch them into a microwave cavity structure for Ramsey interrogation.

III. ACOUSTO-OPTIC MODULATORS

A. Bragg scattering

For most practical applications of acousto-optic modulators, the Bragg description of the modulation process is a good approximation to the behavior of the system.^{4,5} The main features of AOMs are derived in this picture by treating the modulation as a photon-phonon scattering process.

The Bragg matching condition can be derived by treating the acoustic and optical fields as particles with momentum κ and k , respectively, where κ (k) is the phonon (photon) wave vector for the acoustic (optical) field; $\kappa = \Omega/v_s$, where Ω is the rf modulation frequency and v_s is the speed of sound in the crystal. Similarly, $k = \omega/v_L$, where ω is the light frequency and v_L is the speed of light in the crystal.

A scattering process between photons and phonons results in the absorption or emission of acoustic phonons. A first-order scattering process between a photon and a single phonon is described by the energy-momentum relations

$$\begin{aligned}\omega_d &= \omega_i \pm \Omega, \\ k_d &= k_i \pm \kappa.\end{aligned}\quad (1)$$

The subscripts i and d designate whether the corresponding photon is incident or diffracted. The sign depends on whether the phonon is absorbed or emitted, which depends on the relative orientations of the incident photon and phonon wave vectors.

The Bragg matching condition that determines the optimum angles for the incident laser and acoustic beams for peak first-order diffraction efficiency can be derived with energy and momentum conservation arguments. Figure 2


FIG. 2. Momentum conservation diagrams for the absorption of a phonon. The two cases [(a) and (b)] represent two configurations for the incident and diffracted light beams where the conditions are Bragg matched for a given orientation of the phonon momentum. Note that κ , and thus θ_i and θ_d , are made large in the diagrams for clarity. The Bragg angle is normally about 1° or less, and the phonon momentum is roughly two orders of magnitude smaller than the photon momentum.

shows momentum conservation diagrams that describe photon-phonon scattering events in which a phonon is absorbed (+first-order diffraction). The photon momentum of the diffracted light is equal to the sum of the momenta of the phonon and the incident photon. Conservation of energy requires that the frequency of the diffracted beam be shifted upward from ω to $\omega + \Omega$ for a phonon absorption process; but since $\omega \gg \Omega$, the frequency shift can be ignored in the momentum conservation analysis, and $|k_i| = |k_d|$. Adding together the x and y momentum components leads one to Bragg's law

$$\sin \theta_B = \frac{\kappa}{2 \cdot k_i}, \quad (2)$$

where θ_B is the Bragg angle, and $\theta_i = \theta_d = \theta_B$. Note that we have derived Eq. (2) without considering the effect of the boundaries of the acoustic medium. Often it is desirable to calculate the Bragg angle outside of the crystal, $\theta_{B,\text{ext}}$. To close approximation, if the crystal boundaries are parallel to κ , the external Bragg angle is larger than the internal Bragg angle by a factor of n , where n is the refractive index of the acoustic medium.

The cases shown in Figs. 2(a) and 2(b) represent the two geometries where the Bragg matching condition is met for +first-order phonon absorption for a given orientation of κ . Notice that the incident photon wave vector for one case and the diffracted photon wave vector for the other case counter-propagate. This point is relevant to the discussion of the double-pass configuration below. Diagrams similar to those in Fig. 2 can be drawn for the single-phonon emission process, where the diffracted light enters the -first order, and the frequency of the diffracted beam is shifted to $\omega - \Omega$.

B. Acousto-optic devices

A schematic drawing of an acousto-optic modulator is shown in Fig. 3. A rf signal is fed to a strain transducer in contact with the AOM crystal. The rf modulation at frequency Ω causes a traveling density wave to form inside the crystal; the wave propagates at the speed of sound in the crystal, v_s , with the frequency Ω . The refractive index is therefore modulated with a wavelength of $\Lambda = 2\pi \cdot v_s / \Omega$, and the crystal acts like a thick diffraction grating with the rulings traveling away from the transducer with a velocity v_s . The Bragg approximation is valid only when the acoustic wave is describable by a plane wave and all phonons have


FIG. 3. (Color online) Schematic drawing of an aligned AOM. The AOM is shown in the mode where the light is diffracted into the +first order and the laser frequency is up-shifted by the rf modulation frequency Ω . Note that phonons are absorbed in this configuration, and momentum is transferred to the laser beam in the direction away from the strain transducer. The crystal can also be aligned to be Bragg matched for the -first diffraction order (not shown). In that case, corresponding to a phonon emission process, all of the laser beams are mirror symmetric about the horizontal axis with respect to those shown in the figure, and the diffracted laser frequency is down shifted by Ω . Note that the external Bragg angle is larger than the internal Bragg angle because of refraction at the crystal surfaces.

the same wave vector. In practice, this limiting case can be achieved to a close approximation when the strain transducer is long compared to the acoustic wavelength in the direction of laser beam propagation and acoustic diffraction is minimized. It is therefore not uncommon to achieve greater than 80% diffraction efficiency into a single diffraction order with careful engineering.⁵ The AOM that we use (Crystal Technology, Inc. 3080-120)⁶ has a TeO_2 crystal driven with a longitudinal acoustic mode. At a typical rf frequency of 80 MHz the external Bragg angle is 8 mrad for light with a vacuum wavelength of 852 nm.

C. Bandwidth

The bandwidth is generally defined as the frequency range over which the modulated laser power drops by 50% (3 dB) from the peak value. For the AOM model that we use the specified modulation bandwidth is 20 MHz. The diffraction efficiency that we measured versus rf modulation frequency is presented in Fig. 4 for a single pass through the AOM of a 1.1 mm [full width (FW) $1/e^2$] collimated beam. The full width at half maximum (FWHM) of the curve is ~ 50 MHz, which is 2.5 times larger than the manufacturer's specified bandwidth. The diffraction efficiencies are the average of four measurements. The values for the absolute errors in the single-pass diffraction efficiency are up to 4%. The double-pass measurements shown in Fig. 4 are discussed in the next section.

IV. DOUBLE-PASS CONFIGURATION

We send light to our apparatus with a number of single-mode optical fibers that facilitate the delivery of laser light to hard-to-reach points on the apparatus and act as spatial filters. While the spatial filtering improves the beam quality and facilitates optics alignment downstream from a fiber, it makes the coupling efficiency for light into the fiber ultra-sensitive to the stability of the beam pointing upstream of the


FIG. 4. (Color online) Measured single-pass and double-pass diffraction efficiencies. When comparing the data in the graph, note that for the double-pass measurements the light was diffracted twice by the AOM, and the measurements of the diffracted beam power were performed at the output of a single-mode fiber, which had a coupling efficiency of 80%. It is the overall diffraction efficiency through the fiber that is plotted in the figure for the double-pass data. The double-pass diffraction efficiency at the fiber input is 75%. The double-pass diffraction efficiency measurements were performed with and without the cat's eye lens. Note that since the light frequency is shifted twice for the double-pass measurements, the actual frequency tuning bandwidth is twice the FWHM of the curves. For the double-pass measurements performed with the cat's eye, this corresponds to 68 MHz.

fiber. This presents a challenge when the frequency of the light being injected into the fiber is ramped with AOMs, since the Bragg angle depends on the modulation frequency.

The alignment problems caused by AOM frequency sweeps can be practically eliminated by using the double-pass configuration, in which the laser beam travels through the AOM twice and the beam deflection is compensated in the second pass. Modulating the light with the AOM twice also means that the total frequency shift is twice the modulation frequency.

The two momentum conservation diagrams in Figs. 2(a) and 2(b) represent the momentum vectors for the two passes. Notice that the incident and diffracted beams for one case counterpropagate the diffracted and incident beams for the other case. Thus with the alignment optimized, if the +first-order diffracted beam is reflected back onto itself through the AOM for a second pass, the +first-order diffracted beam for the second pass will spatially overlap the original incident beam, independent of the modulation frequency.

Figure 5 is a schematic drawing of our double-pass AOM apparatus, and a photograph is shown in Fig. 6. Light with linear polarization oriented parallel to the page enters the apparatus from the left as shown in Fig. 5. The beam is steered through a polarizing beam splitter and a Galilean telescope by two adjustable mirrors before it enters the AOM. The Galilean telescope reduces the beam diameter from 2.2 to 1.1 mm (FW $1/e^2$) and consists of a plano-convex lens and a plano-concave lens separated by 5 cm. The position of the AOM is constrained by a 3 mm dowel pin, around which the AOM can be rotated about an axis through the center of the crystal. The AOM can be clamped down to the mounting block with two screws after rough


FIG. 5. (Color online) Schematic drawing of the double-pass system. M: mirror; PBS: polarizing beam splitter; PCX: plano-convex lens; PCC: plano-concave lens; AOM: acousto-optic modulator; QWP: quarter-wave plate; BB: beam block; HWP: half wave plate; FC: fiber coupler; PM: polarization maintaining fiber. The circles that are attached to some of the optics represent angular adjustment knobs. The arrows above two of the lenses indicate that they are held in lens-focusing mounts.

angular alignment of the AOM is achieved. Fine angular alignment is accomplished by adjusting the orientation of the two input mirrors. Immediately downstream from the AOM is a quarter-wave plate that circularizes the polarization. The plano-convex lens and mirror downstream from the AOM combine to form the cat's eye retroreflector, which is discussed in detail below. The zeroth-order undiffracted beam is blocked with a beam block before the mirror. The +first-order diffracted beam reflects back from the cat's eye and through the quarter-wave plate, where the polarization is converted back to linear polarization such that the polarization vector is oriented perpendicular to the page on the second pass through the AOM. The light is modulated by the AOM a second time and the resulting +first-order diffracted beam for the second pass overlaps the incident laser beam. The modulated light is then reflected from the adjustable polarizing beam splitter and directed to an adjustable mirror, a half wave plate, and a polarization-maintaining single-mode fiber coupler that contains an 11 mm coupling lens. If beam misalignments of over $1 \mu\text{m}$ for the focused spot on the fiber core are to be avoided, the angular tolerance for the light beam leaving the AOM is ~ 0.2 mrad. Note that the Galilean telescope reduces the angular steering of the beam entering the fiber by a factor of 2. Alignment of the laser beam into the fiber coupler is accomplished by adjusting the angles of the mirror and the polarizing beam splitter. To avoid fluctuations of the polarization of the light emerging from the polarization-maintaining fiber, the half wave plate is


FIG. 6. (Color online) Our double-pass apparatus. The dimensions of the base plate are $25 \text{ cm} \times 46 \text{ cm}$.

needed to align the light polarization with the proper axis of the fiber. The base plate to which the optics are mounted is temperature controlled with an inexpensive commercial foil heater and controller. We describe the function of some key features in the apparatus in detail below.

A. Galilean telescope

The diffraction efficiency for an AOM is a strong function of beam waist diameter, and increases with increasing waist size before leveling off above 90% for a typical diameter of about 0.5 mm. The case enclosing the AOM that we use has laser apertures of 3 mm in diameter, but the active aperture over which the manufacturer guarantees that the device meets specifications is only 1 mm in diameter. The beam size we chose is 1.1 mm ($\text{FW } 1/e^2$) and is a trade-off between performance and convenience.

The Gaussian beam waist diameter d for a plano-convex lens may be estimated by the formula $d = 4\lambda f / \pi D$ in terms of the focal length f , the wavelength λ , and the beam diameter at the lens D . Following the manufacturer's suggestion of using a single plano-convex lens for focusing down the beam while simultaneously optimizing the diffraction efficiency would require the use of a lens with 1 m focal length. If sufficient space is available on the optical table, then we can follow this approach and achieve good diffraction efficiencies with a weakly focused beam.

We can dramatically reduce the size of the apparatus by using a Galilean telescope to reduce the sizes of the laser beams. In principle, by using lenses with very short focal lengths, the length of the telescope can be reduced to a few centimeters. However, there is a trade-off in how small to make the focal lengths. In general, it is much more difficult to fine tune the optical alignment in the apparatus when short focal-length lenses are used to resize the beams. In our apparatus, the input beam has a diameter of 2.2 mm ($\text{FW } 1/e^2$). We reduce the diameter by a factor of 2 with a +100 mm plano-convex and a -50 mm plano-concave lens spaced by ~ 50 mm such that the beam entering the AOM is collimated. We mounted the 100 mm converging lens on a commercial lens-focusing mount in order to fine tune the lens spacing. Note that proper orientation of the planar side of the lenses is required to minimize spherical aberrations⁴ (see Fig. 5). The alignment of a prototype of the apparatus was much more difficult to optimize and maintain when the Galilean telescope consisted of $f = +50$ mm and $f = -15$ mm lenses.

We also measured the single-pass diffraction efficiency of the AOM with a focused beam by removing the diverging lens, without which the beam waist diameter was $\sim 50 \mu\text{m}$ in the AOM. If the incident laser beam is strongly focused into the AOM, then there is effectively a range of angles of incidence for the incoming light, and only a fraction of the converging rays are Bragg matched with the crystal. From the point of view of geometric optics, when the beam is focused down, the light rays enter with an angular spread of ~ 20 mrad. This range is large compared to the external Bragg angle (~ 8 mrad), and in this mode only a component of the zeroth-order beam is diffracted into the first-order beam. Since the component of the zeroth-order beam that diffracts depends on the rf frequency, a dark spot is observed

moving across the zeroth-order beam as the rf frequency is swept over the frequency range of the AOM. The maximum diffraction efficiency peaked at approximately 50% of what we observed with a collimated beam.

B. Cat's eye retroreflector

Having a cat's eye⁷ as opposed to a simple flat-mirror retroreflector dramatically improves the frequency-tuning range (bandwidth) for the double-pass apparatus. Owing to the dependence of the Bragg angle on the rf modulation frequency, the angle of incidence onto the retroreflector depends on Ω . This can cause the frequency-tuning bandwidth for the double-pass apparatus to be much narrower than for single-pass diffraction, since in practice, the optical alignment is optimized at only a single rf frequency. When the beam arrives at the AOM for its second pass, an imperfect retroreflector causes the angle of incidence on the AOM to deviate from the external Bragg angle and the beam spot position to shift on the AOM. This effect is most problematic when a flat-mirror retroreflector is used.

The standard cat's eye consists of a lens and a mirror with their spacing equal to the focal length of the lens, f . For our application, there is a special configuration for the cat's eye. For optimum performance the spacing between the AOM and the lens should equal the lens focal length.

In this special configuration, as Ω is swept, the diffracted angle is also swept, but when a lens is positioned a distance f from the AOM, the rays emanate from the focal point of the lens, and will emerge from the lens parallel to the zeroth-order beam but deviated by a distance $d=f \cdot \tan(2\theta_{B,\text{ext}})$ for the first-order beam. The angle of the mirror can then be adjusted such that the diffracted beam always hits the mirror at normal incidence independent of the modulation frequency.

The advantage of this special cat's eye configuration was previously noted by Boulanger and co-workers,⁸ but they did not present their system in detail and did not inject the modulated light directly into a single-mode fiber. For our apparatus, we find that to efficiently inject the double-passed light into the single-mode fiber without adjusting the beam collimation it is important that the lens-mirror spacing be close to the lens focal length so that the beam is collimated at the single-mode fiber coupler. To be able to optimize the fiber coupling efficiency, we mounted the cat's eye lens onto a commercial lens focusing mount to optimize the lens-mirror spacing. Also, in contrast to what is implied by Boulanger and colleagues, the plano-convex lens in the cat's eye does not have to have a long focal length. We used a lens of half the focal length that they recommended, found the alignment sensitivity unchanged, and still achieved an ideal double-pass bandwidth.

Our cat's eye retroreflector consists of a +75 mm focal-length plano-convex lens and a planar mirror. While the performance is ideal when the AOM-lens spacing is equal to the focal length, the spacing can be reduced by shifting the whole cat's eye toward the AOM to save space if desired without dramatically compromising the bandwidth. For our system, we observe only a 20% decrease in the frequency-

tuning bandwidth when the AOM-cat's eye spacing is reduced by 50%.

In principle, we can achieve the same results using a spherical mirror in place of the cat's eye with the mirror positioned at a distance equal to its radius of curvature from the AOM. However, it is difficult to find an inexpensive concave mirror with high reflectivity.

C. Performance

We typically observe diffraction efficiencies exceeding 85% for a single pass and 75% for a double pass. Our data are summarized in Fig. 4, where single-pass and double-pass measurements are presented. It is important to emphasize that, although the three curves are shown together, the values of the diffraction efficiencies for the single-pass and double-pass measurements have different meanings. Whereas the single-pass diffracted power was directly measured at the AOM output, the light for the double-pass measurements was diffracted by the AOM twice and this diffracted light was injected into a single-mode fiber. The power for the double-pass configurations was measured at the fiber output to have a sensitive test on beam pointing stability. The peak double-pass diffraction efficiency is 75% without the fiber, as compared to 85% for the single-pass efficiency, which indicates that the second-pass diffraction through the AOM can be more efficient than the first-pass diffraction.

The ratio of the FWHM for the single-pass measurements to the FWHM of the double-pass measurements performed with the cat's eye is 1.43(7). Assuming a Gaussian dependence, in theory, if there were no beam steering problems, the diffraction efficiency function for the single-pass should be squared to arrive at the diffraction efficiency frequency dependence for the double pass. Thus, the optimum ratio of the single-pass to the double-pass bandwidth is $\sqrt{2}$, in agreement with what we observe. This agreement shows that the cat's eye acts as an ideal retroreflector. Note that the frequency dependencies shown in Fig. 4 are in terms of the rf modulation frequency, not the actual frequency shift, which is two times larger for a double pass. So, in terms of the frequency tuning range, the bandwidth of the double pass with the cat's eye is larger than that of the single pass by $\sqrt{2}$.

Quite striking is the difference between the double-pass measurements performed with and without the cat's eye lens. Having the cat's eye retroreflector as opposed to a planar mirror increases the frequency-tuning bandwidth of the device by about one order of magnitude.

ACKNOWLEDGMENTS

The authors thank M. A. Lombardi and D. R. Smith for helpful comments on this manuscript. This work is a contribution of NIST, an agency of the U.S. government, and is not subject to copyright.

¹C. Solomon, J. Dalibard, W. D. Phillips, A. Clairon, and S. Guellati, *Europhys. Lett.* **12**, 683 (1990).

²C. G. Townsend *et al.*, *Phys. Rev. A* **52**, 1423 (1995).

³S. R. Jefferts *et al.*, *Metrologia* **39**, 321 (2002).

⁴R. Guenther, *Modern Optics* (Wiley, New York, 1990).

⁵E. H. Young and S. Yao, *Proc. IEEE* **69**, 54 (1981).

⁶Products or companies named here are cited only in the interest of complete scientific description, and neither constitute nor imply endorsement by NIST or by the US government. Other products may be found to serve just as well.

⁷J. J. Snyder, *Appl. Opt.* **14**, 1825 (1975).

⁸J. S. Boulanger, M. C. Gagné, and R. J. Douglas., Proceedings of the 11th European Frequency and Time Forum, Neuchatel, France, 4–6, March, 1997, pp. 567–571.

MT80

AO Modulator/Shifter


TeO2 modulator for 450-1100 nm lasers

- Wide spectral range
- Linear or random Polar

These modulators have been specially designed for general purpose low speed applications. They are proposed with various wavelength ranges or with V-coating.

They can also be used as fixed frequency shifters @80 MHz, as well as variable frequency shifters with a frequency range up to 80 +/- 15 MHz.

With an adapted frequency range, user will be able to operate this device as a high speed low resolution deflector.

Specifications

Material-Acoustic mode	TeO2 [L]
Acoustic Velocity	V=4200 m/s
Optical Wavelength range	VIS: 450-700 nm IR: 700-1100 nm 1064 nm: 980-1100 nm
Transmission	VIS/IR: > 95 % 1064 nm: 98% @1064nm
Optical Input / Output polarizations	Linear
Aperture	1 x 2 mm ² or 1.5 x 2 mm ²
Carrier frequency / Frequency shift	80 MHz
Separation angle	10.1 mrd @532 nm 20.3 mrd @1064 nm
Diffraction efficiency (with TEM00 beam, M² ≤ 1.1)	85 %
Rise time	160 ns/mm (min 24 ns)
Amplitude modulation bandwidth	> 20 MHz (-3 dB, @150µm)
Static extinction ratio	> 2000/1
Max optical power density	VIS : 5 W / mm ² @532 nm, IR : > 10 W / mm ²
Input impedance	Nom 50 Ω
V.S.W.R.	Nom < 1.5/1
RF Power	VIS: ≤ 1 Watts IR/1064 nm: ≤ 2 Watts
Connector	SMA
Size / Weight	(LxIxh) 50.9 x 22.4 x 15.8 mm ³ / 50 g
Operating Temperature	10 to 40 °C


Relative Diffraction Efficiency vs RF Power

→ Separation angle ($\Delta\theta$) is wavelength (λ) sensitive:

$$\Delta\theta = \frac{\lambda F}{V}$$

→ RF power (P) is wavelength (λ) sensitive:

$$\frac{P_1}{P_2} = \frac{\lambda_1^2}{\lambda_2^2}$$

OPTION



Frequency range 80+/-15MHz
 Nominal efficiency over 80+/-15MHz > 70%

MT80-Ax-zz

X (aperture, mm) = **1 or 1.5**
Y = **frequency range** (MHz) if any
ZZ = **VIS** (450-700 nm), **IR** (700-1100 nm), **@1064nm**

Outline Drawing


(sizes in mm)

DRFA10Y

Fast VCO driver

Fast sweeping time VCO

- Random access or raster scan • OEM & Laboratory versions • Sweeping time $\leq 1 \mu\text{s}$

These Voltage Controlled Oscillators are the simplest way to control a deflector or a variable frequency shifter, in random access or raster scan mode.

A voltage between 0 and 10 V is applied on the frequency control input in order to sweep the frequency. The sweeping time from Fmin to Fmax can be as low as $1 \mu\text{s}$. The typical stability of a VCO is 50 to 100 KHz/°C : in option higher stabilities are proposed, and for the sharpest applications, a DDS will be preferred. These drivers are used in combination with AA amplifiers.

Specifications

Frequency range	Up to one octave (40-100, 60-150, 80-200, 140-300, 190-350) Matched to AO device at factory
Frequency stability	Nom +/- 50 to 100kHz/°C
Sweeping time	$\leq 1 \mu\text{s}$
Frequency control ¹	Analog 0-10 V / 1 K Ω
Linearity (full bandwidth)	Nom +/- 5 %
Rise time / Fall time (10-90 %)	< 10 ns
Modulation input control ¹ Video In	Analog 0-5 V / 50 Ω
Extinction ratio ²	> 40 dB
Harmonics	Nom > 20 dBc
Output RF power ³	Nom - 30 to 0 dBm (to be associated with AA Amplifier)
Output impedance	50 Ω
V.S.W.R.	$\leq 1.5 : 1$
Power supply	OEM version : 24 VDC – nom 150 mA Laboratory version ⁴ : 110-230 VAC – 50-60 Hz
Input / Output connectors ²	SMA
Size	OEM version : 84 x 69 x 15.8 mm ³ Laboratory version ⁴ : 310 x 250 x 105 mm ³
Weight	OEM version : 0.2 Kg Laboratory version ⁴ : 4 Kg
Cooling	Conduction through baseplate
Maximum case temperature	50 °C
Operating temperature	10 to 40 °C

¹ On request different Video In are proposed : 0-5 V/50 Ohms, 0-5V/500 Ohms, TTL...

² Other values or models on request

³ Fully adjustable with screw potentiometer from 0 to 1 watt or 0 to 2 Watts


⁴ These versions are complete turn key systems with minimum set up time.

The source and convenient amplifiers are integrated in a rack which is supplied with 110-230 VAC.


A selection switch on front panel allows the user to select both operating modes:

- CW mode: internal CW modulation of RF power with front panel cursor.
- AM mode: external amplitude modulation controlled through external modulation input.
- F CW mode: internal CW modulation of RF Frequency with front panel cursor.
- FM mode: external frequency modulation controlled through external input.

Relative Power/Amplitude vs MOD IN

Frequency versus Input voltage control


DRFA10Y-B-0

B = B (24 VDC), **D** (110-230 VAC) (Power supply)
O = 0 dBm (30, 33...)

Outline Drawing


sizes in mm OEM Version

AMPA-B-30 -33

Large Band RF amplifier

RF amplifiers 1 & 2 Watts

• Large band • General purposes

These amplifiers are made for general purpose applications, covering frequency ranges up to 450-600 MHz.

In association with a VCO driver or a DDS driver, they will provide the necessary RF power to drive an acousto-optic device up to 2 watts.

Specifications

Frequency range¹	1 Watt : 20-450 MHz 2 watts : 20-600 MHz
Gain²	1 Watt : ≥ 33 dB 2 watts : ≥ 40 dB
Gain Flatness	Nom +/- 0.5 dB, ≤ +/1 dB
Noise Figure	1 Watt : nom 5 dB 2 watts : nom 7 dB
Output RF Power (1 dB compression)	≥ 30 dBm (≥ 29.5 dBm @ <40 MHz), 1 Watt ≥ 33 dBm , 2 Watts
Output Impedance	50 Ω
CLASS	A
Power supply³	1 Watt : 24 +/- 0.5 VDC - ≤ 340 mA 2 watts : 24 +/- 0.5 VDC - ≤ 500 mA
Input / Output connectors	SMA female
Size	76 x 40 x 42 mm ³
Weight	0.22 Kg
Heat exchange	Conduction through baseplate
Operating temperature	-10 to +55 °C

¹ Optimized at factory according to AO device if any

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

Item no. 19 09 74 **1/3" black/white**
Item no. 19 09 88 **1/3" colour**

Intended use

The card camera with manual focus adjustment is designed for surveillance and security purposes in uncontrolled and critical areas (such as shop, door intercoms, entrance areas or parking lots etc.). The recorded images can be played back with sound on any suitable display device with video input (such as TV or monitor).

The card camera is not permitted to use without protection (without the corresponding casing). It must be power by 12 V= only. Contact with moisture must be avoided at all times.

Any usage other than described above is not permitted and can damage the product and lead to associated risks such as short-circuit, fire, electric shock, etc. No part of the product may be modified or rebuilt. Please read the operating instructions thoroughly and keep the operating instructions for further reference.

Safety instructions

! We do not resume liability for resulting damages to property or personal injury if the product has been abused in any way or damaged by improper use or failure to observe these operating instructions. The guarantee will then expire!

An exclamation mark in a triangle indicates important information in the operating instructions. Carefully read the whole operating instructions before operating the device, otherwise there is risk of danger.

- Unauthorised conversion and/or modification of the device are inadmissible because of safety and approval reasons (CE).
- The card camera must not be exposed to dust, extreme temperatures, direct sunlight, intense vibration or dampness.
- The card camera must be installed in proper housing in prior to operation.
- The card camera should not be installed in vicinity of strong magnetic or electric fields, e.g. mobile phones, walkie-talkies, electric motors etc.
- The card camera is only permitted to be powered by power pack or a battery/accumulator with screened and stabilized direct voltage. The voltage source must be able to supply sufficient power.
- After rapid temperature changes, the device needs approximately 15 minutes of stabilization to accommodate to the new ambient temperature before use.
- When connecting the camera, make sure that the connection cables are not damaged by sharp edges.
- When using other devices or tools to assist the installation, the corresponding operating instructions must be observed.
- When using a power pack as the voltage source, valid safety regulations must be followed under all circumstances. Car chargers and toy transformers are not suitable as voltage sources and will lead to damage to the components or to the device failing to function.
- In commercial institutions, the accident prevention regulations of the relevant professional insurance association for electrical systems and operating materials have to be observed.
- If safe operation of the card camera is no longer possible, disconnect the appliance immediately and secure it against inadvertent operation. Safe operation is no longer possible if:
 - the device shows visible damages,
 - the device no longer works and
 - the device was stored under unfavourable conditions for a long period of time,
 - the device was subject to considerable transport stress.
- If you have any doubt about the card camera installation or operation, consult adequate qualified personnel.
- Packaging materials must not be lying around carelessly. They could become dangerous toys in hands of children.
- This card camera is not a toy and should be kept out of reach of children!
- Servicing, adjustment or repair works must only be carried out by a specialist/ specialist workshop.
- If any questions arise that are not answered in this operating instruction, please contact our Technical Advisory Service or other experts.

Connection and operation

! Use 12V= voltage source to power the card camera. Never try to operate the camera with another voltage source. The DC source for the camera has to supply a current of 100 mA as minium load capacity. Never overload the voltage source.

If you use a power pack as voltage source, this must be turned off during connecting. Observe the polarity when connecting the voltage supply. If the polarity is not correct, the components of the card camera may be destroyed.

Camera connections

Connect the small white plug of the connection cable with the camera module. This only fits into the socket with the right polarity.

Video connection

Connect the video input of your TV or monitor with the yellow cinch socket of the cameramodule.

 Use a suitable optional cinch connection cable (polarity: inside video, outside ground). For any extensions, only use suitably shielded cinch cables. The use of other, unsuitable cables may lead to interference. Keep the length of the cable as short as possible.

Power supply connection

Connect the red power supply plug of the camera connection cable (VDC) with a matching power pack. The DC low voltage plug must have the following properties:

- Dimensions outside 5.5 mm, inner hole 2.1 mm.
- Polarity inside plus (+), outside minus (-).

Manual focal distance adjustment

- Loosen the screw on the len mount.
- Turn the len clockwise/ anticlockwise to adjust the focal distance.
- Tighten the screw after the desired focal distance is set.

Maintenance and cleaning

Apart from sporadically cleaning the lens, the card camera is maintenance-free.

Use a clean, lint-free, antistatic and dry cloth to clean the len. Do not use any abrasive or chemical agents or detergents containing solvents.

Disposal

 In order to preserve, protect and improve the quality of environment, protect human health and utilise natural resources prudently and rationally, the user should return unserviceable product to relevant facilities in accordance with statutory regulations.

The crossed-out wheeled bin indicates the product needs to be disposed separately and not as municipal waste.

Technical data

Item no.	190974	190988
Operating voltage	12V=	12V=
Current consumption	approx. 100mA	approx. 100mA
Image sensor	1/3" CCD black/white	1/3" CCD colour
Video format	PAL	PAL
Lens (focal distance)	1.5m	1.5m
Video level	1.0Vp-p ± 0.2	1.0Vp-p ± 0.2
Light intensity	0.01 Lux at F1.2	0.5 Lux at F1.2
Shutter speed	1/50 ~ 1/100000 second	1/50 ~ 1/100000 second
Effective pixels	PAL: 512 x 582	PAL: 512 x 582
Resolution (TV lines)	380	420
Video signal-to-noise ratio	> 48dB	> 48dB
Operating temperature	-20°C to +55°C	-20°C to +55°C
Operating humidity	< 90% (non-condensing)	< 90% (non-condensing)
Dimensions (LxWxH)	32 x 32 x 27 mm	32 x 32 x 27 mm
Fastening/ hole pitch	27mm (ø2.2mm)	27mm (ø2.2mm)
Weight (without conection cable)	13.5g	13.5g

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PDA10A Operating Manual - Wideband Amplified Silicon Detector

Description:

The PDA10A is a wideband amplified, silicon detector designed for detection of light signals from DC to 150 MHz. A buffered output drives a 50Ω input impedance up to 5 volts. The PDA10A housing includes a removable threaded coupler that is compatible with any number of Thorlabs 1" and ½" threaded accessories. This allows convenient mounting of external optics, light filters, and apertures, as well as providing an easy mounting mechanism using the Thorlabs cage assembly accessories.

The PDA10A has two 8-32 tapped mounting hole with a 0.25" mounting depth and includes a 120VAC AC/DC power supply. The PDA10A-EC has two M4 tapped mounting holes and includes a 230VAC AC/DC power supply.

Specifications:

Electrical:	
Detector	Silicon
Active Area	0.8mm ² (∅1.0mm)
Response	200 to 1100 nm
Peak Response (typ)	0.45 A/W (750nm)
Small Signal Bandwidth ¹	150MHz (min.)
NEP (960 nm)	5.5x10 ⁻¹¹ W/√Hz (max.)
Noise (RMS)	2.0mV (max.)
Dark Offset	20mV (max.)
Output Voltage (50Ω)	0 to 5V
Output voltage	0 to 10V
Transimpedance Gain	
Hi-Z	1 x 10 ⁴ V/A
50Ω	5 x 10 ³ V/A

General:	
On / Off Switch	Slide
Output	BNC
Optical Head Size	2.8"x1.9" x 0.83" 70mm x 48mm x 21mm
Weight	60 grams
Accessories	SM1T1 Coupler SM1RR Retainer Ring
Storage Temp	-25 to 70°C
Operating Temp	10 to 50°C
AC Power Supply	AC - DC Converter
Input Power	100-120VAC (220-240VAC -EC version) 50-60Hz, 5W

1. The small signal bandwidth was measured with output amplitude of 200mV and a dc offset of 200mV, driving a 50Ω load termination.
2. All measurements performed with 50Ω load unless stated otherwise.


Setup

- Unpack the optical head, install a Thorlabs TR-series ½" diameter post into one of the 8-32 (M4 on -EC version) tapped holes, located on the bottom and side of the sensor, and mount into a PH-series post holder.
- Connect the power supply 3-pin plug into the mating receptacle on the PDA10A.
- Plug the power supply into a 50-60Hz, 100-120VAC outlet (220-240VAC for -EC version).
- Attach a 50Ω coax cable (i.e. RG-58U) to the output of the PDA. When running cable lengths longer than 12" we recommend terminating the opposite end of the coax with a 50Ω resistor (Thorlabs p/n T4119) for maximum performance.

Operation

- The PDA10A is switched on by the 'POWER' Slide switch, located on the top of the optical sensor.
- The light to voltage conversion can be estimated by factoring the wavelength-dependent responsivity of the Silicon detector with the transimpedance gain as shown below:
 [E.g. output in volts / watt = transimpedance gain (V/A) x responsivity (A/W)]
- The maximum output of the PDA10A is 10 volts for high impedance loads (5V for 50Ω loads). The output signal should be below the maximum output voltage to avoid saturation. If necessary, use external neutral density filters to reduce the input light level.
- For maximum linearity performance when measuring focused beams, fiber outputs, or small diameter beams, do not exceed a maximum intensity of 10mW/cm².
- Because of the finite gain-bandwidth performance common to all amplifier circuits, the bandwidth of the PDA10A decreases with increased output signal levels.

Figure 1 - PDA10A Spectral Responsivity Curve

Fiber Adapters and Other Accessories

Thorlabs sells a number of accessories that are compatible with the 1" and ½" threads on the PDA housing including FC, SMA, and ST fiber adapters, stackable lens tubes for mounting optics, and cage assemblies that allow the PDA to be incorporated into elaborate 3-D optical assemblies.

Caution: The PDA10A was designed to allow maximum accessibility to the photodetector by having the front surface of the diode flush with the outside of the PDA housing. When using fiber adapters, make sure that the fiber ferrule does not crash into the detector. Failure to do so may cause damage to the diode and / or the fiber. An easy way to accomplish this is to install a SM1RR retaining ring (included with the PDA10A) inside the 1" threaded coupler *before* installing the fiber adapter.

Also available in the PDA series are InGaAs, switchable gain InGaAs and switchable gain silicon models.

Maintaining the PDA10A

There are no serviceable parts in the PDA10A optical head or power supply. The housing may be cleaned by wiping with a soft damp cloth. The window of the detector should only be cleaned using optical grade wipes. If you suspect a problem with your PDA10A please call Thorlabs and an engineer will be happy to assist you.

Contact

Americas
Thorlabs Inc.
435 Route 206 North
Newton NJ 07860
USA
Ph: (973) 579-7227
Fax: (973) 300-3600
www.thorlabs.com
Email:
techsupport@thorlabs.com

Europe
Thorlabs GmbH
Gaustr. 11
85757 Karlsfeld
Germany
Ph: +49 (0) 8131-59-56-0
Fax: +49 (0) 8131-59-56-99
www.thorlabs.com
Email:
Europe@thorlabs.com

UK and Ireland
Thorlabs, LTD.
1 Saint Thomas Place, Ely
Cambridgeshire CB7 4EX
Great Britton
Ph: +44 (0) 1353-654440
Fax: +44 (0) 1353-654444
www.thorlabs.com
Email:
sales.uk@thorlabs.com

Scandinavia
Thorlabs Sweden AB
Box 141 94
400 20 Göteborg
Sweden
Ph: +46-31-733-30-00
Fax: +46-31-703-40-45
www.thorlabs.com
Email:
Scandinavia@thorlabs.com

Japan
Thorlabs Japan, Inc
5-17-1, Ohtsuka
Bunkyo-ku, Tokyo 112-0012
Japan
Ph: +81-3-5977-8401
Fax: +81-3-5977-8402
www.thorlabs.jp
Email:
sales@thorlabs.jp

WEEE

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

This offer is valid for Thorlabs electrical and electronic equipment

- sold after August 13th 2005
- marked correspondingly with the crossed out “wheelie bin” logo (see fig. 1)
- sold to a company or institute within the EC
- currently owned by a company or institute within the EC
- still complete, not disassembled and not contaminated

As the WEEE directive applies to self contained operational electrical and electronic products, this “end of life” take back service does not refer to other Thorlabs products, such as

- pure OEM products, that means assemblies to be built into a unit by the user (e. g. OEM laser driver cards)
- components
- mechanics and optics
- left over parts of units disassembled by the user (PCB's, housings etc.).

If you wish to return a Thorlabs unit for waste recovery, please contact Thorlabs or your nearest dealer for further information.

Waste treatment on your own responsibility


If you do not return an “end of life” unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

Ecological background

It is well known that WEEE pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS directive is to reduce the content of toxic substances in electronic products in the future. The intent of the WEEE directive is to enforce the recycling of WEEE. A controlled recycling of end of live products will thereby avoid negative impacts on the environment.

Crossed out “wheelie bin” symbol

TOLERANCES	
UNLESS OTHERWISE SPECIFIED:	
DIMENSIONS ARE IN INCHES	
LINEAR TOLERANCES:	
TWO PLACE DECIMAL: ±0.010	
THREE PLACE DECIMAL: ±0.005	
ANGULAR: ±30'	
SURFACE FINISH: 32 MICROINCHES	
PARALLELISM: 0.002	
FLATNESS: 0.002	
STRAIGHTNESS: 0.002	
CONCENTRICITY: 0.002	
PERPENDICULARITY: 0.002	
THREAD: CLASS 2 FIT	

NAME	DATE
DRAWN EC	9/26/05
ENG APPR. EC	9/26/05
MFG APPR. EC	9/26/05
PROPRIETARY AND CONFIDENTIAL	
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF THORLABS, INC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF THORLABS, INC. IS PROHIBITED.	


THORLABS INC.		PO BOX 366 NEWTON NJ	
TITLE: WIDEBAND Si AMPLIFIED DETECTOR			
MATERIAL: PDA10A		SIZE A	REV. A
SCALE: 1:1		SHEET 1 OF 1	
DWG. NO. 13054-E01	PART NO. PDA10A		

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

General information concerning the CE marking

HAMEG instruments fulfill the regulations of the EMC directive. The conformity test made by HAMEG is based on the actual generic- and product standards. In cases where different limit values are applicable, HAMEG applies the severer standard. For emission the limits for residential, commercial and light industry are applied. Regarding the immunity (susceptibility) the limits for industrial environment have been used.

The measuring- and data lines of the instrument have much influence on emission and immunity and therefore on meeting the acceptance limits. For different applications the lines and/or cables used may be different. For measurement operation the following hints and conditions regarding emission and immunity should be observed:

1. Data cables

For the connection between instruments resp. their interfaces and external devices, (computer, printer etc.) sufficiently screened cables must be used. Without a special instruction in the manual for a reduced cable length, the maximum cable length of a dataline must be less than 3 meters and not be used outside buildings. If an interface has several connectors only one connector must have a connection to a cable.

Basically interconnections must have a double screening. For IEEE-bus purposes the double screened cable HZ72 from HAMEG is suitable.

2. Signal cables

Basically test leads for signal interconnection between test point and instrument should be as short as possible. Without instruction in the manual for a shorter length, signal lines must be less than 3 meters and not be used outside buildings.

Signal lines must be screened (coaxial cable - RG58/U). A proper ground connection is required. In combination with signal generators double screened cables (RG223/U, RG214/U) must be used.

3. Influence on measuring instruments

Under the presence of strong high frequency electric or magnetic fields, even with careful setup of the measuring equipment an influence of such signals is unavoidable.

This will not cause damage or put the instrument out of operation. Small deviations of the measuring value (reading) exceeding the instruments specifications may result from such conditions in individual cases.


4. Noise immunity of spectrum analyzers

In the presence of strong electric or magnetic fields it is possible that they may become visible together with the signal to be measured. The methods of intrusion are many: via the mains, via the signal leads, via control or interface leads or by direct radiation. Although the spectrum analyzer has a metal housing there is the large CRT opening in the front panel where it is vulnerable. Parasitic signals may, however, also intrude into the measuring object itself and from there propagate into the spectrum analyzer.


HAMEG Instruments GmbH

Deutsch	3		
English			
General information concerning the CE marking	24		
Spectrum Analyzer: Series HMS	26		
Specifications	27		
1 Installation and safety instructions	28		
1.1 Setting up the instrument	28	9.6 Transient Limiter HZ560	41
1.2 Safety	28	9.7 75/50- Ω -Converter HZ575	41
1.3 Correct operation	28		
1.4 Ambient conditions	28	10 Front panel connections	41
1.5 Warranty and repair	29	10.1 USB connector	41
1.6 Maintenance	29	10.2 PHONE	41
1.7 CAT I	29	10.3 PROBE POWER	41
1.8 Mains voltage	29	10.4 EXTERNAL TRIGGER	41
		10.5 OUTPUT 50 Ω (Tracking Generator)	41
		10.6 INPUT 50 Ω	41
2 Controls and display	30	11 Rear panel connections	42
3 Short description of the series HMS	32	11.1 USB connector	42
		11.2 DVI connector	42
4 Setting of parameters	32	11.3 REF IN / REF OUT	42
4.1 Numerical keyboard	32		
4.2 Knob	32	12 Remote Control	42
4.3 Arrow buttons	32		
4.4 Interactive softkeys	33	13 Appendix	42
4.5 How to enter numerical values	33		
5 Instrument functions	33		
5.1 Setting of the frequency (FREQ)	33		
5.2 Frequency range displayed (SPAN)	33		
5.3 Setting of the amplitude parameters (AMPL)	33		
5.4 Setting of the bandwidth (BANDW)	33		
5.5 Setting of the SWEEP	34		
5.6 Curve display settings (TRACE)	34		
5.7 The use of markers	35		
5.8 Peak-Search	35		
5.9 Operation in the Receiver-Mode	35		
6 Store and recall instrument settings	36		
6.1 Instrument settings	36		
6.2 Reference signals	36		
6.3 Waveforms	37		
6.4 Screenshots	37		
7 Extended operating modes	39		
7.1 Using the help function	39		
7.2 Display settings	39		
8 General instrument settings	39		
8.1 Language settings	39		
8.2 Basic settings	39		
8.3 Interface settings	40		
8.4 Printer settings	40		
8.5 Update (Firmware / Hilfe)	40		
9 Optional Accessories	40		
9.1 Activation of the Preamplifier H03011	40		
9.2 19" Rack mount kit 4HE HZ46	40		
9.3 Carrying case HZ99	40		
9.4 Near field probe HZ530/HZ540	41		
9.5 Measurements of spectra with a HMS3010 and a VSWR bridge HZ547	41		

1GHz [3GHz] Spectrum Analyzer HMS1000 [HMS3000]

3GHz Spectrum Analyzer
HMS3000 without TG

3GHz EMI Near Field Probe
Set HZ550L

VSWR Test Unit HZ547

- ✓ Frequency range 100kHz...1 GHz [3 GHz]
- ✓ Amplitude measurement range -114dBm...+20dBm
DANL - 125 dBm [DANL - 135 dBm] with Preamp. Option H03011
- ✓ Sweep time 20ms...1000s
- ✓ Resolution bandwidth 1 kHz [100 Hz]...1MHz in 1-3 steps,
200kHz (-3dB) additional [200 Hz], 9kHz, 120kHz, 1MHz (-6dB)
- ✓ Spectral purity < -100dBc/Hz (@100kHz)
- ✓ Video bandwidth 10Hz...1MHz in 1-3 steps
- ✓ Tracking Generator HMS1010 [HMS3010] -20dBm/0dBm
- ✓ Integrated AM and FM demodulator (int. speaker)
- ✓ Detectors: Auto-, min-, max-peak, sample, RMS, quasi-peak
- ✓ 8 Marker with delta marker, miscellaneous peak functions
- ✓ Crisp 16.5cm (6.5") TFT VGA display, LED backlight, DVI output
- ✓ 3 x USB for mass-storage, printer and remote control
optional IEEE-488 (GPIB) or Ethernet/USB Interface

**1GHz Spectrum Analyzer HMS1000, HMS1010 (with TG)
3GHz Spectrum Analyzer HMS3000, HMS3010 (with TG)**
All data valid at 23 °C after 30 minute warm-up

Frequency	
Frequency range:	
HMS1000, HMS1010	100kHz...1GHz
HMS3000, HMS3010	100kHz...3GHz
Temperature stability:	± 2ppm (0...30°C)
Aging:	± 1 ppm/year
Frequency counter (from SW 2.0):	
Resolution	1Hz
Accuracy	± (Frequency x tolerance of reference)
Span setting range:	
HMS1000, HMS1010	0Hz (zero span) and 1kHz...1GHz
HMS3000, HMS3010	0Hz (zero span) and 100Hz...3GHz
Spectral purity, SSB phase noise:	
30kHz from carrier (500MHz, +20°C...30°C)	< -85dBc/Hz
100kHz from carrier (500MHz, +20°C...30°C)	< -100dBc/Hz
1MHz from carrier (500MHz, +20°C...30°C)	< -120dBc/Hz
Sweep time:	
Span = 0Hz	20ms...100s
Span > 0Hz	20ms...1000s, min. 20ms/600MHz
Resolution bandwidths (-3 dB):	
HMS1000, HMS1010	1kHz...1MHz in 1-3 steps, 200kHz
HMS3000, HMS3010	100Hz...1MHz in 1-3 steps, 200kHz
Tolerance:	
≤ 300kHz	± 5% typ.
1MHz	± 10% typ.
Resolution bandwidths (-6 dB):	
HMS1000, HMS1010	9kHz, 120kHz, 1MHz
HMS3000, HMS3010	200Hz, 9kHz, 120kHz, 1MHz
Video bandwidths:	10Hz...1MHz in 1-3 steps
Amplitude	
Display range:	Average noise level displayed up to +20 dBm
Amplitude measurement range:	Typ. -114dBm...+20dBm
Max. permissible DC at HF input:	80V
Max. power at HF input:	20dBm, 30dBm for max. 3 Min.
Intermodulation free range:	
IM3 products, 2 x -20dBm (-10dBm ref. level)	66dB typ. (typ. +13dBm third-order intercept)
[at distance between signals ≤ 2MHz]	60dB typ. (+10dBm TOI)
[at distance between signals > 2MHz]	66dB typ. (typ. +13dBm TOI)
DANL (Displayed average noise level):	
(RBW 1kHz, VBW 10Hz, ref. level ≤ -30dBm 10MHz...1GHz resp. 3GHz)	-105dBm, typ. -114dBm
With Preamp.	-135dBm typ. (100Hz RBW)
Inherent spurious:	
[ref. level ≤ -20dBm, f < 30MHz, RBW ≤ 100kHz]	< -80dBm
Input related spurious:	
[Mixer level ≤ -40dBm, carrier offset > 1MHz]	-70dBc typ., -55dBc (2...3GHz)
2nd harmonic receive frequency	
[mixer level -40dBm]:	-60dBc typ.
Level display:	
Reference level	-80dBm...+20dBm in 1dB steps
Display range	100dB, 50dB, 20dB, 10dB, linear
Logarithmic display scaling	dBm, dBμV, dBmV
Linear display scaling	μV, mV, V, nW, μW, mW, W
Measured curves:	1 curve and 1 memory curve
Trace mathematics:	A-B (curve-stored curve), B-A
Detectors:	Auto-, Min-, Max-Peak, Sample, RMS, Average, Quasi-Peak
Failure of level display:	< 1,5dB, typ. 0,5dB
[ref. level to ref. level-50dB, 20°C...30°C]	

Marker / Deltamarker	
Number of marker:	8
Marker functions:	Peak, next peak, minimum, center = marker, frequency, reference level = marker level, all marker on peak
Marker displays:	Normal (level), noise marker, (frequency) counter (from SW 2.0)

Inputs / Outputs	
HF Input	N socket
Input Impedance:	50Ω
VSWR (10MHz...1GHz/3GHz):	< 1,5 typ.
Output tracking generator:	
(HMS1010/HMS3010)	N socket
Output Impedance:	50Ω
Frequency range:	5MHz...1GHz/3GHz
Output level:	-20dBm/0dBm
Trigger and external reference input:	BNC female, selectable
Trigger voltage	TTL
Reference frequency	10MHz
Essential level (50Ω)	10dBm
Supply output for field probes:	6VDC, max. 100mA (2,5mm DIN jack)
Audio output (Phone):	3,5mm DIN jack
Demodulation	AM and FM (internal speaker)

Miscellaneous	
Display:	6,5" TFT Color VGA Display
Save / Recall memory:	10 complete device settings
Trigger:	Free run, Video Trigger (from SW 2.0), external Trigger
Interfaces:	Dual-Interface USB/RS-232 (HO720), USB-Stick (frontside), USB-Printer (rear side) from SW 2.0, DVI-D for ext. monitor
Power supply:	105...253 V, 50/60 Hz, CAT II
Power consumption:	Max. 40Watt at 230V, 50 Hz
Protection class:	Safety class I (EN61010-1)
Operating temperature:	+5°C...+40°C
Storage temperature:	-20°C...+70°C
Rel. humidity:	5%...80% (non condensing)
Dimensions (W x H x D):	285 x 175 x 220 mm
Weight:	3.6 kg

Accessories supplied:
Line cord, Operating manual, Dual-Interface USB/RS-232 (HO720), CD, HZ21 Adapter plug (N plug to BNC socket)
Optional accessories:
HO3011 Preamplifier -135dBm DANL (100Hz RBW)
HO730 Dual-Interface Ethernet/USB
HO740 Interface IEEE-488 (GPIB), galvanically isolated
HZ547 3GHz VSWR Test Unit for HMS1010, HMS3010
HZ520 Plug-in Antenna with BNC connection
HZ530 Near-Field Probe Set 1GHz for EMV diagnostics
HZ540 Near-Field Probe Set 3GHz for EMV diagnostics
HZ560 Transient limiter
HZ575 75/50Ω Converter
HZ46 4RU 19" Rackmount Kit

1 Installation and safety instructions

1.1 Setting up the instrument


As can be seen from the figures, the handle can be set into different positions:

- A and B = carrying
- C = horizontal operating
- D and E = operating at different angles
- F = handle removal
- G = operating using the feet's, batch use and for shipping in original packaging

Attention!

When changing the handle position, the instrument must be placed so that it cannot fall (e.g. placed on a table). Then the handle locking knobs must be simultaneously pulled outwards and rotated to the required position. Without pulling the locking knobs they will latch in into the next locking position.

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Removal/fitting of the handle: The handle can be removed in position F, pulling the side parts outside the housing. Adding the handle works vice versa.

1.2 Safety

The instrument fulfils the VDE 0411 part 1 regulations for electrical measuring, control and laboratory instruments and was manufactured and tested accordingly. It left the factory in perfect safe condition. Hence it also corresponds to European Standard EN 61010-1 and International Standard IEC 1010-1. In order to maintain this condition and to ensure safe operation the user is required to observe the warnings and other directions for use in this manual. Housing, chassis as well as all measuring terminals are connected to safety ground of the mains. All accessible metal parts were tested against the mains with 2200 VDC. The instrument conforms to safety class I. The oscilloscope may only be operated from mains outlets with a safety ground connector. The mains plug has to be installed prior to connecting any signals. It is prohibited to separate the safety ground connection. If suspected that safe operation may not be guaranteed do not use the instrument any more and lock it away in a secure place.

Safe operation may be endangered if any of the following was noticed:

- in case of visible damage.
- in case loose parts were noticed
- if it does not function any more.
- after prolonged storage under unfavourable conditions (e.g. like in the open or in moist atmosphere).
- after any improper transport (e.g. insufficient packing not conforming to the minimum standards of post, rail or transport firm)

1.3 Correct operation

Please note: This instrument is only destined for use by personnel well instructed and familiar with the dangers of electrical measurements. For safety reasons the oscilloscope may only be operated from mains outlets with safety ground connector. It is prohibited to separate the safety ground connection. The plug must be inserted prior to connecting any signals. The oscilloscope is destined for operation in industrial, business, manufacturing, and domestic sites.

1.4 Ambient conditions

Operating ambient temperature: +5 °C to +40 °C. During transport or storage the temperature may be -20 °C to +70°C. Please note that after exposure to such temperatures or in case of condensation, proper time must be allowed until the instrument has reached the permissible temperature, and until the condensation has evaporated before it may be turned on! Ordinarily this will be the case after 2 hours. The oscilloscope is destined for use in clean and dry environments. Do not operate in dusty or chemically aggressive atmosphere or if there is danger of explosion. The any operating position may be used, however, sufficient ventilation must be ensured. Prolonged operation requires the horizontal or inclined position.

Do not obstruct the ventilation holes!

Specifications are valid after a 30 minute warm-up period at 23 degr. C (tolerance ± 2 degr. C). Specifications without tolerances are average values.

1.5 Warranty and repair

HAMEG instruments are subjected to a strict quality control. Prior to leaving the factory, each instrument is burnt in for 10 hours. By intermittent operation during this period almost all defects are detected. Following the burn in, each instrument is tested for function and quality, the specifications are checked in all operating modes; the test gear is calibrated to national standards.

The warranty standards applicable are those of the country in which the instrument was sold. Reclamations should be directed to the dealer.

Only valid in EU countries:

In order to speed claims, customers in EU countries may also contact HAMEG directly. Also, after the warranty expired, the HAMEG service will be at your disposal for any repairs.

Return Material Authorization (RMA):

Prior to returning an instrument to HAMEG, ask for a RMA number either by internet (<http://www.hameg.com>) or fax (+49 (0) 6182 800 501). If you do not have an original shipping carton, you may obtain one by calling the HAMEG service dept (+49 (0) 6182 800 500) or by sending an email to service@hameg.com.

1.6 Maintenance

Clean the outer case using a dust brush at regular intervals. Dirt can be removed from housing, handle, all metal and plastic parts using a cloth moistened with water and 1 % detergent. Greasy dirt may be removed with benzene (petroleum ether) or alcohol. Thereafter wipe the surfaces with a dry cloth. Plastic parts should be treated with a suitable antistatic solution. No fluid may enter the instrument. Do not use other cleansing agents as they may adversely affect the plastic or lacquered surfaces.

1.7 CAT I

This oscilloscope is destined for measurements in circuits not connected to the mains or only indirectly. Direct measurements, i.e. with a galvanic connection to circuits corresponding to the categories II, III, or IV are prohibited! The measuring circuits are considered not connected to the mains if a suitable isolation transformer fulfilling safety class II is used. Measurements on the mains are also possible if suitable probes like current probes are used which fulfill the safety class II. The measurement category of such probes must be checked and observed. The measurement categories were derived corresponding to the distance from the power station and the transients hence to be expected. Transients are short, very fast voltage or current excursions which may be periodic or not.

1.8 Mains voltage

The instrument has a wide range power supply from 105 to 253V, 50 or 60 Hz $\pm 10\%$. There is hence no line voltage selector. The line fuse is accessible on the rear panel and part of the line input connector. Prior to exchanging a fuse, the line cord must be pulled out. Exchange is only allowed if the fuse holder is undamaged. It can be taken out using a screwdriver put into the slot. The fuse can be pushed out of its holder and exchanged. The holder with the new fuse can then be pushed back in place against the spring. It is prohibited to "repair" blown fuses or to bridge the fuse. Any damages incurred by such measures will void the warranty.

Type of fuse:

Size 5 x 20 mm; 250V~, C; IEC 127, Bl. III;
DIN 41 662 (or DIN 41
571, Bl. 3). Cut off: slow blow (T) 2A.

2 Controls and display

Front panel


(HMS1010 differs in frequency range;
HMS3000 / HMS1000 without Tracking Generator)

- 1 **Display** (TFT)
6,5" VGA TFT Display
 - 2 **Interaktive Softkeys** (buttons)
Direct access of all relevant functions
 - 3 **POWER** (button)
Power switch turns the instrument on/off
- Area A:**
This area includes the parameter settings.
- 4 **AMPL** (illuminated button)
Setting of amplitude parameters
 - 5 **SPAN** (illuminated button)
Setting of the Span
 - 6 **FREQ** (illuminated button)
Setting of the frequency
 - 7 **TRACE** (illuminated button)
Configuration of data acquisition and analysis
 - 8 **SWEEP** (illuminated button)
Setting of the Sweep Time and the trigger source
 - 9 **BANDW** (illuminated button)
Setting of the resolution and video bandwidth
 - 10 **LINES** (illuminated button)
Configuration of displayed and limit lines
 - 11 **MEAS** (illuminated button)
Implementation of extended measurements
 - 12 **DISPLAY** (illuminated button)
Setting of the display
 - 13 **PEAK SEARCH** (illuminated button)
Measuring value peak display
 - 14 **MARKER →** (illuminated button)
Search function of marker
 - 15 **MARKER** (illuminated button)
Selection and arrangement of the absolute and relative marker
 - 16 **MODE** (illuminated button)
Switching between SWEEP- and RECEIVER-Mode
 - 17 **PRESET** (button)
Factory reset
 - 18 **AUTO TUNE** (button)
Automatically setting of instrument settings

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Area C (Variation):


This area includes the settings via rotary knob and arrow buttons

- 23 **Rotary knob**
Knob to adjust and activate the values or menu items by pushing
- 24 **Arrow buttons ▲ ▼** (buttons)
Set of signal parameters


Area D (General):

This area includes the general instrument settings

- 25 **FILE/PRINT** (illuminated button)
The key FILE/PRINT on the front panel allows you to store curves, screen shots, screen shots with settings
- 26 **SETUP** (illuminated button)
Display of general instrument settings
- 27 **HELP** (illuminated button)
Including display help
- 28 **SAVE/RECALL** (illuminated button)
Store and restore of instrument settings, references, curves and screenshots

29 **REMOTE** (illuminated button)
 Toggling between front panel and external operation

Area E:

This area includes a series of connectors.

30 **USB port**
 Front USB port for storing parameters

31 **PHONE**
 Headphone connector 3,5 mm jack;
 Impedance > 8 Ω

32 **PROBE POWER**
 Power supply (6 V_{DC}) for field probes
 (2,5 mm jack)

33 **External TRIGGER** (BNC socket)
 BNC input for external trigger signal

34 **OUTPUT 50 Ω**
 Tracking Generator (N connector)
 (HMS3000 / HMS1000 haven't got this connector)

29 **INPUT 50 Ω**
 Input N connector

Rear panel

36 **Mains input connector with fuse**


37 **Interface**
 H0720 Dual-Interface (USB/RS-232) is provided as standard

38 **DVI** (connector)
 Connection of external monitors and projectors

39 **USB port**

40 **REF IN** (BNC socket)
 Reference input

41 **REF OUT** (BNC socket)
 Reference output

3 Short description of the series HMS

Introducing the all new HMS series of 4 TFT spectrum analysers, HAMEG launches a family of instruments featuring frequency ranges up to 1GHz and 3GHz respectively, with or without tracking generator, suiting the budgets of different users.

HMS1000 and HMS3000 provide for frequency ranges between 100kHz and 1GHz/3GHz respectively. Both instruments are also available as HMS1010 and HMS3010 including tracking generators for four-terminal measurements. The amplitude measurement range extends from -114dBm to +20dBm (1kHz resolution bandwidth RBW). Using the optional pre-amplifier, the amplitude range can be extended to -135dBm (100Hz). Spectral purity exceeds -100dBc/Hz at 100kHz carrier frequency spacing.


Fig. 3.1: Menu of the frequency

Special emphasis has been put on a comprehensive and practice oriented set of filter bandwidths. The HMS3000 unit covers the range of 100 Hz to 1 MHz [-3 dB] in 1-3 steps. For pre-compliance measurements, bandwidths of 200 Hz, 9 kHz, 120 kHz and 1 MHz [-6 dB] are additionally provided. Moreover, detectors for auto-, min- and max-peak, sample, rms and quasi-peak as well as built-in AM and FM demodulators (headphone output and built in loudspeaker) are included in the standard package. Eight markers, including delta marker, as well as diverse peak functions considerably shorten analysis and evaluation time.


Fig. 3.2: Measuring functions with up to 8 markers

Despite its large 6.5 Inch VGA-TFT display, the HMS series features a very compact sized case and, due to its stackability, the HMS family easily can be combined with all other HAMEG products requiring only very little footprint. Besides using the DVI output for a projector, a standard TFT monitor can be connected - particularly useful for running manual alignment stations, built into 19" racks. In addition to its three USB ports for mass storage, printer and remote control, optional interfaces for IEEE-488 or Ethernet/USB are available.

4 Setting of parameters

Three methods of setting signal parameters are offered:

- numerical keyboard
- knob
- arrow buttons

Please use the soft menu keys for selecting the respective menu item.

4.1 Numerical keyboard

The simplest method of entering parameters quickly and exactly is the entry via the numerical keyboard. When entering parameters via the keyboard the value will be accepted upon pushing the respective unit key GHz [-dBm], MHz [dBm], kHz [dB..] or Hz [dB..]. Prior to pushing any such key an entry may be deleted by pushing the key BACK. During these operations the window will remain open. The CANCEL key will terminate the entry of parameters and close the window.


Fig. 4.1: Numerical keyboard with unit and command keys

4.2 Knob

It is possible to only use the knob for all settings. Turning the knob CW will increase the value, turning it CCW will decrease it. Pushing the knob or the ENTER key will accept also parameters without units. Such parameters can be modified only by using the knob.

4.3 Arrow buttons

The signal parameters can also be modified by using the arrow buttons: the ▲ button will increase, the ▼ button will decrease the value.

4.4 Interactive softkeys

The grey soft menu keys at the righthand side of the screen are used for the menu field displayed. Use the knob or the numerical keyboard for setting the parameter selected. If a menu field was selected via the soft menu keys, this item will be marked in blue, it is now activated for entering a parameter. If an instrument function should not be available due to a specific setting, the associated soft menu key will be deactivated, the lettering will be shown in grey.

4.5 How to enter numerical values

- Use the grey soft menu keys for the selection of a menu item.
- Enter the value of the parameter using the numerical keyboard or modify it with the knob.
- After a keyboard entry push the respective unit key.

5 Instrument functions

5.1 Setting of the frequency (FREQ)

Pushing the FREQ key will call the menu for setting the frequency. Here, the center and the start and stop frequencies can be entered. The setting is performed as described in chapter 1. The step size of the center frequency can be modified with CF-STEPSIZE. By pushing this soft menu key the settings menu will open.

- **0.1 x SPAN (Basic setting):** the frequency steps will equal 10 % of the displayed frequency range (= 1 vertical division).
- **0.5 x SPAN:** the frequency steps will equal 50 % of the frequency range displayed (= 5 vertical divisions).
- **SET TO CENTER:** The step width of the frequency is equal to the present center frequency. This mode is especially useful for the measurement of harmonics because each step will move the center frequency to the next harmonic.
- **MANUAL:** any step width is available. This allows the easy measurement of spectra with regular frequency steps.

5.2 Frequency range displayed (SPAN)

The frequency range called span is the range on both sides of the center frequency which a spectrum analyzer displays on its screen. The span to be selected depends on the signal to be analyzed, in general, it should be at least twice as wide as the bandwidth of the signal.

The series HMS offers spans from 1 KHz to 1 GHz (HMS1000/1010) resp. from 100 Hz to 3 GHz (HMS3000/3010). If the span is set to 0 Hz, the signal waveform vs. time will be shown. The

spectrum analyzer measures the signal amplitude only at the center frequency set. The setting is performed as described in chapter 4. In order to select the full frequency range of 100 KHz to 1 GHz resp. 100 Hz to 3 GHz by pushing once, the soft menu item FULL is provided. The soft menu key LAST will restore the former setting.

5.3 Setting of the amplitude parameters (AMPL)

The key AMPL is used for all settings of the amplitude displayed. The reference level (soft menu item REF.LEVEL) is identical to the top graticule line of the display. The setting is performed as described in chapter 4. The reference level setting programs the amplifier gain from the input to the display. For low reference level settings the gain will be high in order to make also weak signals well visible. For strong input signals the reference level must be set high in order to prevent overdriving of the signal amplifier chain and in order to keep the signal within the visible display window. For spectra with many signals, the reference level should be so high that all signals remain within the display area.

The receiver input will be overloaded by a disadjusted reference level.

The basic unit (UNIT) of the reference level is the unit dBm. Additionally, the unit dB μ V can be selected by pushing the soft menu key and using the knob. The range (RANGE) defines the resolution of the amplitude axis of the display. The basic scaling is in dB. The standard scaling is 10 dB/DIV. In order to obtain a higher visual resolution, the spectrum analyzer also offers the scalings 5 dB/DIV, 2 dB/DIV, and 1 dB/DIV. A higher resolution does not increase the accuracy, it only improves the readability.

The setting of the reference level will also directly affect the amount of hf attenuation at the input of the spectrum analyzer. The instrument offers 2 different modes of coupling:

- LOW NOISE (highest possible sensitivity)
- LOW DISTORTION (lowest intermodulation products)

These modes are selected via the soft menu key ATT.-SETTING. In the mode LOW DISTORTION 10 dB of additional attenuation is switched in. The preamplifier may be switched in or out with the respective soft menu key (see chapter 9.1 for activate the optional preamplifier).

5.4 Setting of the bandwidth (BANDW)

Spectrum analyzers resolve the spectral content of a signal and display a frequency spectrum. The quality of the resolution is determined by the resolution bandwidth. Additionally, the spectrum analyzers offer a selectable video bandwidth. The instrument will automatically (or, if desired, manually) choose a slower sweep time if the span was set too wide for the RBW (resolution bandwidth) selected (provided the user did not set the span to manual operation).

The video bandwidth affects the smoothing (reduction of noise) of the displayed curve. It is determined by the bandwidth of the low pass filter inserted between the video signal and the display. In contrast to the resolution bandwidth the video bandwidth has no influence on the resolution properties of the spectrum analyzer.

 If the span was set manually too wide or the sweep time to too high, the amplitudes will be displayed with incorrect level; in such cases a red „UNCAL“ message will warn. The span must then be reduced until the „UNCAL“ message disappears.


Fig. 5.1: Hf signal modulated by a sine wave signal and the resultant video signal vs. time.

By pushing the key BANDW you will enter the menu for setting the bandwidths. Both the resolution bandwidth (RBW) and the video bandwidth (VBW) may be set within the specified limits. The following step sizes are available:

RBW	VBW
100 Hz	10 Hz
200 Hz	30 Hz
1 kHz	100 Hz
3 kHz	300 Hz
10 kHz	1 kHz
30 kHz	3 kHz
100 kHz	10 kHz
200 kHz	30 kHz
300 kHz	100 kHz
1MHz	300 kHz
	1 MHz
	3 MHz

Table 5.1: Available RBW and VBW settings

Additionally automatic selection for both bandwidths (AUTO RBW/AUTO VBW) may be chosen with the respective soft menu key. The knob is used for the setting of the parameters.


Fig. 5.2: The selections offered in the RBW menu.

5.5 Setting of the SWEEP

At a frequency spectrum of $f > 0$ the sweep time is the time the spectrum analyzer requires for sweeping the selected frequency range to measure the spectrum. Certain limits have

to be observed (e.g. the resolution bandwidth set) in order to obtain a correct display.

Pushing the key SWEEP will call the selection menu. The SWEEP TIME can be varied within the specified limits. The setting of the parameters is performed as described in chapter 4. In order to assist the user when setting the sweep time, an automatic selection of the sweep time with respect to the RBW and the span settings may be chosen with the soft menu key AUTO. The automatic mode will always set the sweep time to the shortest possible value consistent with the correct display of the spectrum content.

The series HMS will sweep the selected frequency range continuously, i.e., after a sweep was completed, a new one will be started and the display refreshed. If continuous sweeping is not desired (e.g., if a single event shall be recorded upon a trigger), there is also the possibility of selecting SINGLE sweep. If single sweep is selected, the spectrum analyzer will sweep the frequency range once or it displays the video signal vs. time if the span is set to zero. The instrument will only repeat the measurement after the soft key SINGLE was pushed again.

Additionally the soft menu TRIGGER offers diverse trigger functions in order to react to events.

- **FREE RUN:** a new sweep will start after the preceding one was completed (standard setting).
- **EXTERN / and EXTERN \:** the sweep will be started upon the positive or negative slope of an external trigger signal; the external trigger signal is applied via the BNC connector EXTERNAL TRIGGER (TTL logic levels).

Use the respective soft key for the selection of the desired trigger mode.

5.6 Curve display settings (TRACE)

Pushing the key TRACE will call the settings menu.. There are several modes of curve display (TRACE MODE):

- **CLEAR / WRITE (Basic setting):** the previous curve will be erased during a new sweep.
- **MAX HOLD:** the maxima of the curve being measured and all previous ones will be displayed. MAX HOLD allows to easily find intermittent signals in the spectrum or the maximum values of varying signals.
- **MIN HOLD:** the minima of the curve being measured and all previous ones will be displayed. MIN HOLD allows to recover sine wave signals out of the noise floor or to suppress intermittent signals.
- **AVERAGE:** In this mode, the average level of consecutive measurements will be displayed. In the standard setting, averaging will be performed pixel by pixel and over the last measured curves. The average mode is hence suitable for an improved display of sine wave signals close to the noise level (available with a software update V1.2 from www.hameg.com/downloads).
- **VIEW:** freezes the curve being displayed, the measurement will be terminated, this allows to subsequently use the markers for the measurement of spectra.

The sub function TRACE ⇌ MEMORY allows to transfer a curve to the background curve memory; by pushing the soft menu

key SHOW MEMORY it will be displayed and can be compared to the presently displayed curve. The stored curve will always be shown in white and thus can be easily differentiated from the presently displayed curve. In order to let the stored curve disappear, push the SHOW MEMORY key again.

The spectrum analyzer can subtract a stored curve from an active curve and display the difference. If there is a curve stored under TRACE → MEMORY the difference between the stored and the active curves will be displayed by pushing the soft menu key TRACE MATH. In order to let the stored curve disappear push the key TRACE MATH and select OFF. Extended functions (TRACE SETUP) and memory mathematics will be available with a software update V1.2 from www.hameg.com/downloads.


Fig. 5.3: Display of a measured and a stored reference curve.

The curve in the memory (Memory Trace) will be stored in the video memory as a bitmap. The spectrum analyzer hence will not adapt the stored curve if the reference level or the displayed frequency range are changed.

5.6.1 Detector

A detector converts the video signal of a spectrum analyzer before it will be displayed. It functions pixel by pixel, determining how the value of a pixel will be measured. Pushing the soft menu key DETECTOR will call the settings menu for the selection of various types of detectors.

- **AUTO PEAK:** the spectrum analyzer will display the maximum and minimum value of each pixel from the frequency range represented by that pixel, no signal will be lost; if the signal level fluctuates (noise), the width of the curve will indicate the width of the signal fluctuations (Basic setting).
- **SAMPLE:** only displays an arbitrary point within a display pixel. The sample detector should be always used at span = 0 Hz, because this is the only method for a correct display of the video signal vs. time. Can be used for the measurement of noise power. For spans wider than than the resolution bandwidth x 501, signals may be lost.
- **MAX PEAK:** in contrast to the auto peak detector this detector will deliver only the maximum value of the spectrum within a pixel of the curve (e.g. the measurement of pulsed signals or frequency modulated signals).

- **MIN PEAK:** delivers the minimum of a spectrum within a pixel of the curve. Sine wave signals will be displayed with their correct levels while noise-like signals will be suppressed (e.g. for filtering sine wave signals from noise).

5.7 The use of markers

The series HMS offers several markers and delta markers for the evaluation of curves. The markers are always tied to the curve and indicate the frequency and the level at that. The frequency position of the marker is marked by a arrow icon. The numerical values of the frequency and the level are shown at the top of the screen. The unit of the level is the same as the unit selected for the reference level.

The knob allows to choose from 8 different markers. The individual markers can be switched on and off with the respective soft menu key. The soft menu key POSITION is used to set the frequency position of the marker along the curve.

The delta marker level is always relative to the level of the main marker, the unit of level is always dB. The value of frequency of the delta marker is always relative to the reference marker (marker 1) and thus indicates the difference in frequency to that marker. The use of the delta marker will be available with a software update V1.2 from www.hameg.com/downloads.

5.8 Peak-Search

The so-called Peak-Search key will show the user the display of the next maximum value. By pushing the softkey PEAK the next maximum value can be shown. Extended functions will be available with a software update V 1.2 from www.hameg.com/downloads.

5.9 Operation in the Receiver-Mode

By pushing the MODE key the selection menu will be called which allows to switch to the receiver mode. The spectrum analyzer acts as a receiver which measures the level at a preselected frequency. The most important parameters such as e.g. frequency, amplitude, resolution bandwidth may be set using the appropriate keys. In the receiver mode the same bandwidths are available as in the spectrum analyzer mode. Additionally the bandwidths: 200 Hz, 9 KHz, 120 KHz and 1 MHz are available for emi emission measurements according to CISPR. The key FREQ and the softkey DETECTOR selects the detector (Peak, rms, quasi-peak). The measurement time is the time during which the spectrum analyzer collects measurements and combines them according to the detector selected for a display.

With the soft menu item AUDIO the series HMS offers an AM and a FM demodulator allowing listening to modulated signals. The demodulated signal may be listened to with a headphone and an intern speaker. The headphone is connected to the headphone connector (3.5 mm female connector). If the headphon is activated, the intern speaker will be deactivated. The respective soft menu keys allow to switch the demodulator on or off and to set the volume.

6 Store and recall instrument settings

Your spectrum analyzer can store 5 different kinds of data:

- Instrument settings
- Reference signals
- Waveforms
- Screen displays
- Sets of formulae.

Signals and screen displays can only be stored on USB sticks. All other data can be stored either on a USB stick or in the instrument's non-volatile memories.

6.1 Instrument settings

Push the key SAVE/RECALL for calling the main menu for storage and load functions. First a listing is shown of the kinds of data which can be stored and loaded. By pushing the key next to the top menu item DEVICE SETTINGS this menu will open.


Fig. 6.1: Basic menu for instrument settings

In this menu, by pushing the respective key, it is possible to call the menu for storing, the data manager for loading, and the menu for exporting and importing instrument settings. Additionally, the menu item DEFAULT SETT. will reset the instrument to the factory settings. The storing menu is opened by pushing the SAVE key.


Fig. 6.2: Saving instrument settings

Here the storage location (internal memory, front panel USB, rear panel USB) is selected, also a name and a commentary can be added; these will be stored by pushing the soft menu key

next to SAVE. In order to recall stored instrument settings, call the main instrument settings menu and select LOAD by pushing the respective soft menu key. The data manager will open, use the menu keys and the universal knob for navigating.


Fig. 6.3: Loading instrument settings

Here the location is selected from which the settings data are to be loaded. After the selection in the data manager, load the settings by pushing the soft menu key LOAD. The data manager also allows you to erase individual settings in the internal memory. If a USB stick is plugged in and has been selected as the location, it is also possible to change or erase directories. In order to export or import instrument settings, a USB stick must be plugged in, otherwise this menu can not be accessed. Provided this is fulfilled, pushing the key next to IMPORT/EXPORT will open a menu allowing to copy instrument settings between the internal memory and a USB stick.


Fig. 6.4: IMPORT / EXPORT menu for instrument settings

The source is selected by pushing the respective key (e.g. INTERNAL), the selection will be indicated by its blue background. Then the destination is selected (e.g. FRONT). By pushing the key next to IMPORT/EXPORT, the selected settings data will be copied as previously chosen (in this example from the internal memory to a USB stick). It is possible to copy from the internal memory to the external memory and also between two USB sticks.

6.2 Reference signals

References are sets of data which consist of settings information. These may be stored and recalled internally or externally. The main feature of references is the fact that all information is always stored along with the data proper; this allows to compare the reference with live signals. If the key SAVE/RECALL is

pushed and the menu item REFERENCES selected, a changeover into the menu IMPORT/EXPORT is possible; here the standard menu of the data manager will appear which allows you to copy references between the internal memory and an external USB stick (refer to chapter 3.1 for a detailed description).

6.3 Waveforms

In addition to references, the waveform data can be stored only on external USB sticks, not internally.

HAMEG Binary format: A binary data set may contain bytes of any length. The curves will be stored without any time information.

CSV (Comma Separated Values): CSV data sets store the curves in tables, the lines are separated by commas.

HRT (HAMEG Reference Time): Data sets with this code contain data of curves vs. time. If a waveform was stored in this format, it can be used in the reference menu. With the HRT format it is also possible to generate data sets which may be reloaded into the spectrum analyzer via the reference menu.

In order to store waveforms, push the key SAVE/RECALL and select in the main menu the item CURVES by pushing the respective softmenu key.


Fig. 6.5: Menu to save a waveform

In this menu which will open, the top item allows the selection of the front or rear panel USB port. This choice is only possible if the instrument recognized a USB stick at the designated port. If a stick is present and the port selection done by pushing the respective soft key, the first time this happens, the data manager will appear with the associated menu. Here, a listing of destinations can be selected or generated. Confirm the selection of the destination listing by pushing OK, this will recall the menu for storing curves. Pushing the soft key next to the second menu item (CURVE) will activate this function as indicated by the blue background: now the channel can be selected from which the curve shall be taken by turning the universal knob. Only channels which have been activated are eligible.

Pushing the menu key next to DATA NAME will open the menu for entering names: in order to do this first push the CURSOR/SELECT key, then use this menu and the universal knob to enter the desired name which will be stored by pushing ACCEPT. This will recall again the menu for storing waveforms. Now push the soft key FORMAT, this will open a window for selecting the format. The selection is performed again with the universal knob. Additionally, a commentary can be stored along with a curve.

This is done by pushing the menu key next to COMMENTARY, this will open a window for the entry. After entering the commentary and storing it by pushing ACCEPT, again the menu for storing curves will appear. After completion of all these entries, pushing the menu key next to STORE will store the curve according to the selected settings.

6.4 Screenshots

The most important method of storing for documentation purposes is the screen photo. At least one USB stick must be connected, only then will any settings regarding the destination and the format be possible. Push the keys SAVE/RECALL and SCREENSHOTS for opening the appropriate menu.


Fig. 6.6: Menu for screenshots

Also in this menu the destination (according to the USB sticks connected) can be selected with the top menu key. When this is done the first time, the data manager will appear in order to either select or generate a destination listing. After the entry of this information, the SCREENSHOTS storing menu will reappear. The second menu item FILE NAME allows you to enter a name with the respective name entry menu which will open automatically upon selecting this menu item. If FORMAT is selected with the respective menu key, these formats will be offered and can be selected with the universal knob: BMP = Windows Bitmap (uncompressed format) and GIF. Pushing the key next to SAVE will store the actual screen display along with the name and format at the destination selected.

6.4.1 Screenshot example

In order to store data you have to define the kind of data and the destination. First attach a USB stick (refer to 8.1 USB connector) to the front panel connector. Press SAVE/RECALL in order to call the respective menu.


Fig. 6.7: Save/Load menu

Select the kind of data by pressing the respective soft key (in this example SCREENSHOTS) in order to access the settings menu.


Fig. 6.8: Menu with the settings for screenshots

Please verify that the USB connector into which you plugged the USB stick (front or rear) is written in the top softmenu (you can change the destination by opening the respective menu if you press the softkey next to STORAGE). You can now save a Screenshot if you press the softkey next to SAVE using the predefined name written in the menu below FILE NAME. You may name the destination memory with up to 7 characters; in order to do this select the menu item FILE NAME and define the name by using the knob and the CURSOR SELECT key (in this example PRINT).


Fig. 6.9: Defining file names

After the softkey next to ACCEPT was pressed the oscilloscope will have stored the name and return to the settings menu. Here you can now store the actual screen display by pressing the SAVE softkey. Alternatively, you can return to a lower menu level (by pressing the lowest Menu OFF key) and select the menu item key FILE/PRINT. In the following menu press the soft menu key next to SCREENSHOTS: this will assign the function screenshot to the key FILE/PRINT with the settings chosen. This enables you to store a bitmap file on your USB stick by just pressing FILE/PRINT at any time and in any menu.


Fig. 6.10: Settings of the button FILE/PRINT

7 Extended operating modes

7.1 Using the help function

The integrated help function can be activated by pressing the key HELP in the GENERAL area of the control panel. A window will open and the text inside is dynamically updated depending on the key (including softmenu key's) you are pushing or the knob you are turning. If you do not need the help anymore, you can switch off the help window by pushing the "HELP" key.

7.2 Display settings

By pushing the key DISPLAY the display settings menu will be called; here several choices are offered:

- **TRACE:** Adjustment of the trace intensity (0 ... 100 %) of the displayed spectrum.
- **BACKLIGHT:** Adjustment of the backlight intensity (0...100 %).
- **GRID:** Adjustment of the raster intensity (0 ... 100 %). The soft menu item GRID SETUP allows to select a cross, raster lines or no raster with the respective soft menu keys. Also the raster designations (SCALE) can be switched on or off. The soft menu key LED INTENS changes the LED intensity from dark to light, this is effective for all backlighted keys and all other display LED's on the front panel.
- **TRANSPARENCY:** Adjustment of the transparency (0 ... 100 %) of the raster inscriptions.


Fig. 7.1: Display settings menu (DISPLAY)

If a soft menu item is activated, its background will be blue. The setting of the parameters is performed according to chapter 1.

7.3 Selection of the standard (preset) instrument settings (PRESET)

By pushing the key PRESET the spectrum analyzer will resume its preset standard settings. This allows to generate a new configuration, starting out from defined parameters, no parameter from a former setting will be active any more.

- Center frequency: 1.5GHz (HMS3000/3010)
 500MHz (HMS1000/1010)
 Span: 3GHz (HMS3000/3010)
 1GHz (HMS1000/1010)

8 General instrument settings

Basic settings like language for user interface and help, miscellaneous settings and interface settings can be set using the menu which opens after pressing the SETUP key in the GENERAL area of the control panel. Pushing the soft menu key MENU OFF will call the next lower level.

8.1 Language settings

The series HMS provides 4 different languages for the menu and help text:

German, English, French and Spanish

By pushing the soft menu key LANGUAGE the language selection is called, the language selected is active if the menu item's background is blue.

8.2 Basic settings

8.2.1 Clock setting

Pushing the soft menu key SET CLOCK will call the clock and date settings menu. These settings will be used for adding a time and date stamp on print-outs and stored files. The user can modify the time and date with the knob. The respective soft menu item is active if its background is blue. The time and date settings will be accepted by pushing ENTER.

8.2.2 SOUND

The series HMS offers the possibility to sound a warning which can be switched on or off using SOUND. The control resp. warning tone will be active if the respective menu item's background is blue.

8.2.3 Name

In this menu item you can set a name for the series HMS. By pressing the softkey a key panel will show. You can choose the character via the knob. The character will confirm with the enter button (refer to chapter 6.4.1).

8.2.4 About

Choosing this soft menu item will call instrument information such as serial number, software version etc.


Fig. 8.1: Instrument Informations

8.3 Interface settings

Selecting this soft menu item will allow to modify the settings for:

- the Dual Interface H0720 USB/RS-232 (Baud rate, number of stop bits, parity, handshake on/off)
- LAN Interface H0730 (IP address, sub net mask etc., see the manual of the H0730) and
- the IEEE-488 GPIB interface H0740 (GPIB-address)

The interface desired for the communication can be selected with the respective soft menu key. Use the soft menu item PARAMETER to set the necessary interface parameters

8.4 Printer settings

The series HMS supports printing of the screen contents on a connected printer (USB printers with postscript).

8.5 Update (Firmware / Hilfe)

You are invited to download the most recent firmware under www.hameg.com. Firmware and help are packed into one ZIP data packet. After downloading the ZIP data unpack it into an USB stick's basic directory (refer to 6.1 USB connector). Thereupon insert the stick into the USB port of the spectrum analyzer and push the key SETUP in the GENERAL area of the front panel. Here you shall find the menu item UPDATE. After selecting this menu item a window will open which displays the actual firmware version indicating the version number, the date and build information.

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Now choose which to update: the firmware or the help function. If both are to be updated it is recommended to first update the firmware. After you selected firmware updating by pushing the appropriate key the respective date will be searched on the stick, the information of the firmware to be updated from the stick will be displayed below the line NEW. In case the new firmware should be identical to the existing one, the number of the version will be shown in red, otherwise it will be shown in green; only then should you activate the updating by pushing the soft key EXECUTE. If you intend to update the help function or add a help language choose HELP in the updating menu. The information window will now display the languages installed, the date, and the information about the languages available on the stick. With the soft menu, languages may be added, removed or updated.

Please note the format of the date:YYYY-MM-DD according to the multi language norm of ISO 8601.


Fig. 8.3: Info display of help update

9 Optional Accessories

9.1 Activation of the Preamplifier H03011

The Preamplifier for the series HMS (DANL -135dBm typ. / 100 RBW) will activate in the Setup menu with the softkey UPDATE. Please look at chapter 8.5 for the procedure of the update (same procedure like the Firmware-Update). This preamplifier isn't included in the supplied accessories and can be acquired by purchase.

9.2 19" Rack mount kit 4HE HZ46

For the application in rack systems HAMEG provides a kit for the HMS series. Technical details and a description about the mounting you can find in the manual HZ46 on our homepage <http://www.hameg.com/downloads>.

9.3 Carrying case HZ99

The Carrying Case HZ99 is used to transport your spectrum analyzer and is available 'on stock'.

Fig. 9.1: Carrying case HZ99

9.4 Near field probe HZ530/HZ540

The set includes 3 hand-held probes with a built-in preamplifier covering the frequency range from 100kHz to 1GHz resp. < 1MHz to 3GHz. When used in conjunction with a spectrum analyzer or a measuring receiver, the probes can be used to locate and qualify EMI sources, as well as evaluate EMC problems at the breadboard and prototype level. The power can be supplied either from batteries (HZ530) or through a power cord directly connected to a spectrum analyzer (HZ540). Signal feed is via a BNC-cable or SMA/N-cable. They enable the user to evaluate radiated fields and perform shield effectiveness comparisons.

The probes - one magnetic field probe, one electric field probe and one high impedance probe - are all matched to the 50Ω inputs of spectrum analyzers. The technical specifications are described in the HZ530/HZ540 manual on our homepage <http://www.hameg.com/downloads>.

9.5 Measurements of spectra with a VSWR bridge HZ547 (HMS1010/3010)

Fig. 9.2: VSWR bridge HZ547 for HMS1010/3010

The VSWR bridge HZ57 allows the measurement of the voltage standing wave ratio (VSWR) and the reflection coefficient of 50Ω devices. Typical objects are e.g. 50Ω attenuators, load resistors, amplifiers, cables, mixers, frequency selective devices. The frequency range is 150kHz...1GHz (HMS1010) resp. 100kHz...3GHz (HMS3010). The technical specifications and the measurement set-up are described in the HZ547 manual on our homepage <http://www.hameg.com/downloads>.

9.6 Transient Limiter HZ560

The Transient Limiter HZ560 protects the input circuit of spectrum analyzers and measurement receivers, in particular in combination with the use of a Line Impedance Stabilization Network (i.e. LISN HM6050). The technical specifications are described in the HZ560 manual on our homepage <http://www.hameg.com/downloads>.

9.7 75/50-Ω-Converter HZ575

The converter HZ575 has a 75Ω AC coupled input and a 50Ω DC coupled output. Using HZ575, spectrum analyzers with 50Ω input can be used for measurement in 75Ω environment. HZ575 can be used reversed too. A 50Ω signal applied at the 50Ω "output" is present at the 75Ω BNC "input" socket with an impedance of 75Ω. The technical specifications are described in the HZ575 manual on our homepage <http://www.hameg.com/downloads>.

10 Front panel connections

10.1 USB connector

Using the front panel USB connector a software update of the HMS firmware can be performed with a FAT or FAT32 formatted mass memory (chapter 8.5).

10.2 PHONE

The signal available at this connector comes from an AM detector and it helps to identify the sources of interference e.g. when making precompliance measurements. If an antenna is connected to the analyzer input, selecting CENTER and using the knob the analyzer can be tuned to a transmitter (Receiver mode chapter 5.9). The demodulation has to be activated. Please note that this operational mode may be subject to national restrictions!

10.3 PROBE POWER

This connector can be used as a supply (6 V DC) e.g. for HAMEG probes. The inner contact is +6 V, the outer contact is connected to the instrument housing and thus with the measurement inputs' ground potential and also protective earth (PE).

10.4 EXTERNAL TRIGGER

The external trigger input connector is used for the control of measurements by an external signal. (TTL levels.)

10.5 OUTPUT 50Ω (Tracking Generator)

The tracking generator output (only HMS1010 and 3010) has to be connected to the measuring object with a N connector cable. A test signal with a spectrum from 5 MHz to 1 GHz resp. 3 GHz is available.

10.6 INPUT 50Ω

Without attenuation (ATT 0dB) ±25VDC must not be exceeded. With an attenuation of 10 to 50 dB, the maximum level is +20 dBm. Levels or DC voltages above the values mentioned may destruct the input stage. The outer contact is connected to the instrument chassis and thus to safety ground (PE). The maximum input levels resp. voltages must not be exceeded. Danger of destruction!

11 Rear panel connections

11.1 USB connector

The USB interface on the rear panel can be used for a printer (see chapter 8.4).

11.2 DVI connector

The rear panel of the spectrum analyzer also has the standard DVI connector for the connection of external monitors or projectors.

11.3 REF IN / REF OUT

In order to further increase the frequency stability, the internal oscillator may be replaced by an external one which can be connected to the „10 MHz REF IN/REF OUT“ connectors on the rear panel. The external reference frequency signal must comply with the specifications given with respect to frequency accuracy and amplitude.

The switching between internal and external reference frequency can be effected via the button SETUP and the softkey REF. FREQUENCY.

12 Remote Control

The HMS series is basically supplied with an USB/RS-232 interface. The respective drivers are available on the enclosed Product CD or can be downloaded at <http://www.hameg.com>.

To establish a basic communication a serial cable (1:1) as well as a terminal program like Windows HyperTerminal is required. The Windows HyperTerminal program is part of any Windows operating systems. A detailed instruction how to setup a basic communication using HyperTerminal is available at the HAMEG Knowledge Base at <http://www.hameg.com/hyperterminal>.


The HMS series uses SCPI (= Standard Commands for Programmable Instruments) for remote control. Remote control is possible via the built-in dual interface USB/RS-232 (options: Ethernet/USB, IEEE-488). This allow access to nearly all functions which are available on the front panel. A detailed dokument about the provided SCPI commands is available at <http://www.hameg.com>.

13 Appendix


List of figures

Fig. 3.1:	Menu of the frequency	32
Fig. 3.2:	Measuring functions with up to 8 markers	32
Fig. 4.1:	Numerical keyboard with unit and command keys	32
Fig. 5.1:	Hf signal modulated by a sine wave signal and the resultant video signal vs. time.	34
Fig. 5.2:	The selections offered in the RBW menu.	34
Fig. 5.3:	Display of a measured and a stored reference curve.	35
Fig. 6.1:	Basic menu for instrument settings	36
Fig. 6.2:	Saving instrument settings	36
Fig. 6.3:	Loading instrument settings	36
Fig. 6.4:	IMPORT / EXPORT menu for instrument settings	36
Fig. 6.5:	Menu to save a waveform	37
Fig. 6.6:	Menu for screenshots	37
Fig. 6.7:	Save/Load menu	38
Fig. 6.8:	Menu with the settings for screenshots	38
Fig. 6.9:	Defining file names	38
Fig. 6.10:	Settings of the button FILE/PRINT	38
Fig. 7.1:	Display settings menu (DISPLAY)	39
Fig. 8.1:	Instrument Informations	39
Fig. 8.2:	Aktualisierungsmenü Firmware	40
Fig. 8.3:	Info display of help update	40
Fig. 9.1:	Carrying case HZ99	40
Fig. 9.2:	VSWR bridge HZ547 for HMS1010/3010	41

Oscilloscopes

Spectrum Analyzer

Power Supplies


Modular System
Series 8000

Programmable Instruments
Series 8100

authorized dealer

42-3000-0020

www.hameg.com

Subject to change without notice
42-3000-0020 (1) 23092009
© HAMEG Instruments GmbH
A Rohde & Schwarz Company

DQS-Certification: DIN EN ISO 9001:2000
Reg.-Nr.: 071040 QM

HAMEG Instruments GmbH
Industriestraße 6
D-63533 Mainhausen
Tel +49 (0) 61 82 800-0
Fax +49 (0) 61 82 800-100
sales@hameg.de

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© Copyright 2006 by Conrad Electronic Benelux B.V. 2006/06/RC.