Axel Semrau®

MOSH/MOAH Recap and new developments

Starting point

October 2009

- Cantonal Laboratory of Zurich informed that they found mineral oil hydrocarbons (MOH) in food:
- Results up to 100 mg/kg cacao-powder (Σ MOSH +MOAH)
- Rice: 14 mg/kg MOSH (<C25)
 4,0 mg/kg MOAH (<C25)
- High results especially in food with large specific surfaces:
cacao-powder, rice, semolina, noodles, pasta, sauce thickener
- These products are often packaged in cardboard

Eur Food Res Technol (2010) 230:785–796
DOI 10.1007/s00217-010-1223-9

ORIGINAL PAPER

**Is recycled newspaper suitable for food contact materials?
Technical grade mineral oils from printing inks**

Maurus Biedermann · Koni Grob

Migration of mineral oil from packaging materials to foodstuffs

- Due to the large portion of **mineral oil fractions with short chain and aromatic hydrocarbons**, such **contaminations of foodstuffs are adverse**
- **Short chain hydrocarbons** are more **easily taken** up by the body
- **Frequent intake** of such contaminated foodstuffs can thus **lead to exceedances of the toxicological limit values**
- Animal studies have shown that **mineral oil mixtures with low viscosity are stored in the body** and can lead to **accumulations and damage** in the liver, heart valves and lymph nodes

Conclusion: urgent need to minimize the migration of mineral oils to foodstuffs

Possible courses of action

- use of inner bags in order to prevent substance migration from the cardboard box
- use of cardboard made from virgin fibres
- changing the composition of printing inks, mineral oil free systems for newspaper printing


This suggestion is further supported by the fact that the uptake of mineral oils from printing inks can also occur through skin contact

Starting point

The reasons

- The load of recycled cardboard with mineral oil components
- The printing of food-packaging with unsuitable systems

- ❖ Transfer via gas phase into the food

a. Depending on vapor pressure

b. Boiling range up to <C25

Source for the introduction of mineral oil into food

- Contamination of recycled fibres by mineral oil components
- Inks used for the offset printing which contain 20-30% mineral oil as solvent
e.g. newspapers, leaflets, packaging material
- The recycling process cannot get all of the mineral oil out, so some of it ends up in recycled cardboard and paper
- Contamination of virgin fibres via gas phase from outer packaging
- Of special concern are dry foods having a large specific surface, containing fat and with long shelf life

If recycled fibres are used as raw material for food packagings, mineral oil can contaminate the food by transfer via gas phase

Most known Mineral Oil Sources

Environment

Traffic & Industrial Pollution

Lubricating oils
(belt conveyor, agitators)

Pesticide Formulations

Detergent

Cosmetics

Food

Batching Oil for Jute- and
Sisal- Bags

Bakery & Sweets release oils

Dust binding (rice, cereal,
animal feed)

Fruit & Dried Fruit Polishing

Packing materials
Waxes (water repellency), Inks,
Polyolefins, Adhesives

What is mineral oil ?

Mineral oil

Aliphatic fraction

MOSH: Mineral Oil Saturated Hydrocarbons
(ca. 65-90 %)

Paraffinic (open chain, mostly branched)
and naphthenic (cyclic) hydrocarbons

Distribution of chain length in printing inks is
centered at C17 – C19

Aromatic fraction

MOAH: Mineral Oil Aromatic Hydrocarbons
(ca. 10-35 %)

1 – 4 aromatic rings, partly hydrogenated,
mostly highly alkylated

Structures have different degree of alkylation,
alkyl side chains differ in length and branching

Additional approach

➤ June 2010

2 days Workshop of the Federal institute for Risk assessment (BfR) and the Cantonal Laboratory of Zurich

Mineral oil analytic

- ❖ Themes : Methods for the MOSH und MOAH analytic (*GC-FID, on-line HPLC-GC-FID, GCxGC*)
Occurrence and risk assessment

➤ September 2011

Conference “Mineral oils in food packaging - developments and solutions”, Berlin, BfR
International discussion round with 400 participants
Toxicology, Analytic

➤ November 2011

Method development kit for mineral oil, BfR
Different samples: Rapeseed oil; Mixture of Mineral oil; Mixture of n-alkanes;
Extract of recycled cardboard; Mineral oil contaminated rice

Full description of the mineral oil kit

Standard operating procedure (SOP) of the „Manual“ method Mineral oil contaminated

Decision support project (BMELV)

Recycled Paper – Appropriate for Food Packaging?

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BMELV (Federal Ministry of Food, Agriculture and Consumer Protection) Project 09HS012 (2010-2012)

see: <http://download.ble.de/09HS012.pdf>

Confirmation of the results

- **Food from the retail 119 samples study:**

119 dry food samples packed in recycled board were stored at room temperature and analyzed at the best before date or after 1.3 years storage at the latest

Max. amounts in food:	MOSH : 100 mg/kg	MOAH: 16 mg/kg
Main average results :	MOSH: 17 mg/kg	MOAH: 3 mg/kg

- **Migration study – storage experiments**

After 2, 4 and 9 months storage at room temperature the food was analysed for migrated contaminants

Migration in food directly packed in cardboard (% of potential):


MOSH: 35 mg/kg	MOAH: 7 mg/kg (< C24: 65 to 80 %)
----------------	-----------------------------------

Functional barriers foils in this experimental setup :


- PE no barrier for long term storage
- PP no barrier for long term storage
- PP/acrylate no absolute barrier
- PET tight barrier, with and without Al metallization

„Packaging material“

„Dry Food“

Chocolate from advent calendars

November 2012

- The Goods Testing Foundation (Stiftung Warentest) has tested chocolate from advent calendars and found that they contained mineral oils




24 samples	MOAH (mg/kg)	24 samples	MOSH,POSH;PAO (mg/kg)
15 samples	< 0,5	21 samples	> 0,5 - 10
6 samples	> 0,5 - 3	3 samples	> 10
3 samples	> 3 (max. 7)		

8 samples with mineral oil from other sources → lubricating oils

- Federal Institute for Risk Assessment (BfR) has conducted an initial evaluation to establish whether mineral oils in chocolate can pose a health risk
- It is possible to avoid such contamination, because there are calendars where none of these substances was found in the chocolates

„Advent calendars“

- no printing inks with MOH
- no recycled paper
- partially with barriers
- low basic load with MOH

Additional approach (QM based)

Quality management

- Proficiency Test March 2013 organized from Institut Kirchhoff Berlin
 - 17 laboratories
 - 4 different samples (chocolate, recycled cardboard, rice, hazelnut oil)

- German reference office for food proficiency testing and reference materials (DRRR) (commercial provider)
 - a. Mineral oil in cardboard
 - b. Mineral oil in food
 - c. Mineral oil in fatty foods
 - d. Mineral oil in Edible fats & oil
 - e. Migration of mineral oil from cardboard - migration of mineral oil in the simulant: Tenax

More info: www.drrr.de

Additonal approach (QM based)

➤ Collaborative trial, CEN, April 2015 organized from ITERG, France

- Mandat M/243by European Commission to CEN (CEN/TC275/WG13)
- „*Vegetable oils and foodstuff on basis of vegetable oils- Determination of mineral oil saturated hydrocarbons (MOSH) and mineral oil aromatic hydrocarbons (MOAH) with on-line-HPLC-GC-FID analysis*
- 12 laboratories (from 6 European countries)
- 8 oil samples (4 sunflower-, soybean-, olive-, olive pomace-, palmoils)
- 2 food samples (mayonnaise, margarine)
- partly with epoxidation and ALOX-Purification

❖ **Conclusions:**

- palm oil or olive oil need an epoxidation step for the MOAH content determination
- Horrat values and considering as acceptable a value of 25 % for the reproducibility relative standard deviation for this specific determination

LOQ should be fixed at 10 mg/kg for MOSH and MOAH



Cosmetic products

May 2015

- Goods Testing Foundation tested **25 cosmetic products** (creams, baby care and lip care products, body oils, hair waxes and vaseline)
- All samples contained MOAH amount in range of 0.005 % to 9 %
- The Federal Institute for Risk Assessment (BfR) reported by up to 5 % MOAH content
- MOAH are especially problematic in products for lip care because these can be swallowed for the most part

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Basics of the LC-GC Technique

The method was developed from Konrad Grob and Maurus Biedermann at the Kantonales Labor Zürich.

Axel Semrau modified the method and brought it to routine labs

LC-GC coupling

Heart-cut technique

LC as sample preparation for GC

- LC: High sample capacity → Clean-up, separation of the matrix
- GC: High separation efficiency, selective detectors


Mostly normal phase LC or SEC for GC-compatible solvents

Why do you need HPLC-GC for MOSH-MOAH?

Axel Semrau®

No GC separation between MOSH / MOAH on an apolar GC columns possible

(NP)LC-Conditions

2 mm-columns

- Sufficiently high capacity for sample and matrix
- Group sizes manageable
- Flow rates: 200 - 500 $\mu\text{L}/\text{min}$ \rightarrow optimal GC-transfer rate
- Intermediary detector for control
- Time or signal-controlled fractionation

After transfer often back-flush used with polar solvent

- Example: direct injection of oil \rightarrow example: removal of triglycerides

Requirements for the LC-GC coupling used for MOSH/MOAH analysis


- **Has to work with a FID as detector**
- **Must be compatible with at least 2mm columns, smaller diameter is not possible**
- **Solvent Hexan/Dichlormethane must be possible**
- **No discrimination of C10**
- **Up to C40, better C50+**
- **Heart-cut technique**
- **Easy to handle**
- **Robust**
- **2 LC Fractions in one run should be possible**

All this requirements can only be fulfilled with a LCGC solution based on retentiongap technology, Solvent split techniques based on aPTV system are not suitable.

Understanding the principle

Extraction of food / paper with hexane and hexane / ethanol Maybe after concentration, direct injection into NPLC-GC-FID

- 2,1 mm x 250 mm Restek Allure 60 Å
- Column suitable for fatty foods
 - 20 mg Trapping capacity Triglyceride ca. 50% of the column → Trapping function

Pragmatic verification standard:

Cholestane → End of MOSH

TBB → beginning of MOAH


Perylene → end of MOAH

C11: Volatile component (MOSH)

Pentylbenzene: Volatile component (MOAH)

LC-gradient

- Separation of MOSH and MOAH with Hexane/CH₂Cl₂ gradients
- UV-detection (230 nm)
- Gradient tracking
- 2 fractions with each 450 µL


450 μ L-Fractions in Hexane (MOSH) or Hexane/ CH_2Cl_2 (MOAH)

- Different enrichment conditions necessary


C11 and Pentylbenzene are analyzed quantitatively

- 10 m x 0.53 mm retention gap
 - Blank, deactivated guard column,
 - most suitable: Restek MXT
 - Solvent Vapor Exit
 - Partial Concurrent Solvent Evaporation
- 15 m x 0.25 mm x 0.25 μ m separation column
 - No separation of all hydrocarbons isomers possible
 - Quantification on signal hill
 - Short column, high carrier gas flow (H_2)
 - Maximum S/N-ratio
 - most suitable: Restek RXI-5Sil MS


LC-GC technique: Partial Concurrent Solvent Evaporation

LC-GC technique: Partial Concurrent Solvent Evaporation


LC-GC technique: Partial Concurrent Solvent Evaporation


LC-GC technique: Partial Concurrent Solvent Evaporation

SOLVENT TRAPPING ENDS


LC-GC technique: Partial Concurrent Solvent Evaporation

SEPARATION OF THE VOLATILES


LC-GC technique: Partial Concurrent Solvent Evaporation

CONCENTRATION OF THE HIGH BOILERS

LC-GC technique: Partial Concurrent Solvent Evaporation

SEPARATION OF THE HIGH BOILERS


The interface – CHRONECT LC-GC

To be connected with any GC
1 or 2 channel operation

Complete CHRONOS Integration

General System Design

Agilent 1260 with UV-det
Agilent 7890 GC with FID
CTC PAL
Clarity
Chronos

Samplelist using Chronos

Axel Semrau®

Chronos - New Sample List*

Main menu

- Sample list
- Schedules and run control
- Method editor
- Settings

Sample

- Add...
- Remove
- Duplicate
- Move up
- Move down

List

- New
- Load...
- Save
- Save as...
- Print...
- Import CSV...
- Export...

Schedule

- Overlapped
- Priority
- Infinite loop
- Create

	Analysis Method	Comment	Source Tray	Source...	Inj. Volume [µL]	1. Fraction star...	1. Fraction l...	2. Fraction start [...]	2. Fraction length [s]	LC Runtime [s]
1	C:\Use...sh.cam	MOSH/MOAH	Tray1	1	50	90	90	4.5	90	2100
2	C:\Use...sh.cam	MOSH/MOAH	Tray1	2	50	90	90	4.5	90	2100
3	C:\Use...sh.cam	MOSH/MOAH	Tray1	3	50	90	90	4.5	90	2100
4	C:\Use...sh.cam	MOSH/MOAH	Tray1	4	50	90	90	4.5	90	2100
5	C:\Use...sh.cam	MOSH/MOAH	Tray1	5	50	90	90	4.5	90	2100
6	C:\Use...sh.cam	MOSH/MOAH	Tray1	6	50	90	90	4.5	90	2100


Main menu

Activities


CombIPAL: not connected

100 %

Schedule for samples

Method in Chronos

The screenshot displays the 'Chronos - New Sample List*' application window. The interface is divided into several sections:

- Main menu:** Includes icons for 'Sample list', 'Schedules and run control', 'Method editor', and 'Settings'.
- Method editor:** Contains buttons for 'New', 'Load...', 'Save', 'Save as...', and 'Print...'. Below these are 'CDS method' options: 'ChemStation', 'ChromPerfect', 'Clarity', 'Eksigent NanoLC', and 'Xcalibur'.
- General tab:**
 - Information and settings:** Fields for 'File name' (C:\Users\Public\Documents\Chronos\Methods\Demo\LC-G), 'Method name' (LC-GC Dual Channel Method), and 'Description'. Includes checkboxes for 'Error method', 'Overwrite device settings', and 'Show invisible task properties'.
 - Columns:** A table defining method parameters.
 - Tasks:** A list of tasks with descriptions, such as 'Acquisition', 'WaitOverlapped', and 'Injection'.
 - Properties:** A table of method properties and their values.

Name	Token	Cell type	Variable type	Visible	Editable	Sortable	Standard value
Comment	%COMMENT%		Text	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Source Tray	%TRAY%	Tray with vial		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Tray1:1
Inj. Volume [µL]	%INJ_VOL%		Decimal number	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	10
1. Fraction start [min]	%1FR_START%		Decimal number	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2
1. Fraction length [s]	%1FR_1FN%		Positive number or zero	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0n

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
- **MOSH + MOAH in the same GC-run**

Status instrumentation 2017

- **50 systems installed in routine labs in Germany, Austria, Benelux, France**
- **Service and support by Axel Semrau; SIM GmbH; Instrument Solution; Sugelabor**
- **Development of FAT and SAT procedure by Axel Semrau**
- **Further automation options and improvement of the methodology**
- **First european regulation**

DIN EN 16995:2016-05 Draft :

Foodstuffs - Vegetable oils and foodstuff on basis of vegetable oils -
Determination of mineral oil saturated hydrocarbons (MOSH) and mineral oil
aromatic hydrocarbons (MOAH) with on-line HPLC-GC-FID analysis; German
and English version prEN 16995:2016

will be finished in 2017

New Developments for MOSH/MOAH

Axel Semrau®

- **Fract and Collect option for further identification of MOAH fraction to example by GCxGCMS**
- **Automated expoxidation**

Axel Semrau®

Fract and Collect Option

- **Requires additional valve on the system**
- **Requires PAL RTC**
- **Complete integration in samplelist**
- **Collects a complete Fraction (MOSH or MOAH; selectable by user)**
- **Fraction can be used for further identification of compounds using other techniques like GCxGC TOFMS or NMR**



GCxGC/TOFMS



jutebag

coffee

m/z 184 198 212 – Dibenzothiophene, 1-2C alkylated Dibenzothiophenes, 4C alkylated Naphthalenes

m/z 152 154 168 182 196 Biphenyl, 1-2C Biphenyl, Acenaphthylene, Acenaphthene, Dibenzofuran, Dibenzothiophene, Phenanthrene

Automated Epoxidation


Axel Semrau®

- **Requires additional modules like centrifuge; solvent reservoir; tools**
- **Requires PAL RTC**
- **Modification of manual method for better automation**

Why Epoxidation ?



Axel Semrau®


Selective derivatisation of matrix compounds

- Olefins are derivatised to be removed during LC
- Epoxidation using mCPBA

Automated Epoxidation



- No solvent exchange necessary
 - Epoxidation direct in Hexane instead of DCM
 - Ethanol as solvent for mCPBA
- Complete automation is possible!

→ Quenching using Na₂S₂O₃

Example Palmoil

Axel Semrau®

Summary

- **LCGC-FID is a routine method for the analysis of MOSH/MOAH**
- **Runs in a lot of routine labs 24/7**
- **Further automation steps were developed to reduce manual interaction**
- **Service and support infrastructure is available by a strong network of different partners**

Axel Semrau®

Thank you for your attention!

 Instrument
Solutions
Technology for Your Success!