



# TrisKem International

vUGM

part of the NPL vCARM conference

News Radiopharmacy and other on-going R&D

S. Happel - 24/11/2020

Radiopharmacy  
and  
Nuclear Medicine


Environment and  
Bioassay

Geochemistry  
and  
Metals Separation


Decommissioning



- R&D , TechSupport & QC team: 9 persons
  - 3 radiochemistry PhD, 2 organic chemistry PhD, 4 technicians
- Two R&D labs:
  - Synthesis Lab (new resins and extractants)
    - Incl. grafted resins (silica or polymers), macrocycles,...
    - Extractants
  - Application Lab
    - Preparation of extraction chromatographic resins
    - Resin characterisation and method development
- Counting equipment:
  - ICP-MS, IC, TOC, TGA, IR, moisture analyzer, surface area analyser, particle size and shape analyser, pycnometer
  - Production and packing lab with four 20L reactors
- No handling of radioactivity => R&D cooperation
  - Resin and method development “cold” => R&D partner




- Radionuclide production
  - Separation of radionuclides from irradiated targets
    - » Diagnostics: Zr-89, Cu-64, Ga-68, Ge-68, Ti-44/5, Tc-99m, Sc-43/4, ...
      - **ZR Resin, CU Resin, TK200 Resin, TK400, TK201, TK202, ...**
    - » Therapy: alpha emitters, Lu-177, Tb-161, Cu-67, Sn-117m, Sc-47, ...
      - **TK400, TK200, TBP Resin, CU Resin, TK211/2/3, TK221, ...**
- Quality control
  - **DGA sheets** (functionalized TLC, Ra-223, Ga-68, Pb-212, ...  
=> CVUT Prague)
  - Cartridge based methods
- Decontamination of effluents
- Purification of generator eluates => under development



- ZR Resin
  - Hydroxamate based resin => different from Holland publication!
  - Standard for Zr separation from Y targets
  - Ready to use / no activation
  - Facile Zr elution (avoid 1M oxalic acid)
- Zr-89 production via (p,n) reaction from  $^{nat}Y$  targets
  - High Zr/Y selectivity necessary
  - Alternative e.g. TBP Resin (=> Graves et al.)
- Application for other separations: **Ti/Sc, Ga/Zn, Ge/Ga**

# Zr-89 separation on ZR Resin

Zr separation on 1 mL ZR Resin column



- No activation with acetonitrile!
- Load preferably from 2M HCl
- Rinsing described by Holland may be employed (2M HCl / water)
- Clean Fe removal
- Quantitative Zr elution in 1.5 - 2 mL  $\geq 0.05$ M oxalic acid
- Oxalic acid conc. may be higher (e.g. 1M) => smaller elution volume
- For  $\leq 100$ mg Y often 0.3mL ZR Resin

# Zr-89 separation on TBP Resin




- Method published by Graves et al.
- 400mg Y foils irradiated at 14 MeV (50  $\mu$ A)
- Dissolution in 10 mL conc. HCl
- Separation on 220 mg TBP Resin
- Load from 9.6M HCl, rinse with 20 mL 9.6M HCl
- Zr elution with 1 mL 0.1M HCl
- Zr yield:  $89 \pm 3\%$ , Y decontamination:  $1.5 \times 10^5$
- Zr elution should also be possible with oxalate, citrate, phosphate...
- Other applications of TBP Resin:
  - Sc isotope production from Ca targets (=> presentation EANM'18, Polatom)
  - Sn-117m from Cd targets



Nuclear Medicine and Biology  
Volumes 64–65, September–October 2018, Pages 1-7




Evaluation of a chloride-based  $^{89}\text{Zr}$  isolation strategy using a tributyl phosphate (TBP)-functionalized extraction resin

Stephen A. Graves <sup>a</sup>, Christopher Kutyreff <sup>b</sup>, Kendall E. Barrett <sup>b</sup>, Reinier Hernandez <sup>c</sup>, Paul A. Ellison <sup>b</sup>, Steffen Happel <sup>d</sup>, Eduardo Aluicio-Sarduy <sup>b</sup>, Todd E. Barnhart <sup>b</sup>, Robert J. Nickles <sup>b</sup>, Jonathan W. Engle <sup>b</sup>,   


[Show more](#)

<https://doi.org/10.1016/j.nucmedbio.2018.06.003>


[Get rights and content](#)



# Use of TK400 for Fe/Nb removal





- On-going work following questions about Fe removal
  - On TBP only: Fe and Nb follow Zr
- Removal of Fe & Nb upfront possible using TK400 Resin
- Test with stacked 2 mL TK400/TBP cartridges
  - Load and Rinse: TK400 stacked above TBP
  - Elution: splitting of cartridges and separate elution
    - TBP => ZR
    - TK400 => Fe & Nb
- Best results:
  - TK400 run at 9M HCl => no Zr retention
  - Adjustment of eluate to  $\geq 10M$  HCl
  - Loading of adjusted eluate onto TBP for Zr purification









## ➤ Ti-44 production

- 4g irradiated Sc
- 5 mL Zr Resin
- Ti-44 yield >95%
- 65.2 MBq Ti-44
- $D_f(\text{Sc}): 10^5$

Ti-45 => K. Olguin presentation

Fig. 3. HCl concentration dependency of  $K_d$  for  $^{44}\text{Ti}/^{46}\text{Sc}$  on ZR hydroxamate resin. Fig. 5.  $^{44}\text{Ti}/^{46}\text{Sc}$  elution profile using ZR hydroxamate resin with a load of 4 g of scandium.

## Use of ZR Resin as support in Ti-44/Sc-44 generators




Fig. 1. Schematic concept of a forward/reverse flow radionuclide generator.

- Direct (1 mL ZR) and reverse elution (2 mL ZR)
  - 65 column volumes tested up until publication
  - High Sc yields, max. Ti-44 breakthrough:  $4.1 \cdot 10^{-4}\%$
  - Obtained Sc gave labelling yields > 94%
- Generator been set-up at BNL/SBU => Poster S. Houclier ISRS 2019

# Ge-68 separation from GaNi or GaCo



- Loading from  $\text{HNO}_3$ ,  $\text{HCl}$  or  $\text{H}_2\text{SO}_4$ 
  - Target dissolution in  $\text{HNO}_3$  or  $\text{H}_2\text{SO}_4$  often preferred  $\Rightarrow$   $\text{GeCl}_4$  volatile
- On-going: Cold test on  $>5\text{g}$  GaNi
- First cycle on ZR Resin
  - 2 mL ZR Resin cartridge
  - Load/rinse from 5M  $\text{H}_2\text{SO}_4$
  - High Ge retention/purification from Ga, Ni & Co
  - Elution: 0.1M citric acid (pH 3)
  - Preferably at  $\leq 1$  mL/min
- Second cycle on ZR Resin
  - Preferably 1 mL cartridge
    - Alternative 2 mL cartridge
  - Adjustment of eluate to 5M  $\text{H}_2\text{SO}_4$
  - Load of adjusted eluate through ZR
    - Rinse with 5M  $\text{H}_2\text{SO}_4$
    - Elution with 0.1M citric acid (pH 3)
    - Preferably at  $\leq 1$  mL/min
- Conversion step (citric acid to HCl):
  - Acidification to 9M  $\text{HCl}$ , load onto Guard Resin
    - Alternatives: Prefilter or TK400
    - Advantage Guard Resin: further Ga removal
    - Advantage Prefilter/TK400: works at 6M  $\text{HCl}$  - but no selectivity for Ge over Ga
  - Rinse with 9M  $\text{HCl}$
  - Elution with 0.05M  $\text{HCl}$   $\Rightarrow$  pH control!
- Overall very high Df & chemical yield



# Ga-68(/Ga-67) separation from Zn targets

- **ZR Resin**

- Loading possible from:
  - dilute HNO<sub>3</sub> (**liquid targets**)
  - > 6M HCl (**solid targets**)
- Rinse under loading condition
- Ga separation on ZR Resin
- Elution with ~1.5M HCl

- **Ga conversion step on TK200**

- TK200 load from 1.5M HCl
- Rinse with 1.5M HCl
  - Better pH control of eluate via rinse with NaCl/HCl before elution=> Gagnon et al.
- Elution in 2 – 3 BV water



⇒ **New IAEA TechDoc:**

<https://www-pub.iaea.org/books/IAEABooks/13484/Gallium-68-Cyclotron-Production>

# Cyclotron production of Ga-68

Rodnick et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:25  
https://doi.org/10.1186/s41181-020-00106-9

EJNMMI Radiopharmacy  
and Chemistry

RESEARCH ARTICLE

Open Access



## Cyclotron-based production of $^{68}\text{Ga}$ , $[^{68}\text{Ga}]\text{GaCl}_3$ , and $[^{68}\text{Ga}]\text{Ga-PSMA-11}$ from a liquid target

Melissa E. Rodnick<sup>1</sup>, Carina Sollert<sup>2</sup>, Daniela Stark<sup>3</sup>, Mara Clark<sup>1</sup>, Andrew Katsifis<sup>3</sup>, Brian G. Hockley<sup>1</sup>, D. Christian Parr<sup>2</sup>, Jens Frigell<sup>2</sup>, Bradford D. Henderson<sup>1</sup>, Monica Abghari-Gerst<sup>1</sup>, Morand R. Piert<sup>1</sup>, Michael J. Fulham<sup>4</sup>, Stefan Eberl<sup>5</sup>, Katherine Gagnon<sup>2</sup> and Peter J. H. Scott<sup>1\*</sup>




Fig. 4 Two-column approach for  $^{68}\text{Ga}$  chemical separation

**Table 1** High level schemes of  $[^{68}\text{Ga}]\text{GaCl}_3$  purifications

	Scheme A*	Scheme B
① ZR Load	< 0.1 M $\text{HNO}_3$	
② ZR Wash	15 mL 0.1 M $\text{HNO}_3$	
③ ZR Elution / Trapping on TK200	5–6 mL ~ 1.75 M HCl	
④ TK Wash	–	3.5 mL 2.0M NaCl in 0.13 M HCl
⑤ TK Elution	$\text{H}_2\text{O}$	1–2 mL $\text{H}_2\text{O}$ followed by dilute HCl to formulate

\*Process as reported previously (Nair et al. 2017)

- J. Kumlin et al. (preprint):
  - ZR, LN & TK200 for solid targets

ORIGINAL RESEARCH

## Multi-Curie Production of Gallium-68 on a Biomedical Cyclotron and Automated Radiolabelling of PSMA-11 and DOTATATE

> Helge Thisgaard, Joel Kumlin, Niels Langkjær, Jansen Chua, Brian Hook, Mikael Jensen, Amir Kassaian, Stefan Zeisler, Sogol Borjian, Michael Cross, Paul Schaffer, Johan Hygum Dam

DOI: 10.21203/rs.3.rs-70698/v1 [Download PDF](#)

- High Ga-68 activities
- ARTMS/Odense: 10 Ci production: <https://physicsworld.com/a/cyclotron-based-gallium-68-generator-breaks-production-records/>
- W. Tieu et al.: Use of single TK400 cartridge for solid Zn targets



Nuclear Medicine and Biology  
Volumes 74–75, July–August 2019, Pages 12–18




## Rapid and automated production of $[^{68}\text{Ga}]\text{gallium chloride}$ and $[^{68}\text{Ga}]\text{Ga-DOTA-TATE}$ on a medical cyclotron

William Tieu<sup>a,\*,</sup>, Courtney A. Hollis<sup>a</sup>, Kevin K.W. Kuan<sup>a</sup>, Prab Takhar<sup>a</sup>, Mick Stuckings<sup>b</sup>, Nigel Spooner<sup>b,c</sup>, Mario Malinconico<sup>d</sup>


# Other examples for separations on TK200 (TOPO based)

- Zn separation from Cu




- Zn/Cu separation. Elution study, ICP-MS measurement

- Pt separation from Ir



- Pt/Ir separation. Elution study, ICP-MS measurement

<sup>65</sup>Zn/Cu (spiked with <sup>64</sup>Cu) separation on TK200 resin




- Zn-65 separation. Data kindly provided by Fedor Zhuravlev, DTU

- Other systems under testing
- Radioanalysis: actinide separation
- Tri and tetravalents generally more difficult to elute (Sc,...)

# Cu-61/4 separation on TK201

## ➤ Cu-64 separation from solid Ni targets on TK201:

- Load and rinse at 6M HCl => Ni removal => recovery/recycling
- Co elution with 4 – 5M HCl
- Gagnon et al. use of NaCl/HCl for better pH control of eluate
- Cu elution with 0.5M HCl => Fe and Zn remain retained
- Preferred alternative: Use of TBP (or TK400) upfront for Fe/Ga removal => allows for Cu elution in 0.05M HCl



Svedjehed et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:21  
<https://doi.org/10.1186/s41181-020-00108-7>

(2020) 5:21


EJNMMI Radiopharmacy  
and Chemistry

### RESEARCH ARTICLE

### Open Access

Automated, cassette-based isolation and formulation of high-purity [ $^{61}\text{Cu}$ ]CuCl<sub>2</sub> from solid Ni targets

Johan Svedjehed<sup>1</sup>, Christopher J. Kutyręff<sup>2</sup>, Jonathan W. Engle<sup>2,3</sup> and Katherine Gagnon<sup>1\*</sup>



## Purification of <sup>67</sup>Cu and Recovery of its Irradiated Zn Target

A.J. DeGraffenreid<sup>a</sup>, R. Nidzyn<sup>a</sup>, B. Jenkins<sup>a</sup>, D.E. Wycoff<sup>b</sup>, T.E. Phelps<sup>b</sup>, A. Goldberg<sup>a</sup>, D.G. Medvedev<sup>a</sup>, S.S. Jurisson<sup>b</sup>, C.S. Cutler<sup>a</sup>

<sup>a</sup>Brookhaven National Laboratory, C-AD/MIRP—Upton, NY (USA)

<sup>b</sup>University of Missouri, Department of Chemistry—Columbia, MO (USA)

Poster  
presented at  
ISRS 2017

- 13.7g Zn metal dissolved to give 312 mg ZnCl<sub>2</sub>/mL solution at pH 2
- Loading of 60,6 mL => 18.9g ZnCl<sub>2</sub> onto 2.4g CU Resin column => 8 mL
- Rinse with 80 mL pH2 HCl
- Elution in 2 x 20 mL 6M HCl
- Evaporation to dryness
- Chemical yield ~100%
- Single column D<sub>f</sub> for Zn ~10 000
  - Additional removal indicated
- Ideally further Zn and Co removal
- Original suggestion: AIX

Nuclide	EOB Activity (mCi ± 1σ)	Cu Resin Recovery (%)			
		Load w/ Quant. Transfer	pH 2 HCl Rinse	Acid #1	Acid #2
<sup>64</sup> Cu	4700 ± 200	ND	ND	102	ND
<sup>65</sup> Zn	41.0 ± 0.8	103	ND	0.04	ND
<sup>58</sup> Co	63 ± 1	104	0.04	0.1	0.01

- Produced 143 mCi <sup>67</sup>Cu
- Quantitative recovery of radiocopper
- 99.5% radionuclidic purity—single column
- ICP-OES: 132.9 μg Cu and 1.3 mg Zn
  - Anion exchange column still needed to remove trace Zn
- Specific activity <sup>67</sup>Cu at EOB: 1.07 mCi/μg

### Cu Resin

Robust separation that could shorten the overall processing time to separate co-produced radionuclides and large quantities of Zn from radiocopper  
Cation and anion exchange columns still needed to suitably purify radiocopper

### Alternatives to AIX:

- TK201: preferred option. Cu elution from CU Resin with 6M HCl directly onto TK201, followed by Cu elution from TK201 in dilute acid
- TK200: Cu eluted from CU Resin in 1 - 2M HCl, direct load through TK200 (Zn retained, Cu passes)


- On-going developments radiolanthanides (nca Lu-177, Tb-161)

- Separation of nca Lu-177 from Yb-176 targets (500 – 2000 mg)
- Partial simplification via sequential separation steps
  - ‘Sequential separations’ approach also applicable e.g. to Tb separation
- New resins: TK211/2/3
  - On-bead mix of different extractants for improved selectivity
  - Higher extractant load
  - Small amount of long-chained alcohol and use of inert support containing aromatic groups => aim: improved radiolysis stability
  - 20 - 50  $\mu\text{m}$  beads
    - Originally developed as 15 $\mu\text{m}$  beads => too small for large scale separations
  - Resins also applicable to Horwitz method
    - TK212 & TK221 instead of LN2 and DGA
- Prepacked PP columns under development
  - 29 mL, 53 mL, 150 mL, 375 mL and 750 mL
  - With CoA
- Upscale to 1g and 2g Yb targets under finalisation






- Simplified method for Lu separation from 500 mg Yb – TK211/2 & TK221




- Partial simplification via sequential separation step (direct load from TK212 onto TK211 for polish)

# Lu separation from 500 mg Yb - TK212/TK221/TK212/TK211




- Large tailing due to high Yb content
- Improved separation through use of 1.25M HNO<sub>3</sub> / 10% EtOH (v/v) instead of 1.3M
- Additional benefit from use of EtOH => improved radiolysis stability
- Online separation: switch at start of 'Lu fraction' => ideally radiation detector driven
- General aim: ~95% Lu recovery, ~80% Yb removal

# Lu separation from 500 mg Yb - TK212/TK221/TK212/TK211




- 2<sup>nd</sup> separation step on smaller TK212 (53 mL) after TK221 for conversion from high HNO<sub>3</sub> to dilute HCl => 5g TK221
- Separation with e.g. 1.25M HNO<sub>3</sub> (with or without 10% EtOH)
- Generally aim: Lu >90% recovery, ~99% Yb removal
- Direct loading of obtained Lu fraction onto TK211 Resin
  - Alternatively TK221/TK212 according to Horwitz et al.

# Lu separation from 500 mg Yb - TK212/TK221/TK212/TK211




- Overall Lu recovery of process in the order of 85% => higher if more Yb allowed in final Lu
- Low remaining Yb
- Flow rates may still be optimized
- Final step: concentration/conversion to  $\leq 0.05\text{M}$  HCl on TK221, nitrate removal via A8

# Simplified method for Tb separation from 500 mg Gd – TK211/2 & TK221




# Tb separation from 500 mg Gd targets

- « Sequential separation » TK212 => TK211
- Initial separation on TK212 (150 mL column)
- Fine tuning ongoing (e.g. adjustment of eluents to 10% EtOH)
- Separation easier than Lu/Yb
- Polishing via direct load onto TK211 (29 mL)



# Tb polish on TK211





- Gd breakthrough during load & rinse with 0.5M HNO<sub>3</sub>
- Method optimisation on-going => yields and purity very high
- Alternatively: Tb elution with 1M HNO<sub>3</sub>  
=> smaller elution volume / better Tb retention on TK221
- Conversion to dilute HCl via TK221, A8 for nitrate removal
- Next steps => Upscale to 1g and 2g of Gd

# On-going developments Lu-177


## – new TK221 Resin

- DGA well suited for ‘conversion’ and purification (Ca, Al, Fe,... removal)
  - Convert Lu from high nitric acid to dilute HCl
- Elution of heavy lanthanides in 0.05M HCl needs elevated volumes
  - smaller volume preferred => higher activity concentration
- Optimisation of DGA Resin => new TK221 Resin
  - Lu eluted in smaller volume
  - Radioanalytical applications: higher U retention compared to DGA





- TK202 Resin
  - Based on Polyethylene Glycol (PEG) grafted on inert support
    - Presentation Izabela Cieszykowska
  - Poster presentation at ISTR2019, Vienna, 28/10/19 – 01/11/19
    - Tc recovery > 90% for 6 – 8g Mo per g of TK202
    - Tc recovery > 80% for 12g Mo per g of TK202
  - On-going tests:
    - Tc separation from larger Mo targets
    - Use of carbonate rinse for better pH control of eluate
    - Use of an aluminium oxide 'guard column' for trace Mo removal
    - Decommissioning => Tc determination in decommissioning samples after alkaline fusion => presentation A. Bombard
- SE Resin
  - Functional group covalently bound onto polymer support
  - Piazselenol chemistry (retention as Se(IV))
  - Load from e.g. 2M HCl, elution with 0.1M NaOH/0.1 H<sub>2</sub>O<sub>2</sub>
  - Se-72 separation and generator support, Se-79, Se isotope ratios => First beta testing on-going



# Some other on-going projects

- Decontamination
  - PAN based resins,...
- Cs/Rb separation (TK300)
  - Beta testing on-going
  - Calixarene based
  - Probably 2 resins
- Ra Resin
  - Macrocycle based
  - Probably 2 resins
- Li Resin
  - Macrocycle based
- Rapid tests
  - Test sticks => Uni Southampton
  - DGA Sheets (2D TLC)
  - Extractive membranes
- At separation (TK400,...)
- Range of PSm based resins
  - TK-ElScint (p.ex. TK-TcScint)
- DGT (Diffusive Gradients in Thin Films) => 'bio-availability'
- Functionalised polymers & silicates, ...
  - e.g. DO-DGA, DE-DGA, macrocycles,...
- Microfluidics
- Improvement of radiolysis stability
- Other types of supports

- If any of these R&D topics are of interest – we are always looking for research collaborations or common projects so please contact us! => [shappel@triskem.fr](mailto:shappel@triskem.fr)
- In case you have developed a new product that might be of interest please get in touch to discuss
- Please never hesitate to contact us if you have any technical questions or facing the need to develop new separations

Thank you for your attention!



## SUBSCRIBE TO OUR NEWSLETTER

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Infos [here](#)

