## ROTI®Store Hefe-Cryoröhrchen

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For preparing cryogenic cultures

#### A. Product description

ROTI®Store is a product for preparing cryogenic cultures. It enables safe storage and recultivation of microorganisms. ROTI®Store yeast cryo-vials have been specially developed for the preparation of yeast cryogenic cultures. The special freezing medium facilitates freezing and recultivation of different yeast stocks such as *Saccharomyces cerevisiae* or *Candida albicans*. They may be stored at -20 °C for several weeks. For long term-storage we recommend a minimum of -80 °C. In order to recultivate simply remove a bead and bring it into contact with desired culture medium.


Cryo-conversation has been successfully tested for *S. cerevisiae*, *C. albicans*, *C. sake*, *C. tropicalis*, *Zygosaccharomyces rouxii*, *Dekkera bruxellensis* and *Malassezia furfur*

#### B. Pack contents and storage

ROTI®Store contains 50 cryo-vials with approx. 25 beads per vial in a special yeast cryo-medium. The non-inoculated ROTI®Store yeast cryo-vials can be stored between +4 °C to +25 °C (optimal storage: +4 to +8 °C). They should not be constantly exposed to sunlight.

#### C. Inoculation

- You will require a fresh overnight culture or colonies resuspended in culture medium with an OD<sub>660</sub> of at least 2.\*
- Check the vial prior to use for possible contamination (cloudiness). Please label the cryo-vial with a waterproof marker (e.g. ROTILABO®-LabMarker, Roth Art. No. 3083.2) and add approx. 0.5 ml culture solution.
- Close the vial tightly and shake it carefully several times.
- Remove the supernatant above the pellet with a sterile pipette (Roth Art.No. EA69.1) and dispose of it correctly.
- The inoculated and tightly closed ROTI®Store yeast cryo-vial can now be frozen and stored at -20 °C for several weeks or at -80 °C (or lower for years).

\*Direct inoculation of ROTI®Store cryo-vials using picked colonies is possible, it does result, however, in much less efficient recultivation rates. In these cases, we recommend to resuspend the colony well in the cryo medium. Slightly agitate the tube 15-20 mins. at 4 °C and remove the supernatant only then.

#### D. Recultivation

- Remove the ROTI®Store yeast cryo-vial from the deep-freezer and transport it in a cold cryo-block (e.g. Cooling box ROTILABO® blue, Roth Art.No. K088.1). The vial should not be allowed to thaw.
- Remove a bead under sterile conditions with a sterile pincette or inoculation needle (e.g. Roth Art.No. EA88.1) and then close the cryo-vial.
- The bead may now be added directly to the liquid medium or rolled out onto a solid culture medium.

Ø: 12 mm  
Length: 49 mm

#### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com jh 12/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

#### ROTI®Store yeast cryo-vials

**X983.1** 50 vials with 25 beads per vial