

ABG II VERSUS ANATO

Stryker Confidential. Do not distribute.



ANATO STEM

- ONE global anatomic stem
- Create a power brand
- Combine the benefits of all stems in Stryker's Anatomic Stem History
 - PCA
 - Citation
 - ABG II




Introduction


- Anatomic stems have been used for over 30 years with clinical success as seen with the PCA and Citation families¹ as well as the ABG and ABG II monolithic families².
- The anatomic stem philosophy is intended to transfer loads more circumferentially through implant geometries which may better match the metaphyseal bone.

1-Bojescul, et. al, JBJS-A, 2003 2-Nourisat, et. al, JOA, 2003



ABG Body Rationale





Anatomic load transfer occurs in the metaphysic

The ABG design aims to re-establish anatomic load transfer

Orthopaedics

History of ABG

ABG I

1985-1989: ABG development 1989: 1st ABG THP implantation ABG II

1993-1995: ABG II development

1997: 1st ABG II stem implantation

ANATO

2010: Anatomic stem

2014: Release next generation anatomic stem







Orthopaedics


PCA and Citation stems

1983 PCA 1990 PCA 1996 **Citation AT E-Series**

Orthopaedics

COMPARE

ANATO



Orthopaedics

ANATO is a new concept in stem design, combining the benefits of a fit-and-fill device with the procedural flexibility of a shorter stem

ANATO incorporates the best features of Stryker's successful ABGII and Citation hips, and its revolutionary design has been validated by SOMA:

- 8 sizes
- Anteverted and neutral neck options
- Anatomic Shaping

Pronounced as in <u>ANATO</u>my

– /əˈnato/





Instrumentation

- Works in conjunction with existing ABG and Accolade TMZF and II instruments
- Works with DA instruments (new and old)
- New Instrumentation
 - Shortened broaches- size 6, 7 and 8
 - Stem Extractor
 - MIS Neck Trials
 - Case and Trays



THANK YOU