

## **FEATURES AND BENEFITS**

- Easy-to-read contents indicator Visible through the clear window and color coded
- Audible alarms Sounds at 325 PSI / 120 liters
- Visual alarms LED lights flash red to indicate low content
- Lightweight construction 40% lighter than steel for increased mobility
- Improve patient satisfaction Reduce potential therapy interruption
- Minimize potential run-outs Eliminate emergency inventory requests
- Ready to use CGA compliant for use with appropriate delivery device or regulator
- Convenient Content indicator is permanently attached and system is maintained by Praxair

## PRAXAIR HEALTHCARE SERVICES


## *Medipure*™ Alert N<sub>2</sub>O

## Nitrous Oxide Cylinder with Content Indicator

As a caregiver, you know that if nitrous oxide runs out unexpectedly, a medical procedure could be interrupted.

The physical properties of nitrous oxide are different than most other medical gases. In a cylinder, nitrous oxide under pressure and at room temperature, exists as a liquid. When in use, the pressure remains constant until the last of the liquid vaporizes. The pressure then drops quickly, giving you only minutes to change to a full cylinder.

With *Medipure* Alert  $N_2O$ , Praxair's portable nitrous oxide cylinder, advance notice is given. When the contents are low, it will be visible through the contents indicator and you will hear an alert.



\* approximate based on withdrawal rate of 2 LPM

Call Praxair to Schedule an Evaluation of Your Cylinder Needs Today at 1.800.PRAXAIR

© Copyright 2016 Praxair Technology Inc. All rights reserved.

Praxair, Medipure, and the Flowing Airstream design are trademarks or registered trademarks of Praxair Technology, Inc. in the United States and/or other countries. \*Refer to Operating Instructions and Product Label for details

The information contained herein is offered for use by technically qualified personnel at their discretion and risk without warranty of any kind.